
## Epistaxis Catheter Cathéter à épistaxis Epistaxe-Katheter Catéter para epistaxis Catetere per epistassi Epistaxiskatheter Kateter for neseblødning

## Instructions for Use Mode d'emploi Gebrauchsanweisung Modo de empleo Istruzioni per l'uso Gebruiksaanwijzing Bruksanvisning

Read instructions before use.

Rx Only

For US Audiences Only

### DESCRIPTION

The EpiStax is a two-balloon silicone catheter with two independently inflatable silicone balloons designed to control intra-nasal hemorrhage. The smaller balloon marked "10cc" is designed to control bleeding in the posterior chamber, while the larger balloon marked "30cc", controls bleeding in the anterior chamber. The airway tube permits nasal breathing and access for suctioning.

### INTENDED USE

The balloon nasal catheters are intended to be used to control nasal epistaxis (bleeding).

### CAUTIONS

- This device is for **SINGLE USE ONLY**. Do **NOT** re-sterilize or re-use.
- Inspect for damage prior to use. Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- Do not use petroleum based lubricants with this device as it will seriously affect the performance of the balloons potentially causing them to burst.
- Sterile saline is not be used to inflate the device as the salts used may crystallize and cause a blockage.
- This device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

**CAUTION: US Federal Law restricts sales and use to or on the order of a Physician.**

### INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA or the competent authority of the Member State in which the user and/or patient is established.

### INSTRUCTIONS FOR USE

#### Test Inflation

- Attach a syringe to the valve marked "10cc", and test-inflate the distal balloon with 5 to 10ml of air. Deflate balloon. Repeat for valve marked "30cc", test-inflate with 15 to 25 ml of air. Deflate balloon. If either balloon does not inflate, use a new device and call Summit Medical Customer Service.

#### Insertion

- Lubricate catheter using sterile water or water-based lubricant.
- After removal of nasal blood clots, insert the catheter so that the beveled tip is in the nasopharynx.
- Fill a 20cc syringe with sterile water.
- Insert tip of syringe into valve marked "10cc" to inflate posterior balloon. Inflate balloon with fluid from syringe. Inflation volumes of 4cc to 8cc are usually adequate.

**CAUTION: Care must be taken not to over-inflate balloons.**

- Remove the syringe and gently pull the catheter anteriorly to seat the posterior balloon.
- Insert tip of syringe into valve marked "30cc" to inflate anterior balloon. Inflate balloon with fluid from syringe. Inflation volumes of 10cc to 25cc are usually adequate.
- Check balloon pressures frequently and relieve pressure when indicated. Suction as necessary through airway tube using lubricated suction catheter.

#### Removal

- After 2-5 days remove catheter.
- Insert syringe and evacuate each balloon completely. After balloons are completely deflated, slowly remove catheter and discard per facility protocol.

### STERILIZATION

The device is a **SINGLE USE ONLY** device supplied sterile and ready for use. Sterilization is by Ethylene Oxide (EO).

Single-use, one-patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of the patient.

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