The Denver Splint®

**SYMBOL REFERENCE KEY / RÉFÉRENCE DE SYMBOLE / SYMBOLLEGENDE /
CHAVE DE REFERÊNCIA DE SÍMBOLO / POUZITÉ SYMBOL / CLAVE DE REFERENCIA DE
LOS SÍMBOLOS / LEGENDA DEI SIMBOLI / SYMBOLREFERENSNSYKEL / LEGENDA**

i Consult instructions for use / Consulter le mode d'emploi / Gebrauchsanweisung
beachten / Consultar as instruções de utilização / Prečítete si návod k použití /
Consultar las instrucciones de uso / Consultare le istruzioni per l'uso / Brugsanvisning/
Raadpleg de gebruiksaanwijzing

Quantity / Quantité / Meng / Quantitate / Množství / Cantidad / Quantità / Antal / Aantal

REF Catalogue number / Numéro de référence / Bestellnummer / Número de catálogo /
Katalogový číslo / Número de catálogo / Numero catalogo / Katalognummer /
Catalogusnummer

LOT Date of manufacture / Date de fabrication / Herstellungsdatum / Data de fabrico /
Datum výroby / Fecha de fabricación / Data di fabbricazione / Tillverkningsdatum /
Datum van vervaardiging

NON STERILE Non-sterile / Non stérile / Nicht steril / Não estéril / Nesteril /
Non sterile / Ej steril / Niet-steriel /

**Do not use if package is damaged / Ne pas utiliser si l'emballage est endommagé /
Beschädigte Verpackung nicht verwenden / Não utilizar se a embalagem estiver
danificada / Nepoužívajte, pokud je obal poškozen / No utilzar si el envase está dañado /
Non utilizzare se il pacchetto è danneggiato / Far inte användas om förpackningen är
skadad / Niet gebruiken als de verpakking beschadigd is**

**Caution / Attention / Vorsicht / Atenção / Upozornění / Atención /
Försiktighet / Opgelet**

X Do not reuse / Ne pas réutiliser / Nicht wiederverwenden / Não reutilizar /
Nepoužívajte opakované / No reutilizar / Ateranvänd inte /
Niet hergebruiken

**Use by date / Date limite d'utilisation / Verwendbar bis / Prazo de validade /
Použijte do data / Fecha de caducidad / Data di scadenza / Anvärds före datum /
Houdbaarheidsdatum**

**Manufacturer / Fabricant / Hersteller / Fabricante / Výrobce / Fabricante /
Produttore / Tillverkare / Fabrikant**

**Country of Manufacture (Made in the US) / Pays de fabrication (Fabriqué aux États-Unis) /
Herstellungsland (Hergestellt in den USA) / País de fabricación (Fabricado nos EUA) /
Země výroby (vyrobeno v USA) / País de fabricación (Hecho en los EE. UU.) /
Paese di produzione (fabbricato negli Stati Uniti) / Tillverkningsland (Tillverkad i USA) /
Land van vervaardiging (gemaakt in de Verenigde Staten) /**

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NOTIFICATION DE INCIDENTES GRAVES

INTENDED USE/INDICATIONS FOR USE
The Denver Splint™ is intended to provide external support after rhinoplasty or nasal fractures.

The Denver Splint is indicated to be used postoperatively on rhinoplasty and nasal fracture patients to provide cushioned pressure to reduce edema and to provide external stabilization of bony fragments.

INTENDED USER
The Denver Splint is intended to be used by healthcare professionals, such as, plastic and cosmetic surgeons, facial trauma surgeons, and otolaryngologists.

WARNINGS AND PRECAUTIONS

- Do not apply splint directly to skin - may cause mild sensitization in certain sensitive individuals. Apply paper tape first.
- For Single Use Only. Do not reuse. Reuse of device could result in infection, contamination or device failure which could lead to patient harm.
- Do not resterilize.
- Do not use if package is damaged or unintentionally opened prior to use, or if dorsal pad is compressed.
- Do not use in known cases of severe soft tissue injuries that compromise blood supply to the skin.
- Remove splint prior to X-ray, CT or MRI imaging. Splints may interfere with image quality.

INSTRUCTIONS OF APPLICATION

PREPARATION

- *For Thermoplastic Splints: Preheat unit or heat water to 160-180°F (71-82°C)
- Cleanse and dry skin of nose.
 - Use alcohol sponge to wipe nose. Dry skin of nose after wiping.
 - Express edema manually from skin.
 - Apply a uniform coating of SKIN-PREP® to nasal skin and allow to dry (approximately 30 seconds). If an area intended to be covered is missed, wait until SKIN-PREP has dried then reapply.
 - Express edema again.
 - Apply overlapping paper tape to nose.

APPLICATION

Malleable Splint

- Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. **DO NOT TOUCH SKIN WITH DORSAL PAD.**
 - Remove paper backing from adhesive side of splint.
 - Shape metal. Apply adhesive side of splint to tape on the nose. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**
- Thermoplastic Splint (Pre-Cut)**
- Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. **DO NOT TOUCH SKIN WITH DORSAL PAD.**
 - Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds.
 - Remove splint from water using forceps or a tongue depressor to avoid burning fingers.
 - Blot briefly on towel. The splint will cool but remain workable for approximately one minute.
 - Apply softened splint and mold over the nose.
 - If further molding is required, remove splint, resoften in hot water and reapply.
 - When in place, the splint should touch only the underlying tape and should not touch the skin. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**

Thermoplastic Splint (Custom Cut)

- Design desired shape using gauze or paper towel. Check fit to patient.
- Place shape over Thermoplastic Sheet/Flat and cut splint to shape with bandage scissors.
- Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds.
- Remove splint from water using forceps or a tongue depressor to avoid burning fingers.
- Blot briefly on towel. The splint will cool but remain workable for approximately one minute.
- Apply softened splint and mold over the nose.
- If further molding is required, remove splint, resoften in hot water and reapply.
- When in place, the splint should touch only the underlying tape and should not touch the skin **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**

REMOVAL

Remove splint in seven to ten days. Recommend removing by spreading bayonet forceps (dull) between skin and tape. Do not tent up skin with removal. Dispose of device as clinical waste in line with local protocol.

SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the Member State in which the user and/or patient is established.

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