SYMBOL REFERENCE TABLE / TABLEAU DE RÉFÉRENCE DES SYMBOLES / TABELLE: SYMBOLLEGENDE / TABELA DE REFERÊNCIA DE SÍMBOLOS / CLAVE DE REFERENCIA DE LOS SÍMBOLOS / LEGENDA DEI SIMBOLI / LEGENDA	
	Consult instructions for use / Consulter le mode d'emploi / Gebrauchsanweisung beachten / Consultar as instruções de utilização / Consultar las instrucciones de uso / Consultare le istruzioni per l'uso / Raadpleeg de gebruiksaanwijzing/
QTY 	Quantity / Quantité / Menge / Quantidade / Cantidad / Quantità / Aantal
REF 	Catalogue number / Numéro de référence / Bestellnummer / Número de catálogo / Número de referencia / Numero di catalogo / Catalogusnummer
	Date of manufacture / Date de fabrication / Herstellungsdatum / Data de fabrica / Fecha de fabricación / Data di fabbricazione / Datum van vervaardiging
LOT 	Batch code / Code de lot / Chargennummer / Código de lote / Código de lote / Codice lotto/ Batchcode
	Do not resterilize / Ne pas restériliser / Nicht erneut sterilisieren / Não reesterilizar / No reesterilizar / Non ri-sterilizzare / Niet opnieuw steriliseren
STERILE R 	Sterilized using irradiation / Sterilisé par irradiation / Durch Bestrahlung sterilisiert / Esterilizado por irradiação / Esterilizado mediante irradiação / Sterilizzato mediante irraggiamento/ Gesteriliseerd door bestraling
STERILEEO 	Sterilized by Ethylene Oxide / Sterilisé à l'oxyde d'éthylène / Mit Ethylenoxid sterilisiert / Esterilizado por óxido de etileno / Esterilizado con óxido de etileno / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide
	Non-sterile / Non stérile / Nicht steril / Não estéril / No estéril / Non sterile / Niet-steriel
	Do not use if package is damaged / Ne pas utiliser si le conditionnement est endommagé / Beschädigte Verpackung nicht verwenden / Não utilizar se a embalagem estiver danificada / No utilizar si el envase está dañado / Non utilizzare in caso di imballaggio danneggiato / Niet gebruiken als de verpakking beschadigd is
	Caution / Attention / Vorsicht / Atenção / Atención / Attenzione / Opgelet
	Do not reuse / Ne pas réutiliser / Nicht wiederverwenden / Não reutilizar / No reutilizar / Non riutilizzare / Niet hergebruiken
	Use by date / Date limite d'utilisation / Verwendbar bis / Prazo de validade / Fecha de caducidad / Data di scadenza / Houdbaarheidsdatum
	Manufacturer / Fabricant / Hersteller / Fabricante / Fabricante / Produttore / Fabrikant
	Keep dry / Garder au sec / Vor Feuchtigkeit schützen / Mantar seco / Mantener seco / Mantere asciutto / Droog houden
	Keep away from sunlight / Garder à l'abri du soleil / Vor Sonneneinstrahlung schützen / Mantar afastado da luz solar / Mantener alejado de la luz solar / Tenere al riparo dalla luce / Uit de buurt van zonlicht houden
	Upper limit of temperature / Température maximale autorisée / Maximal zulässige Temperatur / Limite superior de temperatura / Limite superior de temperatura / Limite superior della temperatura / Bovengrens van temperatuur
	Flammable / Inflammable / Entzündbar / Inflammêl / Inflammable / Inflammabile / Brandbaar
	Health hazard/hazardous to the ozone layer / Dangereux pour la santé/Dangereux pour la couche d'ozone / Gesundheitsrisiko/Schädlich für die Ozonschicht / Perigoso para a saúde / perigoso para a camada de ozono / Peligroso para la salud/peligroso para la capa de ozono / Pericoloso per la salute e lo strato di ozono / Gezondheidsrisico/schadelijk voor de ozonlaag
	Caution: US Federal law restricts this device to sale by or on the order of a physician / Attention: La loi fédérale américaine n'autorise la vente de ce dispositif qu'à un médecin ou sur ordonnance d'un médecin / Vorsicht: Das US-Bundesgesetz beschränkt den Verkauf dieses Produkts auf Ärzte oder auf Anordnung eines Arztes. / Atenção: A lei federal dos EUA restringe a venda deste dispositivo à classe médica ou por indicação de um médico / Atención: De acuerdo con la legislación Federal de los EE.UU., este dispositivo solo puede ser vendido por médicos o por prescripción médica / Attenzione: La legge federale degli Stati Uniti limita la vendita di questo dispositivo solo su prescrizione medica / Opgelet: De federale wetgeving van de Verenigde Staten beperkt de verkoop van dit hulpmiddel tot de verkoper door of op voorschrift van een arts
EC REP 	Authorized representative in the European Community / European Union / Représentant autorisé établi au sein de la Communauté européenne/Union européenne / Bevollmächtigte Vertretung in der Europäischen Gemeinschaft/ Europäischen Union / Representante autorizado na Comunidade Europeia / União Europeia / Representante autorizado en la Comunidad Europea / Unión Europea / Rappresentante autorizzato nella Comunità europea/Unione europea / Officiële vertegenwoordiger in de Europese Gemeenschap/Europese Unie
MD 	Medical Device / Dispositif médical / Medizinprodukt / Dispositivo Médico / Dispositivo Médico / Dispositivo medico / Medisch hulpmiddel
UDI 	Unique Device Identifier / Identifiant unique du dispositif (UDI)/ Eindeutige Geräteknennung / Identificador Único do Dispositivo / Identificador único de dispositivo / ID univoco del dispositivo / Unieke apparaat-id
	Not made with natural rubber latex / Ne contient pas de latex de caoutchouc naturel / Ohne Naturkautschuklatex hergestellt / Não é fabricado com látex de borracha natural / No está fabricado con látex de caucho natural. / Non fabbricato con lattice di gomma naturale/ Niet vervaardigd met natuurlijk rubber latex

		AN INNOVIA MEDICAL COMPANY
EC REP		

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



*SIZES AVAILABLE: Large, Small/Medium, Petite & Mini
Patent #’s 4,153,051 4,213,452
Denver Splint® is a registered trademark of Summit Medical, an Innovia Medical Company. SKIN-PREP® is a registered trademark of Smith & Nephew. VELCRO® is a registered trademark of Velcro Industries B.V. B. A. **SINGLE-USE**

INTENDED USE/INDICATIONS FOR USE
The Denver Splint® is intended to provide external support after rhinoplasty or nasal fractures. The Denver Splint is indicated to be used postoperatively on rhinoplasty and nasal fracture patients to provide cushioned pressure to reduce edema and to provide external stabilization of bony fragments.

INTENDED USER
The Denver Splint is intended to be used by healthcare professionals, such as, plastic and cosmetic surgeons, facial trauma surgeons, and otolaryngologists.

CONTRAINDICATIONS
SKIN-PREP® for use on intact skin only.

WARNINGS AND PRECAUTIONS

- Do not apply splint directly to skin- may cause mild sensitization in certain sensitive individuals. Apply paper tape first.
- For Single Use Only. Do not reuse. Reuse of device could result in infection, contamination or device failure which could lead to patient harm.
- Do not resterilize.
- Do not use if package is damaged or unintentionally opened prior to use, or if dorsal pad is compressed.
- Do not use in known cases of severe soft tissue injuries that compromise blood supply to the skin.
- Remove splint prior to X-ray, CT or MRI imaging. Splints may interfere with image quality.
- SKIN-PREP – Warning: causes serious eye irritation. Avoid contact with eyes. In the case of accidental contact, flush eyes well with water. Should redness or other signs of irritation appear, discontinue use. Warning: may cause drowsiness and dizziness. Danger: Highly flammable liquid and vapor. Storage: Keep dry and away from sunlight, upper temperature limit of 25°C/77°F.
- Protecto Tape - Skin should be clean, dry, and free of oils to assure good adhesion.

APPLICATION INSTRUCTIONS

PREPARATION
**For series 4000 and Thermoplastic Sheets/Flats: Preheat unit or heat water to 160-180°F (71-82°C)*

- Cleanse and dry skin of nose.
- Use alcohol sponge to wipe nose. Dry skin of nose after wiping.
- Express edema manually from skin.
 - If an area intended to be covered is missed, wait until SKIN-PREP has dried then reapply.
- Express edema again.
- Apply overlapping paper tape to nose.

APPLICATION
Series 1500 Kit & Non-Kit
<ol style="list-style-type: none">Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. DO NOT TOUCH SKIN WITH DORSAL PAD. Remove paper backing from Velcro® (i.e. adhesive segment) of splint and apply to paper tape on nose. DO NOT TOUCH SKIN WITH SPLINT. Narrow end of splint should point toward nasofrontal angle. Shape metal segment of splint and attach to Velcro. Mold by bending wings of splint toward face to make Velcro contact. DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.
Series 2000 Kit & Non-Kit
<ol style="list-style-type: none">Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. DO NOT TOUCH SKIN WITH DORSAL PAD. Remove paper backing from adhesive side of splint. Shape metal. Apply adhesive side of splint to tape on the nose. DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.
Series 4000 Kit & Non-Kit
<ol style="list-style-type: none">Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. DO NOT TOUCH SKIN WITH DORSAL PAD. Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds. Remove splint from water using forceps or a tongue depressor to avoid burning fingers. Blot briefly on towel. The splint will cool but remain workable for approximately one minute. Apply softened splint and mold over the nose. If further molding is required, remove splint, resoften in hot water and reapply. When in place, the splint should touch only the underlying tape and should not touch the skin. DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.
Series 5000 Kit & Non-Kit
<ol style="list-style-type: none">Remove paper backing from adhesive segment of Velcro® Loop Splint Segment and apply to paper tape on nose. Narrow end of Velcro Loop Splint Segment should point toward nasofrontal angle. Shape metal segment of splint and attach hook Velcro to loop splint segment. Attached dorsal pad will be vertical on dorsum of the nose. Mold by bending wings of splint toward the face to make Velcro contact. DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.
Thermoplastic
<ol style="list-style-type: none">Design desired shape using gauze or paper towel. Check fit to patient. Place shape over Thermoplastic Sheet/Flat and cut splint to shape with bandage scissors. Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds. Remove splint from water using forceps or a tongue depressor to avoid burning fingers. Blot briefly on towel. The splint will cool but remain workable for approximately one minute. Apply softened splint and mold over the nose. If further molding is required, remove splint, resoften in hot water and reapply. When in place, the splint should touch only the underlying tape and should not touch the skin. DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.

REMOVAL
Remove splint in seven to ten days. Recommend removing by spreading bayonet forceps (dull) between skin and tape. Do not tent up skin with removal. Dispose of device as clinical waste in line with local protocol.

SERIOUS INCIDENT REPORTING
Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the Member State in which the user and/or patient is established.

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