**INTENDED USE/INDICATIONS FOR USE**  
The Denver Splint® is intended to provide external support after rhinoplasty or nasal fractures.  
The Denver Splint is indicated to be used postoperatively on rhinoplasty and nasal fracture patients to provide cushioned pressure to reduce edema and to provide external stabilization of bony fragments.

**INTENDED USER**  
The Denver Splint is intended to be used by healthcare professionals, such as, plastic and cosmetic surgeons, facial trauma surgeons, and otolaryngologists.

**CONTRAINdications**  
SKIN-PREP® for use on intact skin only.

#### WARNINGS AND PRECAUTIONS

- Do not apply splint directly to skin- may cause mild sensitization in certain sensitive individuals. Apply paper tape first.
- For Single Use Only. Do not reuse. Reuse of device could result in infection, contamination or device failure which could lead to patient harm.
- Do not sterilize.
- Do not use if package is damaged or unintentionally opened prior to use, or if dorsal pad is compressed.
- Do not use in known cases of severe soft tissue injuries that compromise blood supply to the skin.
- Remove splint prior to X-ray, CT or MRI imaging. Splints may interfere with image quality.
- SKIN-PREP - Warning: causes serious eye irritation. Avoid contact with eyes. In the case of accidental contact, flush eyes well with water. Should redness or other signs of irritation appear, discontinue use. Warning: may cause drowsiness and dizziness. Danger: Highly flammable liquid and vapor. Storage: Keep dry and away from sunlight, upper temperature limit of 25°C/77°F.
- Protecto Tape - Skin should be clean, dry, and free of oils to assure good adhesion.

#### APPLICATION INSTRUCTIONS

##### PREPARATION

- \*For series 4000 and Thermoplastic Sheets/Flats: Preheat unit or heat water to 160-180°F (71-82°C)
- Cleanse and dry skin of nose.
  - Use alcohol sponge to wipe nose. Dry skin of nose after wiping.
  - Express edema manually from skin.
  - Apply a uniform coating of SKIN-PREP® to nasal skin and allow to dry (approximately 30 seconds). If an area intended to be covered is missed, wait until SKIN-PREP has dried then reapply.
  - Express edema again.
  - Apply overlapping paper tape to nose.

##### APPLICATION

###### Series 1500 Kit & Non-Kit

- Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. **DO NOT TOUCH SKIN WITH DORSAL PAD.**
- Remove paper backing from Velcro® (i.e. adhesive segment) of splint and apply to paper tape on nose. **DO NOT TOUCH SKIN WITH SPLINT.** Narrow end of splint should point toward nasofrontal angle.
- Shape metal segment of splint and attach to Velcro.
- Mold by bending wings of splint toward face to make Velcro contact. **DO NOT PINCH.** **DO NOT TOUCH SKIN WITH SPLINT.**

###### Series 2000 Kit & Non-Kit

- Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. **DO NOT TOUCH SKIN WITH DORSAL PAD.**
- Remove paper backing from adhesive side of splint.
- Shape metal. Apply adhesive side of splint to tape on the nose. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**

###### Series 4000 Kit & Non-Kit

- Trim dorsal pad length to fit splint size. Apply center dorsal pad vertically to tape on the dorsum of the nose. **DO NOT TOUCH SKIN WITH DORSAL PAD.**
- Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds.
- Remove splint from water using forceps or a tongue depressor to avoid burning fingers.
- Blot briefly on towel. The splint will cool but remain workable for approximately one minute.
- Apply softened splint and mold over the nose.
- If further molding is required, remove splint; resofen in hot water and reapply.
- When in place, the splint should touch only the underlying tape and should not touch the skin. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**

###### Series 5000 Kit & Non-Kit

- Remove paper backing from adhesive segment of Velcro® Loop Split Segment and apply to paper tape on nose. Narrow end of Velcro Loop Split Segment should point toward nasofrontal angle.
- Shape metal segment of splint and attach hook Velcro to loop split segment. Attached dorsal pad will be vertical on dorsum of the nose. Mold by bending wings of splint toward the face to make Velcro contact. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**
- Design desired shape using gauze or paper towel. Check fit to patient.
- Place shape over Thermoplastic Sheet/Flat and cut splint to shape with bandage scissors.
- Soften splint in 160-180°F (71-82°C) water for 30 to 60 seconds.
- Remove splint from water using forceps or a tongue depressor to avoid burning fingers.
- Blot briefly on towel. The splint will cool but remain workable for approximately one minute.
- If further molding is required, remove splint; resofen in hot water and reapply.
- When in place, the splint should touch only the underlying tape and should not touch the skin. **DO NOT PINCH. DO NOT TOUCH SKIN WITH SPLINT.**

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