Symbol Reference Key / Définition des symboles / Symbollegende / Definición de los símbolos / Legenda dei simboli / Legenda / Symbolreferansenøkkel / Ključ referenčnih symbolů / Legendas dos símbolos

	Consult instructions for use / Consulter le mode d'emploi / Gebrauchsanweisung beachten / Consulte el modo de empleo / Consultare le istruzioni per l'uso / Raadpleeg de gebruiksaanwijzing / Se bruksanvisning / Прочтите si návod k použití / Consultar as instruções de utilização
	Quantity / Quantité / Stückzahl / Cantidad / Quantità / Aantal / Menge / Množství / Quantidade
MD	Medical Device / Dispositif médical / Medizinprodukt / Producto sanitario / Dispositivo médico / Medisch hulpmiddel / Medisinsk utstyr / Zdravotnický prostředek / Dispositivo médico
REF	Catalog Number / Numéro de catalogue / Artikelnummer / Número de catálogo / Numero di catalogo / Catalogusnummer / Katalognummer / Katalogové číslo / Número de catálogo
LOT	Lot Number / Numéro de lot / Chargennummer / Número de lote / Numero di lotto / Partijnummer / Partinummer / Číslo šarže / Número do Lote
	Use By / Utiliser avant / Verwendbar bis / Fecha de caducidad / Data di scadenza / Te gebruiken voor / Brukes av / Datum spotřeby / Usar até
STERILE EO	Sterilized by Ethylene Oxide / Stérilisé à l'oxyde d'éthylène / Mit Ethylenoxid sterilisiert / Esterilizado con óxido de etileno / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Steriliseret med etylenoksid / Sterilizováno ethylénovým oxidem / Esterilizado por óxido de etileno
STERILE R	Sterilized using irradiation / Stérilisé par irradiation / Durch Bestrahlung sterilisiert / Esterilizado mediante irradiación / Sterilizzato mediante irraggiamento / Gesteriliseerd door bestraling / Steriliseret ved bestråling / Sterilizováno ozářením / Esterilizado por irradiação
	Manufacturer / Fabricant / Hersteller / Fabricante / Produttore / Fabrikant / Produzent / Výrobce / Fabricante
	Date of manufacture / Date de fabrication / Herstellungsdatum / Fecha de fabricación / Data di fabbricazione / Datum van vervaardiging / Produktionsdato / Datum výroby / Data de fabrico
	Do not use if package is damaged / Ne pas utiliser si l'emballage est endommagé / Nicht verwenden, wenn die Verpackung beschädigt ist / No usar si el envase está dañado / Non utilizzare in caso di imballaggio danneggiato / Niet gebruiken als de verpakking beschadigd is / Ikke bruk hvis pakningen er skadet / Nepoužívejte, pokud je obal poškozen / Não utilizar se a embalagem estiver danificada
	Do not re-use / Ne pas réutiliser / Nicht wiederverwenden / No reutilizar / Non riutilizzare / Niet opnieuw gebruiken / Engangsbruk / Nepoužívejte opakovaně / Não reutilizar
EC REP	Authorized representative in the European Community / Représentant autorisé établi au sein de la Communauté européenne / Bevollmächtigter Vertreter in der Europäischen Gemeinschaft / Representante autorizado en la Unión Europea / Rappresentante autorizzato nella Comunità Europea / Oficiële vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autorizovaný zástupce v Evropském společenství / Mandatário na Comunidade Europeia
CH REP	Authorized representative in Switzerland / Représentant autorisé en Suisse / Bevollmächtigte Vertretung in der Schweiz / Representante autorizado en Suiza / Rappresentante autorizzato in Svizzera / Oficiële vertegenwoordiger in Zwitserland / Autoriseret representant i Sveits / Autorizovaný zástupce ve Švýcarsku / Mandatário na Suíça

	Single sterile barrier system / Système de barrière stérile unique / Einzelnes steriles Barriersystem / Sistema de barrera estéril individual / Sistema com barrera sterile singola / Enkelvoudig steriel barrièresysteem / System med én steril barriere / Jednoduchý sterilní barierový systém / Sistema de barreira estéril única
UDI	Unique Device Identifier / Identifiant unique du dispositif (UDI) / Eindeutige Gerätekennung / Identificador único de dispositivo / ID univoco del dispositivo / Uniek apparaat-id / UDI (Unique Device Identifier) / Jediný identifikační prostředek / Identificador Único do Dispositivo
	Do not resterilize / Ne pas restériliser / Nicht erneut sterilisieren / No reesterilizar / Non ri-sterilizzare / Niet opnieuw steriliseren / Må ikke resteriliseres / Nesterilizujzte opakovaně / Não reesterilizar
	Caution / Attention / Vorsicht / Atención / Attenzione / Opgelet / Forsiktig / Upozornění / Atenção
	MR Unsafe / Non compatible IRM / MR-Unsicher / No es seguro para RM / Non sicuro per RM / MR onveilig / Ikke MR-sikker / Nebezpečné pro MR / Não seguro para RM
	Caution: US Federal law restricts this device to sale by or on the order of a physician / Attention : La loi fédérale américaine n'autorise la vente de ce dispositif qu'à un médecin ou sur ordonnance d'un médecin / Vorsicht: Das US-Bundesgesetz beschränkt den Verkauf dieses Produkts auf Ärzte oder auf Anordnung eines Arztes. / Atención: De acuerdo con la legislación federal de los EE. UU., este dispositivo solo puede ser vendido por médicos o por prescripción facultativa / Attenzione: La legge federale degli Stati Uniti consente la vendita di questo dispositivo esclusivamente su prescrizione medica / Opgelet: de federale wetgeving van de Verenigde Staten beperkt de verkoop van dit hulpmiddel tot de verkoop door of op voorschift van een arts / Forsiktig: Federal (amerikansk) lovgivning begrenser salg av denne enheten til kun å selges av, eller på resept fra lege / Upozornění: Federální zákony USA omezují prodej tohoto prostředku na lékaře nebo na jeho objednávku / Atenção: a lei federal dos EUA restringe a venda deste dispositivo a um médico ou mediante prescrição médica

Not made with natural rubber latex. / Ne contient pas de latex de caoutchouc naturel. / Ohne Naturkautschuklatex hergestellt / No está fabricado con látex de caucho natural. / Non fabbricato con lattice di gomma naturale. / Niet vervaardigd met natuurlijk rubber latex. / Ikke laget av naturlig gummilatteks / Pfi výrobě nebyl použit přírodní kaučukový latex. / Não fabricado com látex de borracha natural.

Summit Medical LLC
815 Vikings Parkway, Suite 100
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939
F: 1-888-229-1941 | +1 651-789-3979
www.summitmedicalusa.com

 0459

EC REP

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
CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Epistax™ Balloon Catheter

 Only

Instructions for Use Mode d'emploi Gebrauchsanweisung Modo de empleo Istruzioni per l'uso Gebruiksaanwijzing Bruksanvisning Návod k použití Instruções de utilização

DESCRIPTION

Epistax™ is a two-balloon silicone catheter with two independently inflatable silicone balloons designed to manage internal nasal bleeding. The smaller balloon marked "10cc" is designed to manage bleeding in the posterior chamber, while the larger balloon marked "30cc", manages bleeding in the anterior chamber. The catheter allows the patient to breathe through their nose with Epistax™ in place.

INTENDED USE / INDICATION FOR USE
Epistax™ is an inflatable nasal catheter intended to be used to manage internal nasal bleeding.

INTENDED USER & PATIENT TARGET GROUP
Epistax™ is intended to be used by suitably qualified healthcare professionals in a clinic or hospital environment. The patient target group is any adult in the human population as determined by the healthcare professional.

CONTRAINDICATIONS
Epistax™ is made of silicone which may cause an allergic reaction in patients with sensitivity to this material. Surgeons should risk assess the use of Epistax™ in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to silicone.

WARNINGS & PRECAUTIONS

- Epistax™ is for SINGLE USE ONLY. Do NOT re-sterilize or re-use. Re-use of device could result in infection, contamination, or device failure.
- Do NOT use and discard if the package has been opened or damaged or the date of expiration has passed. If the sterile barrier is opened or damaged, please contact the manufacturer with the product code and lot number listed on the Epistax™ device.
- Do NOT use petroleum-based lubricants with Epistax™ as it will affect the performance of the balloons.
- Sterile saline is not to be used to inflate Epistax™ as the salts used may crystallize and cause a blockage.
- Care must be taken not to over-inflate the balloon.
- Epistax™ has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Epistax™ in the MR environment is unknown and scanning a patient who has Epistax™ may result in patient injury.

INSTRUCTIONS FOR USE
Test Inflation

- Attach a 6% tapered Luer connector syringe to the valve marked "10cc", and test-inflate the distal balloon with 5 to 10ml of air. Deflate balloon. Repeat for valve marked "30cc", test-inflate with 15 to 25 ml of air. Deflate balloon. If either balloon does not inflate, use a new Epistax™ and call Summit Medical Customer Service.

Insertion

- Lubricate Epistax™ using sterile water or water-based lubricant.
- After removal of nasal blood clots, insert the catheter so that the beveled tip is in the nasopharynx.
- Fill a 6% tapered Luer connector syringe with sterile water.
- Insert tip of the syringe into valve marked "10cc" to inflate posterior balloon. Inflate balloon with fluid from the syringe. Inflation volumes of 4cc to 8cc are usually adequate.

CAUTION: Do not exceed the maximum fill volume of the balloon.

- Remove the syringe and gently pull the catheter anteriorly to seat the posterior balloon.
- Insert the tip of the syringe into valve marked "30cc" to inflate anterior balloon. Inflate balloon with fluid from syringe. Inflation volumes of 10cc to 25cc are usually adequate.

CAUTION: Do not exceed the maximum fill volume of the balloon.

- Check fit and position of the balloons as required. Adjust pressure accordingly using the syringe.

Removal

- Epistax™ may be removed as indicated by a suitably qualified healthcare professional. Maximum use is 5 days.
- To remove Epistax™, insert the syringe and evacuate each balloon completely. After balloons are completely deflated, slowly remove catheter, and dispose as per IFU disposal instructions.

Disposal
After use, Epistax™ should be disposed of in accordance with the healthcare settings disposal procedure. Epistax™ comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the Epistax™ after use to prevent contamination. If Epistax™ includes a syringe, the syringe is considered a sharp and must be placed in a secure puncture-resistant bin (conforming to BS EN ISO 23907) which is suitable for incineration.

PERFORMANCE CHARACTERISTICS

- Balloons will inflate to pressures required to manage internal nasal bleeding
- Balloons will maintain pressure required to manage internal nasal bleeding while Epistax™ is in situ

POSSIBLE ADVERSE EFFECTS

- Headache and pain
- Structural nostril damage

CLINICAL BENEFITS

- The balloons can be inflated to different sizes to allow for different anatomies
- Balloons will maintain the required pressure to manage internal nasal bleeding
- The pressure of the balloons is easy to adjust using the syringe
- Can be used for up to 5 days to manage internal nasal bleeding

STERILITY
Epistax™ is provided terminally sterilized by ethylene oxide (EO) gas. The product is sterile in an undamaged, unopened package. For product codes RH-4220-01 & RH-4420-01, the syringe is provided gamma irradiation sterilized.

INCIDENT REPORTING
Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the country in which the user and/or patient is established.

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