# EarPopper®

REF

## EP-2100A/EP-2100

- (EN) MIDDLE EAR INFLATION DEVICE INSTRUCTIONS FOR USE
- (FR) APPAREIL DE GONFLAGE DE L'OREILLE MOYENNE MODE D'EMPLOI
- (DE) MITTELOHR-INFLATIONSGERÄT GEBRAUCHSANWEISUNG
- (ES) DISPOSITIVO DE HINCHADO PARA EL OÍDO MEDIO MODO DE EMPLEO
- (IT) DISPOSITIVO DI GONFIAGGIO PER L'ORECCHIO MEDIO ISTRUZIONI PER L'USO
- (NL) MIDDENOOR BLAASHULPMIDDEL GEBRUIKSAANWIJZING
- (NO) INFLASJONSENHET FOR MELLOMØRE INSTRUKSJONER FOR BRUK AV
- (KO) 중이 확장 장치 사용 지침
- (HE) מכשיר לניפוח האוזן התיכונה הוראות שימוש
- (AR) تعليمات الاستخدام جهاز علاج انتفاخ الأذن الوسطى
- (JA) 中耳加圧装置 使用説明書
- (ZH) 中耳充氣裝置 使用説明書
- (ZH-CN) 中耳充气装置 使用说明
- (LT) VIDURINĖS AUSIES SLĖGIO REGULIAVIMO PRIETAISAS NAUDOJIMO INSTRUKCIJA
- (PL) URZĄDZENIE DO LECZENIA UJEMNEGO CIŚNIENIA W UCHU ŚRODKOWYM INSTRUKCJA OBSŁUGI
- (FI) VÄLIKORVAN PAINEENTASAUSLAITE KÄYTTÖOHJEET

**TABLE OF CONTENTS**

<b>Language</b>	<b>PAGE</b>
<b>EN</b> English	3
<b>FR</b> Français	8
<b>DE</b> Deutsch	13
<b>ES</b> Español	18
<b>IT</b> Italiano	23
<b>NL</b> Nederlands	28
<b>NO</b> Norsk	33
<b>KO</b> 한국어	38
<b>HE</b> עברית	43
<b>AR</b> الإنجليزية	48
<b>JA</b> 日本語	53
<b>ZH</b> 繁體中文	58
<b>ZH-CN</b> 简体中文	63
<b>LT</b> Lietuvių k.	68
<b>PL</b> Polski	73
<b>FI</b> Suomi	78

## EP-2100A/EP-2100 | MIDDLE EAR INFLATION DEVICE

### INTENDED USE/INDICATIONS FOR USE

The EarPopper is intended to be used for the treatment of Otitis Media with Effusion, Eustachian Tube Dysfunction, and ear fullness caused by air travel and scuba diving.

The EarPopper is indicated for use in the treatment of negative middle ear pressure.

Negative middle ear pressure can lead to fluid accumulation in the middle ear, impaired hearing and hearing loss. The EarPopper provides a method for ventilating the middle ear by introducing air through the nose and Eustachian tube into the middle ear. Equalizing middle ear pressure can prevent the accumulation of fluid and prevent hearing loss.

### INTENDED USER


The EarPopper is intended to be used by the patient or parent/guardian after consulting a physician.

### PRECAUTIONS

#### • CAUTION: A LICENSED PHYSICIAN SHOULD BE CONSULTED PRIOR TO USE.

- To ensure compliance, the physician should demonstrate proper use of the EarPopper.
- **DO NOT** force the instrument into the nasal cavity.
- **DO NOT** use without a clean, intact nose piece attached.
- Adult supervision is required when using this product on children.
- Submerging device in water may result in product failure.
- Stop using EarPopper if ear discomfort occurs during use.
- Patients with pre-existing conditions should seek additional advice from a healthcare professional prior to use. Use of the device in these circumstances may lead to perforation of the eardrum.
- Discontinue use if device is damaged.

### PARTS OF THE EARPOPPER

## CONTRAINDICATIONS

**DO NOT** use the EarPopper in the presence of ear infection, upper respiratory infection, nasal congestion, cold symptoms or a perforated eardrum.

The device is not intended for multiple users.

## OPERATING INSTRUCTIONS

1. Pull and remove battery strip before initial use.
2. Hold the EarPopper nosepiece firmly against one nostril. Press the other nostril closed using the index finger of the other hand. (**See Fig. 1**) It is very important to establish a good seal. Press the power button to begin airflow and swallow with mouth closed while the device is running.

### Note:

- Holding a small amount of liquid in mouth prior to the procedure will aid in swallowing.
- The treatment only occurs during the swallow.

Fig. 1

If the user does not feel any change in the ear, repeat the procedure making sure the EarPopper remains firmly against the nostril, while the user swallows with mouth closed.

3. Release power button after swallowing.
4. Repeat steps 2-3 for the other nostril.

Wait five minutes and repeat steps 2-4. This will complete one treatment.

## TREATMENT RECOMMENDATIONS

Treatment recommendations are guidelines based on clinical studies. Always consult your physician before beginning any treatment.


- For fluid behind the eardrum (Otitis Media with Effusion), use twice daily for 7-11 weeks, or until hearing returns to normal.
- For negative middle ear pressure related to air travel or diving, use as needed to equalize pressure. If condition does not improve, consult your physician.

## CARE, CLEANING, & STORAGE

Remove the nosepiece (EP-1010) for periodic cleaning. (See Fig. 2) Clean with mild soap and water. After cleaning, rinse with warm water and allow to thoroughly dry before using. To reattach the nosepiece align the ribs with the slots in the tip and push firmly. Removable nosepiece is single patient use only.

Fig. 2

Removable Nosepiece (EP-1010)

The outside of the device may be cleaned by wiping with a soft, damp cloth.

**DO NOT AUTOCLAVE, USE ABRASIVE CLEANERS, CHEMICALS, OR DISINFECTANTS, AND DO NOT SUBMERGE THE EARPOPPER IN LIQUID.**

Operating environment is 5° C to 40° C (41° F to 104° F), 15 - 93% RH, 700 hPa - 1060 hPa.

Store the EarPopper in a clean dry location at room temperature with the clean nosepiece in place to protect tip from damage and dirt. Remove battery for prolonged storage (6+ months).

## BATTERY REPLACEMENT

Use the following steps to check for a dead battery:

1. Place finger over hole in the end of the nosepiece.
2. Press the Power Button to activate device.
3. If motor runs intermittently or stops, replace the battery.

Open the battery compartment by depressing the thumb tab on the back of the device, then rock the battery cover forward and remove. Install 4 new AAA batteries. Replace the cover by placing the front tab in the device and then rock the thumb tab into place. (See Fig. 3)

Battery life: 4+ hours of continuous operation.

Fig. 3

## **DEVICE DISPOSAL FOR REST OF WORLD**

The nosepiece should be disposed of as normal household waste. Dispose of device and battery according to national or local regulations. Contains no lead, cadmium, mercury or bromide flame retardants.

## **DEVICE DISPOSAL IN THE EU**

As the EarPopper is an electrical device it cannot be treated as normal household waste. It should be handed to the applicable collection point for the recycling. For more information about recycling this product or disposing of it through a take-back program contact [customerservices@innoviamedical.com](mailto:customerservices@innoviamedical.com).

## **WARNINGS**

- No modification of this equipment is allowed.
- Risk of choking hazard - product contains small parts (nosepiece and batteries). Keep away from children and pets.

## **WARRANTY**

This device is warranted to be free from manufacturing defects for a period of one year from the date of purchase under normal use. This excludes damage to the device resulting from accident or misuse. 5 year service life.

## **CUSTOMER SERVICE**

For questions concerning the operation of the EarPopper, please visit [www.earpopper.com](http://www.earpopper.com) or call or write using the contact information on the back of this manual.

**Note:** Product reviews through the Amazon website cannot be investigated as contact details are not provided to Summit Medical.

EarPopper complaints should be reported directly to [customerservices@innoviamedical.com](mailto:customerservices@innoviamedical.com).

## **SERIOUS INCIDENT REPORTING**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the state in which the user and/or patient is established.

## **TECHNICAL INFORMATION**

**PRODUCT SAFETY COMPLIANCE TESTING:** (IEC 60601-1)  
IEC 60601-1-1

**CLASSIFICATIONS:** IEC 60601-1: protection against electrical shock: internally powered equipment, type BF applied parts protection against harmful ingress of water or particulate matter: IP22 Mode of operations: continuous.  
IEC 60601-1-2

EN61000-4-2 Electrostatic Discharge	PASS	15 kV air 8 kV contact
EN61000-4-3 Radiated Immunity	PASS	10 V/m
EN61000-4-4 Electrical Fast Transients	N/A	————
EN61000-4-5 Surge Immunity	N/A	————
EN61000-4-6 Conducted Immunity	N/A	————
EN61000-4-8 Power Frequency Magnetic Field	PASS	30 amp/m
EN61000-4-11 Voltage Dips, Interrupts, and Fluctuations	N/A	————

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Summit Medical LLC
MDSS GmbH