
- (EN) INSTRUMENT STERILIZATION/PROTECTION TRAYS
- (FR) PLATEAUX DE PROTECTION/STÉRILISATION DES INSTRUMENTS
- (DE) STERILISATIONS-/SCHUTZSCHALEN FÜR INSTRUMENTE
- (ES) BANDEJAS DE ESTERILIZACIÓN/PROTECCIÓN DE INSTRUMENTOS
- (IT) VASSOI PER LA PROTEZIONE/STERILIZZAZIONE DEGLI STRUMENTI
- (PT-BR) BANDEJAS DE ESTERILIZAÇÃO/PROTEÇÃO DE INSTRUMENTOS
- (JA) 医療器具滅菌/プロテクショントレイ
- (DA) BAKKER TIL STERILISATION/BESKYTTELSE AF INSTRUMENTER
- (FI) INSTRUMENTIN STERILOINTI-/SUOJATELINEET
- (CS) TÁČY PRO STERILIZACI/OCHRANU NÁSTROJŮ
- (BG) ТАБЛИ ЗА СТЕРИЛИЗАЦИЯ/ЗАЩИТА НА ИНСТРУМЕНТИ
- (SL) PLADNJI ZA STERILIZACIJO/ZAŠČITO INSTRUMENTOV
- (SK) PODNOSY NA STERILIZÁCIU/OCHRANU NÁSTROJOV
- (ET) INSTRUMENDI STERILISEERIMIS-/KAITSEALUSED
- (RO) TĂVI DE STERILIZARE/PROTECȚIE A INSTRUMENTELOR
- (PL) TACE DO STERYLIZACJI / OCHRONY NARZĘDZI
- (NL) STERILISATIE/BESCHERMINGSSCHALEN VOOR INSTRUMENTEN
- (LV) INSTRUMENTU STERILIZĀCIJAS/AIZSARDZĪBAS PĀPLĀTES
- (NO) INSTRUMENTBRETT FOR STERILISERING/BESKYTTELSE
- (SV) INSTRUMENTSTERILISERING/SKYDDSTRÅG
- (HR) PLITICE ZA STERILIZACIJU/ZAŠTITU
- (EL) ΔΙΑΚΟΙ ΑΠΟΣΤΕΙΡΩΣΗΣ/ΠΡΟΣΤΑΣΙΑΣ ΟΡΓΑΝΩΝ
- (HU) MŰSZER-STERILIZÁLÓ/-VÉDŐ TÁLCÁK
- (LT) INSTRUMENTŲ STERILIZAVIMO / APSAUGINIAI PADĖKLAI
- (TR) ALET STERİLİZASYON/KORUMA TEPŞİLERİ

INSTRUCTIONS FOR USE / MODE D'EMPLOI / GEBRAUCHSANWEISUNG / INSTRUCCIONES DE USO / ISTRUZIONI PER L'USO / INSTRUÇÕES DE USO / 使用説明書 / BRUGSANVISNING / KÄYTTÖOHJEET / NÁVOD K POUŽITÍ / ИНСТРУКЦИИ ЗА УПОТРЕБА / NAVODILA ZA UPORABO / POKYNY NA POUŽITIE / KASUTUSJUHEND / INSTRUCȚIUNI DE UTILIZARE / INSTRUKCJA UŻYTKOWANIA / GEBRUIKSAANWIJZING / LIETOŠANAS INSTRUKCIJA / BRUKSANVISNING / ANVÄNDARINSTRUKTIONER / UPUTE ZA UPORABU / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / HASZNÁLATI ÚTMUTATÓ / NAUDOJIMO INSTRUKCIJA / KULLANMA TALIMATLARI

LANGUAGE	PAGE
English	1
Français	8
Deutsch	15
Español	22
Italiano	29
Português	36
日本語	43
Dansk	50
Suomi	57
Čeština	64
Български	71
Slovenski	78

LANGUAGE	PAGE
Slovenčina	85
Eesti	92
Română	99
Polski	106
Nederlands	113
Latviski	120
Norsk	127
Svenska	134
Hrvatski	141
Ελληνικά	148
Magyar	155
Lietuvių k.	162
Türkçe	169

TABLE OF CONTENTS | ENGLISH

SUBJECT	PAGE
Intended Use/Indications for Use Summary	1
Device Description	1
Intended User	1
Sterilization	1-2
Limitation on Processing	2
Warnings	2
Universal Precautions	2
Point of Use	2
Cleaning	2
Disinfection	2
Maintenance, Inspection, and Testing	3
Indications for Use	3-5
Storage	6
Disposal	6
Serious Incident Reporting	6
Warranty	6
Appendix A	7

The instructions provided within have been validated by the tray manufacturer as being capable of reprocessing reusable medical devices.

Individual sterilizers, instrument cleanliness, specific loading of instrument trays, types and geometry of instruments, sterilization containers, filters and wrappings vary at each location.

These instructions only describe possibilities and validations of the instrument sterilization/protection trays and not the BK Medical devices to be protected in these trays. For reprocessing instructions and recommendations for the BK Medical devices please refer to BK Medical "Care and Cleaning" at: www.bkmedical.com/care-and-cleaning/

READ THIS SECTION BEFORE PLACING PRODUCT INTO SERVICE

INTENDED USE/INDICATIONS FOR USE SUMMARY

Instrument Sterilization/Protection Trays are intended to contain and protect reusable medical devices during transport, sterilization, and storage. Instrument Sterilization/Protection Trays are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instrument Sterilization/Protection Trays are intended to allow sterilization of the enclosed medical devices during these sterilization cycles:

- Pre-Vacuum Steam (EN 285)
- Ethylene Oxide (ISO 11135)
- Validated Gas Plasma Sterilization i.e.
 - STERRAD® 100S Standard
 - STERRAD® 100NX® Standard
 - STERRAD® 100NX® Express
 - STERRAD® 100NX® Flex
 - STERRAD® NX® Standard
 - STERIS® AMSCO® V-PRO 1®, V-PRO 1 Plus, V-PRO maX and V-PRO maX 2

The Instrument Sterilization/Protection Trays are intended to be used in conjunction with a legally marketed wrap, Aesculap® rigid containers, or Genesis™ rigid containers. The System trays are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

DEVICE DESCRIPTION

Instrument Sterilization/Protection Trays are used to enclose and hold surgical instruments and instrument accessories in an organized manner during the sterilization process and subsequent storage and transportation. The trays do not have direct patient contact. The trays by themselves do not maintain sterility. The trays are different sizes of the same basic configuration: a rectangular base with latchable cover. The trays have perforations to allow sterilant penetration. The trays contain silicone inserts in the base and/or cover to hold, organize, and protect the surgical instruments within the tray during the sterilization process and subsequent storage and transportation.

INTENDED USER

Instrument Sterilization/Protection Trays are intended to be used by healthcare professionals in the operating room and sterile processing department for transportation, sterilization, and storage of medical devices.

STERILIZATION

See **Table 1** for typical sterilization parameters that have been qualified for the sterilization of Instrument Sterilization/Protection Trays. These are only valid with sterilization equipment that is properly maintained and calibrated. Minimum dry times are required to prevent wet packs.

Table 1

STERILIZATION METHOD	CYCLE (times)
Pre-Vacuum Steam	Parameter:
	Temperature 270°F (132°C)
	Expose Time 4 minutes Dry Time 30 minutes
Ethylene Oxide (EO)	Preconditioning:
	Temperature 131°F (55°C)
	Relative humidity 70±15%
	Precondition time 1 hour
	Sterilization:
Exposure time 120 minutes	
Temperature 131°F (55°C)	
Aeration time 12 hours	
STERRAD 100S	Standard
STERRAD 100NX	Standard
STERRAD 100NX	Express
STERRAD 100NX	Flex
STERRAD NX	Standard
STERIS AMSCO V-PRO 1	Standard
STERIS AMSCO V-PRO 1 PLUS	Lumen, Non-Lumen
STERIS AMSCO V-PRO maX	Lumen, Non-Lumen
STERIS AMSCO V-PRO maX 2	Lumen, Non-Lumen

The following sterilization methods are validated:

- The 4 minute autoclave sterilization cycle in legally marketed wrap, Aesculap rigid container or Genesis rigid container cleared by the FDA.
- The Ethylene Oxide (EO) sterilization cycle in legally marketed wrap and a Genesis sterile container cleared by the FDA.
- The STERRAD 100S Standard sterilization cycle in legally marketed wrap and an Aesculap rigid container cleared by the FDA.
- The STERRAD 100NX Standard sterilization cycle in legally marketed wrap and an Aesculap rigid container cleared by the FDA.

STERILIZATION (continued)

- The STERRAD 100NX Express sterilization cycle in legally marketed wrap cleared by the FDA.
- The STERRAD 100NX Flex sterilization cycle in legally marketed wrap and an Aesculap rigid container cleared by the FDA.
- The STERRAD NX Standard sterilization cycle in legally marketed wrap and an Aesculap rigid container cleared by the FDA.
- The STERIS AMSCO V-PRO, V-PRO 1 PLUS, V-PRO maX and V-PRO maX 2 sterilization cycles in legally marketed wrap and an Aesculap rigid container cleared by the FDA.

DO NOT exceed the load capacity of the sterile container as specified by the manufacturer.

Use an FDA cleared accessory to maintain sterility.

Please consult the sterilizer instruction manual to ensure intended loads are compatible with the intended sterilization cycle.

Please consult the container instructions to ensure that the intended load is compatible with the FDA cleared loads for the container.

Please consult BK Medical's Care and Cleaning Guide to insure intended devices/loads are compatible with the intended sterilization cycle

LIMITATION ON PROCESSING

- 1) The end of useful life on the Instrument Sterilization/Protection Trays is a minimum of 25 sterilization cycles. Inspect the tray before use for wear and damage caused by use. Discontinue use if visible signs of wear are present, including corrosion, mechanical failures, cracking, peeling, flaking, broken welds, damaged feet, damaged latches, damaged Hold-Its®/ Hold-Downs®, discoloration, etc.
- 2) See **Table 1** for sterilization methods and configurations.
- 3) **DO NOT OVERLOAD** trays or components.
- 4) **DO NOT OVERLOAD** individual Hold-Its® slots. Load only one instrument per Hold-Its® slot.
- 5) For rigid container users, **DO NOT WRAP** trays or components and place inside of container for sterilization.
- 6) Inside of sterilizers, **DO NOT STACK** individually wrapped or containerized trays or components. Separate wrapped or containerized trays or components from each other or any other items on separate shelves of the sterilizer to allow for maximum sterilant flow.
- 7) The use of non absorbent tray liners (e.g. silicone fingered organizing mat) can cause condensate to pool. If visible moisture is present, re-sterilize with a longer dry time.
- 8) The total weight of the container system (e.g. container, tray and instrument load) must not exceed 25 pounds (11.34 kg).
- 9) Instruments (e.g. endoscopes and instruments with lumens or channels) should be prepared and sterilized according to the instrument manufacturer's instructions.
- 10) It is the responsibility of the processor to maintain specific validations for the terminal sterilization process being applied to the configurations of instruments and containers being presented to the sterilization process.

REPROCESSING OF THE TRAYS

WARNINGS

For trays of aluminum construction it is recommended to use only neutral pH (6.0 – 8.5) detergents. A detergent with a highly acidic or highly alkaline pH could permanently damage the anodized aluminum finish of the tray. This warning does not apply to trays of stainless steel construction.

Note: Clean and inspect trays according to the instructions provided prior to placing into service.

DO NOT use if package is damaged or unintentionally opened prior to use.

UNIVERSAL PRECAUTIONS

- Personnel should wear all personal protective clothing and equipment as required by their employer's/department's operating procedures for the contamination level they will be exposed to
- Keep dissimilar metals separated during sterilization to prevent corrosion

POINT OF USE

Remove gross soil with disposable cloth/paper wipe.

Contaminated components should be kept moist until qualified cleaning processes can be applied.

CLEANING

Listed below are recommended cleaning options to clean the cassette/tray. Refer to the instrument manufacturer's instructions for use for cleaning the instrument. Use only neutral pH (6.0 – 8.5) solutions, mildly alkaline and free of sodium carbonate to avoid damaging finish for aluminum cassettes.

DO NOT use scouring pads or abrasive cleaners. **DO NOT** store cassette in liquid.

1. Manual Gross Decontamination:

- A. Equipment: Neutral pH (6.0 – 8.5) enzymatic detergent, soft bristle brush, and running water.
- B. Remove all visible soil and contaminants using a soft bristle brush. The entire tray should be immersed while cleaning to aid in the removal of contaminants and to reduce splashing of detergent on personnel for a minimum of 2 minutes.
- C. Rinse thoroughly for a minimum of 1 minute with clean water to remove all detergent. See rinsing instructions on the detergent label.

2. Ultrasonic Clean:

- A. Prepare enzyme wash in an ultrasonic cleaning unit.
- B. Place a single, tray in the detergent and run for a minimum of ten minutes.
- C. Rinse for a minimum of 2 minutes with cold tap water.
- D. Visually inspect tray for contaminants. Repeat the cycle if necessary to remove visible contamination.

3. Automated Washer:

The Instrument Sterilization/Protection Trays have been validated for the automatic wash system cycle listed in **Table 2**. Qualification of specific parameters will need to be conducted by the processor.

Table 2

CYCLE	WATER TEMP.	CLEANING PROCESS
Pre-wash 1 & 2	Cold Tap Water	Re-circulation Time: 2 minutes
Enzyme Wash	Hot Tap Water	Soaking Time: 4 minutes
Wash 1	Heated at 150°F (65.5°C)	Re-Circulation Time: 15 minutes
Rinse 1 & 2	Hot Tap Water	Re-Circulation Time: 5 minutes

Note: After completion of an automatic wash cycle, visually inspect the tray and accessories for any remaining visible soil. All visible soil must be removed by hand cleaning, brushing, ultrasonic, or additional automatic cycles prior to sending to sterilization.

DISINFECTION

Trays are intended to be terminally sterilized.

MAINTENANCE, INSPECTION, AND TESTING

Trays may be reused until unacceptable deterioration such as corrosion, cracking, rust, peeling, flaking, discoloration, or mechanical failure occurs.

Signs of Mechanical Failure Include:

- Broken or cracked corners or welds
- Broken or non-working latches
- Missing, torn, or cut silicone inserts
- Missing or damaged feet

INDICATIONS FOR USE

4 Minute Steam, Wrap & Aesculap® Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a pre-vacuum steam sterilization cycle. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap or Aesculap rigid container. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility.

Autoclave Sterilization Parameter	
Cycle	Pre-Vacuum
Temperature	270°F (132°C)
Exposure Time	4 minutes
Minimum Dry Time	30 minutes
Summit Cassette Model	Aesculap Container Model
IN-8823-AE	*JN444
IN-2880	*JK444
IN-6105	*JN742
*Validated by Summit Medical for use in Pre-Vacuum Steam sterilizers ONLY operating at 270°F (132°C) for 4 minutes exposure time. Consult container instructions to ensure that contents do not exceed the sterilization container's intended load claims.	

Lumen Claims for 4 Minute Pre-Vacuum Steam Sterilization Cycle			
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens
IN-8823-CF	3mm	400mm	4
IN-8823-CF	3mm	200mm	2
IN-2880	1mm	76mm	2
IN-2880	3mm	177mm	1
IN-6105	5mm	241mm	1
IN-2681	1mm	65mm	1
IN-2681	3mm	200mm	1
IN-7823	1mm	400mm	17

4 Minute Steam, Genesis™ Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a pre-vacuum steam sterilization cycle. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed Genesis rigid containers. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

Lumen Claims for 4 Minute Pre-Vacuum Steam Sterilization Cycle			
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens
IN-2681	1mm	65mm	1
IN-2681	3mm	200mm	1
IN-6105	5mm	241mm	1
IN-0000	1mm	400mm	5
	3mm	400mm	1
	5mm	400mm	1

8 Minute Steam & Ethylene Oxide (EO), Wrap & Genesis™ Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe System cassettes are intended to allow sterilization of the enclosed medical devices during Pre-Vacuum steam or ethylene oxide sterilization cycles. The InstruSafe System cassettes are intended to be used in conjunction with central legally marketed wrap or with a Genesis rigid container. The InstruSafe System cassettes are not intended on the own to maintain sterility.

Sterilization methods and configurations

Steam 8 Minute Preconditioning at 132°C	
140 count woven wrap	10 minute dry time
Non-woven wrap (Kimberly Clark)	50 minute dry time
Genesis Container (reference Table 1 for filter paper to use)	30 minute dry time

Ethylene Oxide (EO)	
<ul style="list-style-type: none"> ▪ 1 hour preconditioning at 131°F (55°C) with relative Humidity of 70 ± 15% ▪ 2 hours exposure at 131°F (55°C) ▪ 12 hours aeration 	

Table 1

Genesis Container	Genesis Container Filter Paper
CD2-4B	FX3-1: 9x9"
CD2-5B	FX3-1: 9x9"
CD3-4B	FX3-1: 9x9"
CD3-5B	FX3-1: 9x9"
CD3-6B	FX3-1: 9x9"
CD3-7B	FX3-1: 9x9"
CD4-5B	FO3-2: 9x6"
CD5-61B	FO3-2: 9x6"
CD6-6B	FX3-1: 9x9"

See V. Mueller Genesis Sterilizer Container system User's Operating Instructions when using any Genesis Rigid Container.

INDICATIONS FOR USE (continued)

STERRAD® 100S & STERRAD® 100NX® Standard Cycles, Wrap & Aesculap® Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a Sterrad 100S Standard and Sterrad 100NX Standard sterilization cycles. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with legally marketed wrap or Aesculap rigid container. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

The following sterilization trays were validated with the corresponding rigid containers:

STERRAD 100S Standard Cycle	
Summit Cassette Model	Aesculap Container Model
IN-8823-AE	*JM444
IN-6105	*JM440

STERRAD 100NX Standard Cycle	
Summit Cassette Model	Aesculap Container Model
IN-8823-AE	*JM444
IN-6105	*JM440

**Validated by Summit Medical for use in STERRAD 100S Standard Cycle and STERRAD 100NX Standard Cycle ONLY. Consult container instructions to ensure that contents do not exceed the sterilization container's intended load claims.*

Lumen Claims for STERRAD 100S Standard and STERRAD 100NX Standard Cycles	
4 Stainless steel lumens with 3mm inner diameter x 400mm length	All appropriately sized models are listed in Appendix A with the exception of IN-2681.
2 Stainless steel lumens with 3mm inner diameter x 200mm length	
1 Stainless steel lumens with 3mm inner diameter x 200mm length	IN-2681
1 Stainless steel lumens with 1mm inner diameter x 65mm length	

STERRAD® NX® Standard Cycle, Wrap & Aesculap® Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a Sterrad NX Standard Sterilization Cycle. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap or Aesculap rigid container. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

Sterilization Methods and configurations
- Sterrad NX Standard Sterilization Cycle

Lumen Claims for STERRAD NX Standard Sterilization Cycle			
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens
IN-2681	1mm	65mm	1
IN-2681	3mm	200mm	1
IN-8987-S	1mm	500mm	5
IN-8615	2.3mm	210mm	5
IN-6105	4mm	235mm	1

Note: The worst case validated load based on vent-to-volume calculation is the IN-2681 tray.

STERRAD® 100NX® Express Cycle, Wrap

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a Sterrad® 100NX Express Cycle. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. The InstruSafe Instrument Protection System has no lumen claims for the Sterrad 100NX Express Cycle.

INDICATIONS FOR USE (continued)

STERRAD® 100NX® Flex Sterilization Cycle, Wrap & Aesculap® Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. InstruSafe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during a Sterrad 100NX Flex Sterilization Cycle. The InstruSafe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap or Aesculap rigid container. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

STERRAD 100NX Flex Sterilization Cycle	
Summit Cassette Model	Aesculap Container Model
IN-0000	*JM440
IN-6105	*JM440
*Validated by Summit Medical for use in STERRAD 100NX Flex Sterilization Cycle ONLY. Consult container instructions to ensure that contents do not exceed the sterilization container's intended load claims.	

Lumen Claims for STERRAD 100NX Flex Sterilization Cycle				
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens	Wrap/Rigid Container
IN-0000	1mm	850mm	1	Wrap + Rigid Container
IN-8823	1mm	850mm	1	Wrap + Rigid Container
IN-7344	1mm	850mm	1	Wrap
IN-6105	4mm	235mm	1	Wrap + Rigid Container

The worst case validated load based on vent-to-volume calculation is the IN-0000 tray.

Note: The IN-0000 tray is for testing purposes only.

AMSCO® V-PRO® Low Temperature Sterilization Cycles, Wrap & Aesculap® Rigid Container

InstruSafe® Instrument Protection System cassettes are used to organize and protect other medical devices that are sterilized by a healthcare provider. Instru-Safe Instrument Protection System cassettes are intended to allow sterilization of the enclosed medical devices during Amsco V- PRO Low Temperature Sterilization Cycles. The Instru-Safe Instrument Protection System cassettes are intended to be used in conjunction with a legally marketed wrap or Aesculap rigid container. The InstruSafe Instrument Protection System cassettes are not intended on their own to maintain sterility. A full list of device models is provided in Appendix A.

AMSCO V-PRO Low Temperature Sterilization Systems			
Sterilizer	Standard Cycle	Lumen Cycle	Non Lumen Cycle
V-PRO 1	X	N/A	N/A
V-PRO 1 PLUS	N/A	X	X
V-PRO maX	N/A	X	X
V-PRO maX 2	N/A	X	X
Summit Cassette Model		Aesculap Container Model	
IN-8823		*JM444	
IN-6105		*JM742	
*Validated by Summit Medical for use in AMSCO V-PRO Low Temperature Sterilization Systems ONLY. When using Aesculap container as sterile barrier, the load (Summit tray and contents), should not exceed the load claims for the container's weight or load type.			

Lumen size of instrumentation validated includes:			
Summit Cassette Model	Minimum Inside Diameter	Maximum Length	Number of Lumens
IN-8823	3mm	400mm	2
IN-6105	3mm	200mm	1
IN-2681	1mm	64mm	1
Note: The worst case validated load based on vent-to-volume calculation is the IN-2681 tray.			

STORAGE

Store terminally sterile trays that are wrapped or containerized on storage shelf in a horizontal position. Consult wrap or container manufacturer for shelf life information.

DISPOSAL

In the event the Instrument Sterilization/Protection Trays do not pass inspection prior to use or have otherwise been deemed no longer fit for purpose, the devices shall be disposed of in line with local protocol. The method of disposal shall depend on the potential risks of cross-contamination and infection when the need for disposal is identified.

SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the Member State in which the user and/or patient is established.

WARRANTY

LIMITED WARRANTY FOR INSTRUMENT STERILIZATION/ PROTECTION TRAYS.

THIS LIMITED WARRANTY AND THE REMEDY PROVIDED HEREIN ARE EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS WARRANTIES AND, UNLESS STATED HERE-IN, ANY STATEMENTS OR REPRESENTATIONS MADE BY ANY OTHER PERSON OR FIRM ARE VOID. THE DURATION OF ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL BE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY. NEITHER SUMMIT MEDICAL NOR ITS AFFILIATES SHALL BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL LOSSES OR DAMAGES, RESULTING FROM THE USE OR INABILITY TO USE THE TRAY, WHETHER RESULTING FROM BREACH OF WARRANTY OR ANY OTHER LEGAL THEORY.

This Limited Warranty gives you specific legal rights, and you may also have other rights which vary from State to State.

Some States do not allow limitations on how long an implied warranty lasts, or do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusions may not apply to you.

What Is Covered. Summit Medical warrants the original purchaser that the tray enclosed with this Limited Warranty conforms to the manufacturer's specifications and is free from defects in workmanship and material for a period of 12 months from the date of original purchase. If the original purchaser transfers the tray to another party, this Limited Warranty will not be enforceable by the party to whom the product is transferred.

What We Will Do To Correct Problems. Should your tray prove defective during this period, you must notify Summit Medical or an authorized distributor or dealer of Summit Medical. You must permit Summit Medical or its representatives to make such investigation, examination and tests as Summit Medical deems appropriate and, if requested to do so, you will return the product to the factory at the address set forth below. Summit Medical's sole obligation under this Limited Warranty is, at its option, to repair or replace the defective product or products, without charge for parts or labor. Postage, insurance or shipping costs incurred in presenting your tray product for warranty service are your responsibility.

What Is Not Covered. This Limited Warranty is contingent upon proper use and maintenance of the product; it does not cover products that have been improperly shipped, or that have been misused, abused, neglected, or improperly maintained, cleaned or stored, or that have been serviced other than by Summit Medical or an authorized distributor or dealer of Summit Medical, or that have been modified without the express approval of Summit Medical. Failure to follow the directions in the owner's manual may constitute improper use or maintenance of the product and causes this Limited Warranty not to apply. This Warranty does not extend to normal wear or to replacement items.

If you have questions or claims related to this warranty, contact distributor:

BK Medical
8 Centennial Drive
Peabody, MA 01960 USA
Tel: +1 978 326 1300
bkmedical.com

North America
Sales & Service:
BK Medical
8 Centennial Drive
Peabody, MA 01960 USA
Tel: +1 978 326 1300
Fax: +1 978 326 1399
bkmedical.com

Europe and Rest of World
Sales, Service & Design Center
BK Medical
Mileparken 34
2730 Herlev Denmark
Tel: +45 4452 8100
Fax: +45 4452 8199
bkmedical.com

Asia Sales & Service:
BK Medical Technology
Shanghai Co., Ltd.
Room 602, Building 2, No. 1, Alley 777
Wanrong Road, Jingan District,
Shanghai, China 20072
Tel: +86 21 3603 6378
bkmedical.com

BK Medical Aps, Mileparken 34, 2730 Herlev Denmark
Tel: +45 4452 8100 Fax: +45 4452 8199

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
Asya Satış ve Servis:
BK Medical Technology
Shanghai Co., Ltd.
Room 602, Building 2, No. 1, Alley 777
Wanrong Road, Jingan District,
Shanghai, China 20072
Tel: +86 21 3603 6378
bkmedical.com

BK Medical Aps, Mileparken 34, 2730 Herlev Denmark
Tel: +45 4452 8100 Faks: +45 4452 8199

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CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

InstruSafe is a registered trademark of Summit Medical, LLC. Genesis is a trademark of CareFusion Corporation and Aesculap is a registered trademark of B. Braun Company. All rights reserved.

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