## Rhino Rocket® Nasal Packs (with applicator)

**i<sup>1</sup>** Before using, read the following information:

### DESCRIPTION

The Rhino Rocket® Nasal Packs are made from absorbent polyvinyl alcohol (PVA) foam sponge with open, interconnected cells and a plastic applicator which come in eight (8) sizes as well as a starter pack, which has one (1) of each size. The Rhino Rocket® Nasal Packs have a nylon string that allows for fixation of the device to the cheek or to loop around the ear and removal of the pack from the patient's nasal cavity.

### INTENDED USE

Rhino Rocket® Nasal Packs are intended to be used in the nasal or sinus cavity.

**INDICATIONS FOR USE**
Rhino Rocket® Nasal Packs are indicated to absorb exudate, manage epistaxis and stabilize the nasal cavity.

### TARGET PATIENT GROUP

Anyone of any age in the human population as determined by a suitably trained healthcare professional.

### INTENDED USER

Healthcare professionals

### CLINICAL BENEFITS

- Exudate management
- Epistaxis management
- Stabilisation of the nasal cavity
- Prevention of adhesions
- Prevention of haematoma formation.
- Devices can be left in-situ for up to 24 hours to aid haemostasis

### CONTRAINDICATIONS

No contra-indications associated with the use of the Rhino Rocket® Nasal Packing devices have been identified.

### PRECAUTIONS

Intended for single use only, do not resterilize. Do not use if pouch is damaged. The plastic applicator must not be inserted into the nasal cavity as it may damage the soft tissue.

### POSSIBLE ADVERSE EFFECTS

- Tissue damage
- Infection
- Toxic Shock Syndrome

## ⚠

### WARNINGS

• The device is for SINGLE USE ONLY. DO NOT re-sterilize or re-use.

• **DO NOT** use if the package is opened or damaged or the date of expiration has passed. **DO NOT** use if the device is damaged.

• **DO NOT** leave packing in place longer than 24 hours.

- Removal must be conducted by a healthcare professional
- Always rehydrate with sterile saline or sterile water before removal.
- To minimise the risk of displacement, ensure device is attached to the patient using the nylon string


### TOXIC SHOCK SYNDROME (TSS)

In rare instances, the physicochemical conditions associated with nasal surgery, both with and without nasal packings, may present a risk of toxic shock syndrome (TSS). Warning signs of TSS include: sudden fever (usually 102°F (38.9°C) or more), vomiting, diarrhea, dizziness, fainting (or near fainting when standing up), and/or a rash that looks like sunburn. If these or other signs of TSS appear, remove the packing at once and seek medical attention immediately.

Rhino Rocket Nasal Packs exhibit no antimicrobial properties, they are not bacteriostatic towards pre-existing infections, nor do they prevent the occurrence of new infections. In the case of pre-existing infections, it is recommended that appropriate treatment be instituted.

## EN

### INSTRUCTIONS FOR USE

**NOTE:** Before opening, visually inspect the outer pouch and pouch seal for voids or openings.

**DO NOT USE** if the outer pouch or pouch seal is not intact.

#### 1. Before Insertion:

Verify the correct style and size for the desired procedure. Open the product only in a sterile field. Open by pulling and separating the pouch seal.

#### 2. Insertion:

Rhino Rocket® Nasal Packs expel the nasal pack from the applicator along the nasal floor until the entire packing is within the nose See Fig 1) without inserting the applicator into the nose.

**3. After Insertion:** If necessary, hydrate the outer edge of the packing using a syringe containing sterile saline or sterile water until the packing begins to expand (See Figure 2). Secure the string to prevent packing migration (See Figure 3).

**4. Before Removal:** Rehydrate the packing using a syringe with sterile saline or sterile water until the packing is fully saturated. Allow the packing to remain in place for 5-10 minutes after rehydration to facilitate removal.

**5. Removal:** Use forceps to grasp the pack and gently withdraw.

## Rx Only

### INTENDED LIFETIME OF THE DEVICE

Rhino Rocket® Nasal Packs are short term use devices that are in continuous contact with the nasal and/or sinus structures for up to 24 hours.

### SHELF LIFE

Rhino Rocket® Nasal Packs have a declared shelf life of 3 years.

### STERILITY ESTERILIDAD

Rhino Rocket is provided terminally sterilized by gamma irradiation. The product is sterile in an undamaged, unopened package.

### DISPOSAL

After use the device should be disposed of in accordance with the healthcare settings disposal procedure. The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

### INCIDENT REPORTING

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the FDA or the competent authority of the Member State in which the user and/or patient is established.

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