TISSU-TRANS® FILTRON®
100cc, 250cc, 500cc, 1000cc & 2000cc

Sterile Fat Collection and Filtration Instructions for Use

3-TT-FILTRON 100, 250, 500, 1000 & 2000

DEVICE DESCRIPTION

The Tissu-Trans® FILTRON® devices are designed to collect and filter autologous fat tissue in a sterile, closed system that can be re-injected into patients.

The Tissu-Trans® FILTRON® devices are manufactured in 5 different sizes.

Figure 1

STERILE R Tissu-Trans® FILTRON® is provided in protected packaging and is provided sterile.

The Tissu-Trans® FILTRON® comprises of:

- Harvest Canister with Filter and Waste Extraction Port
- Harvest Tubing 6 feet / 1.83m
- Transfer Tube 6 inches / 152mm with Luer Lock tip

The Tissu-Trans® FILTRON® is a single use device should not be re-sterilized.

Standard Liposuction Equipment Also Needed:

- Universal Bracket on tall canister stand or IV pole clamp
- Sterile mayo-stand drape or equivalent to cover canister stand
- Standard suction tubing (1/4" / 6.35mm or 5/16" / 7.94mm ID)
- Waste canister
- Suction source (set at approx. 15" Hg. (380mmHg) / 51 kPa or 50% power)
- Cannula (3mm cannula recommended)
- Luer Lock syringes for transfer and re-injection

Figure 2

The FILTRON is placed in-line between the cannula handle and waste canister

INTENDED USE

The Tissu-Trans® FILTRON® is intended to be used with appropriate vacuum systems and/or cleared pumps, waste tubing and cannulas, for the collection and filtration of autologous fat that can be re-injected for aesthetic and medical procedures, as deemed appropriate by the Healthcare Professional.

INDICATION FOR USE

The Tissu-Trans® FILTRON® is intended to be used with vacuum and/or cleared pumps, waste tubing and cannulas, for the collection and filtration of aspirated fat, for aesthetic and

medical procedures. The untreated fat, may be re-injected via a cleared injection apparatus.

CONTRAINDICATIONS

There are no known contraindications for the Tissu-Trans® FILTRON® device.

INTENDED USER

Healthcare professionals

PATIENT POPULATION

The intended patient population is adults.

CLINICAL BENEFITS

The clinical benefit of the Tissu-Trans® FILTRON® is to collect and filter autologous fat in a sterile, closed system that then can be re-injected.

Virtually pure fat retained in the filter as demonstrated in a study(1) to produce virtually pure fat containing trace amounts (1%) of residual oil and fluid.

(1)Fisher, C., Plast Reconstr Surg (2013) 132:351

POTENTIAL ADVERSE EFFECTS

There are no known potential adverse effects associated with the use of the Tissu-Trans® FILTRON® device.

CLINICAL PERFORMANCE CLAIMS

- Enables the collection, filtration and transfer of autologous fat for medical and aesthetic purposes
- Separates fat and fluids by filtration
- Filters and removes waste, non-viable cells and excess fluid during harvest
- No decant time required - fat is ready for re-injection

CLINICAL SAFETY CLAIMS

- Closed, sterile, single-use system
- Compatible with most aspirators and instruments
- Eliminates need for centrifuge(2)

(2)Torres, S., Aesthetic Medicine (2016) 2:33.

INTENDED LIFETIME OF THE DEVICE

Devices are single use, transient devices. Therefore, the device lifetime in use is less than 60 minutes continuous use.

SHELF LIFE OF THE DEVICE

The Tissu-Trans® FILTRON® devices have a shelf-life of three years. Use this product before the expiration date indicated on the label.

DISPOSAL

The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

All components should be disposed of in accordance with hospital and/or clinic protocol.

WARNINGS

1. Do not re-use or re-sterilize. Single patient use. Re-use may lead to cross contamination and patient infection.
2. This device will not, and of itself, produce significant weight reduction.
3. This device should be used with extreme caution in patients with chronic conditions.
4. The volume of blood loss and endogenous

body fluid loss may adversely affect intra and/or post-operative hemodynamic stability and patient safety. The capability of providing adequate, timely replacement is essential for patient safety.

- The device should not be exposed to suction pressure higher than 25 inHg /85 kPa.
- Use the appropriate size canister for harvesting the required amount of fat. Over filling the device can result in harvest tubing blockage sub-optimal filtration and a poorer quality fat.
- DO NOT Use if the device or its packaging has been damaged.
- Use of this device is limited to those physicians who have attained proficiency in suction lipoplasty.

- Do not use if package is damaged. Please contact the manufacturer with the product code and lot number listed on the device and dispose of the package in accordance with the "Disposal" section of this instruction for use.


SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the FDA/competent authority of the Member State in which the user and/or patient is established.

TALLOW-DERIVED SUBSTANCES

The harvest tubing may contain a tallow-derived substance that has been processed under conditions at least as vigorous as those specified in Annex 1, Section 3 of (EU) 722/2012. The use of these tallow derivatives in the harvest tubing do not pose a residual risk to patients.

Figure 3


INSTRUCTIONS

- Drape canister stand with sterile Mayo stand cover (or equivalent). NOTE: The outer pouch is the sterile barrier. All inner pouches and packaging materials are protective packaging.
- Remove sterile components and place on sterile surface. Remove and discard protective packaging.
- Secure Tissu-Trans® FILTRON® to canister stand in an upright position by inserting mounting bracket on canister into universal slots on stand - being careful to maintain sterility.
- Connect standard suction tubing to the VAC port on the canister lid. Connect the other end to waste canister and suction source.
- Connect L-Connector end of included harvest tubing to PATIENT port on the canister lid. Connect other end of tubing to cannula.

- Turn on suction. Check all connections. **Adjust vacuum pressure to 15"Hg (380mmHg) / 51 kPa or 50% power.** System may be tested with sterile saline to ensure proper flow.
- Harvest desired volume of fat. Virtually pure* fat will remain in filter and waste will be suctioned out to waste canister. Fat remaining in filter will be ready for re-injection (if applicable). NOTE: The volumetric markings on the canister are not calibrated and are intended for reference purposes only.

Filtered fat may still have some liquid, which is necessary to ensure proper flow during transfer. After fat is transferred to syringes, it may be decanted further and excess liquid expelled. (Fisher, C., Plast Reconstr Surg (2013) 132:351)

Figure 4

REMOVAL / TRANSFER

OPTION 1: DEVICE IS INVERTED.

- Turn off suction. Disconnect tubes from Tissu-Trans® FILTRON® lid, maintaining sterility.
- On lid, close VAC port with attached cap. All ports except PATIENT should now be closed.
- Attach the Transfer Tube to PATIENT transfer port. See Figure 5. Close the white pinch clamp on the Transfer Tube. All ports are now closed.
- Carefully invert the Tissu-Trans® FILTRON® upside down.
- Reattach inverted canister to the canister stand using the same bracket.
- Attach Luer Lock syringe to Transfer Tube. Open pinch clamp and draw out fat to be re-injected. See Figure 6. Repeat with additional syringes.
- When desired volume has been transferred, discard entire system.
- Attach re-injection cannula. Fat is ready to be re-injected to recipient site (if applicable).

Figure 5



Figure 6

REMOVAL / TRANSFER

OPTION 2 (500CC, 1000CC, & 2000CC ONLY): DEVICE REMAINS UPRIGHT.

1. Turn off suction. Disconnect tubes from Tissu-Trans® FILTRON® lid, maintaining sterility.
2. On lid, close both VAC and PATIENT ports with attached caps. All ports should now be closed.
3. Attach transfer or re-injection cannula to luer lock syringe.
4. Remove cap from large CAP1 Port.
5. Insert cannula into CAP1 Port and draw out fat to be re-injected. Repeat with additional syringes.
6. When desired volume has been transferred, discard entire system.
7. Attach re-injection cannula. Fat is ready to be re-injected to recipient site.

Figure 7

PERFORMANCE CHARACTERISTICS

- The Tissu-Trans® FILTRON® devices filter lipoaspirate to provide virtually pure fat* for autologous fat grafting
- Approximate volume of fat which can be filtered within the device
 - 100cm³ (3-TT-FILTRON® 100)
 - 250cm³ (3-TT-FILTRON® 250)
 - 500cm³ (3-TT-FILTRON® 500)
 - 1000cm³ (3-TT-FILTRON® 1000)
 - 2000cm³ (3-TT-FILTRON® 2000)
- Maximum suction pressure: 25"Hg / 85 kPa

*Demonstrated in a study to produce virtually pure fat containing trace amounts (1%) of residual oil and fluid. (Fisher, C., *Plast Reconstr Surg* (2013) 132:351)

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Summit Medical LLC
MDSS GmbH