Instructions for Use

TISSU-TRANS® FILTRON® 500cc, 1000cc, and 2000cc

Sterile Fat Collection and Filtration

3-TT-FILTRON 500, 1000 and 2000

Product is packaged in clear pouches and are provided sterile. DO NOT re-sterilize.

Product part number is provided on pouch label and indicates size.

Shelf Life

Use this product before the expiration date indicated on the label.

Intended Use

The Tissu-Trans® Filtron® is intended to be used with vacuum and/or cleared pumps, tubing and cannulas, for the collection of aspirated fat, for aesthetic body contouring. If the fat is untreated, it may be re-injected via a cleared injection apparatus.

Tissu-Trans Filtron is used in the aspiration, harvesting, filtering and transferring of autologous tissue.

Tissu-Trans Filtron is intended for use in the following surgical specialties when the aspiration of soft tissue is desired:

- Plastic and Reconstructive Surgery
- Neurosurgery
- Gastrointestinal and Affiliated Organ Surgery
- Urological Surgery
- General Surgery
- Orthopedic Surgery
- Gynecological Surgery
- Thoracic Surgery
- Laparoscopic Surgery

Indication for Use

The Tissu-Trans Filtron is intended to be used with vacuum and/or cleared pumps, tubing and cannulas, for the collection of aspirated fat, for aesthetic body contouring. If the fat is untreated, it may be re-injected via a cleared injection apparatus.

STERILE Sterile Items Included:

- Harvest Canister with Filter and Waste Extraction Port
- Harvest Tubing 6 feet
- Transfer Tube 6 inches with Luer Lock tip

Standard Items Also Needed:

- Universal Bracket on tall canister stand or IV pole clamp
- Sterile mayo-stand drape or equivalent to cover canister stand
- Standard suction tubing (1/4" or 5/16" ID)
- Waste canister
- Suction source (set at approx. 15" Hg. (380mm) or 50% power)
- Cannula
- Luer Lock syringes for transfer and re-injection

Notes Before Use:

1. The Tissu-Trans FILTRON is sterile and disposable. One patient use.

2. Do Not Re-sterilize or Re-use.

3. Product is ready for use. Do not open canister or lid.

3. Be certain that the tissue is well hydrated - as in superwet tumescent anesthesia.

4. **Must use** a standard waste canister between the FILTRON and suction source.

During harvest, waste is filtered out to waste canister.

5. **Recommended to use** a 3mm Harvest Cannula.

6. **Recommended to use** the medical grade tubing provided to assure non-toxic pathway.

7. Do not use if sterile pouch is torn or punctured.

8. If Filtron does not hold a vacuum, use a new device and call Customer Service.

Instructions:

1. Drape canister stand with sterile Mayo stand cover (or equivalent).
2. Remove sterile components and place on sterile surface. Remove and discard protective packaging.
3. Secure FILTRON to canister stand by inserting mounting bracket on canister into universal slots on stand - being careful to maintain sterility.
4. Connect standard suction tubing to the VAC port on the canister lid. Connect the other end to waste canister and suction.
5. Connect L-Connector end of included harvest tubing to PATIENT port on the cannister lid. Connect other end of tubing to cannula.
6. Turn on suction. Check all connections. **Adjust vacuum pressure to 15" or 380mm Hg.**
System may be tested with sterile saline to ensure proper flow.
7. Harvest desired volume of fat. Quality fat will remain in filter and waste will be suctioned out to waste canister. Fat remaining in filter will be ready for re-injection (if applicable).
Filtered fat may still have some liquid, which is necessary to ensure proper flow during transfer. After fat is transferred to syringes, it may be decanted further and excess liquid expelled.


Figure 2

Figure 3

Figure 4

Removal / Transfer Option 1: Device is inverted.

1. Turn off suction. Disconnect tubes from FILTRON lid, maintaining sterility.
2. On lid, close VAC port with attached cap. All ports except PATIENT should now be closed.
3. Attach the Transfer Tube to PATIENT transfer port. **See Figure 5.**
Close the white pinch clamp on the Transfer Tube. **All ports are now closed.**
4. Carefully invert the FILTRON upside down.
5. Reattach inverted canister to the canister stand using the same bracket.
6. Attach Luer Lock syringe to Transfer Tube. Open pinch clamp and draw out fat to be re-injected. **See Figure 6.**
Repeat with additional syringes.
7. When desired volume has been transferred, discard entire system.
8. Attach re-injection cannula. Fat is ready to be re-injected to recipient site (if applicable).

Figure 5

Figure 6

Removal / Transfer Option 2: Device remains upright.

1. Turn off suction. Disconnect tubes from FILTRON lid, maintaining sterility.
2. On lid, close both VAC and PATIENT ports with attached caps. All ports should now be closed.
3. Attach transfer or re-injection cannula to luer lock syringe.
4. Remove cap from large CAP1 Port.
5. Insert cannula into CAP1 Port and draw out fat to be re-injected. Repeat with additional syringes.
6. When desired volume has been transferred, discard entire system.
7. Attach re-injection cannula. Fat is ready to be re-injected to recipient site.

Figure 7

Dispose in standard medical biohazard waste receptacle.

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Figura 7

Warnings: 1. This device will not, in and of itself, produce significant weight reduction. 2. This device should be used with extreme caution in patients with chronic medical conditions. 3. The volume of blood loss and endogenous body fluid loss may adversely affect intra and/or postoperative hemodynamic stability.

Precautions: 1. This device is designed to remove localized deposits of excess fat through small incisions. 2. Use of this device is limited to those physicians who have attained proficiency in suction lipoplasty. 3. Results of this procedure will vary depending on patient age, surgical site, and experience of the physician. 4. Results of this procedure may or may not be permanent. 5. The amount of fat removed should be limited to that necessary to achieve a desired medical effect. 6. All disposable components must be discarded.

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Summit Medical LLC
MDSS GmbH