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
# MYRINGOTOMY BLADES AND HANDLE

**BEFORE USING, READ THE FOLLOWING INFORMATION:**

### INDICATIONS FOR USE

Where an incision in the tympanic membrane is required.

### CONTRAINDICATIONS

None known.

### PRECAUTIONS

The tip protector should remain in place until blade is ready to be used to prevent damage to the cutting edge or injury.

### POSSIBLE ADVERSE REACTIONS

- Minor scar tissue from incision.
- Infection due to airborne or water contaminants.

### INSTRUCTIONS FOR USE

#### 1. ASSEMBLY

##### A. Model BL-2024 Flat Blade Handle.

1. Check to insure that the handle and blade are compatible (See Chart Below).
2. Remove blade from packaging using sterile technique.
3. Loosen (but do not remove) locking nut.
4. Insert the blade into handle.
5. Tighten locking nut securely to retain blade.
6. Remove tip protector from end of blade and the Myringotomy Incision System is ready for use.

The Summit Medical Model BL-2024 Flat Blade Handle, is intended for use with the following models of myringotomy blades:

BL-2001	BL-2006
BL-2004	BL-2008

**Note:** No assembly required for BL-2002, BL-2005 and BL-2007.

#### 2. CLEANING (HANDLE ONLY)

Remove blade from handle before cleaning. The handles should be cleaned by scrubbing with a non-oily soap and soft brush. **DO NOT** use substances containing chlorides.

#### 3. STERILIZATION (HANDLE ONLY)

It is recommended that the handle be sterilized by steam autoclave, (prevacuum) 4 minutes 270°F (132°C).

(BL-2001, BL-2002, BL-2004, BL-2005, BL-2006, BL-2007, BL-2008) Single-use, one-patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of the patient.

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Manufacturer  
Fabricant  
Fabricante

Non-Sterile (For BL-2024 Only)  
Non stérile (pour BL-2024  
uniquement)  
No estéril (Para BL-2024  
únicamente)

Not made with natural rubber latex.  
Sans latex naturel.  
No está fabricado con látex de caucho  
natural.

Summit Medical LLC  
815 Vikings Parkway, Suite 100  
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939  
F: 1-888-229-1941 | +1 651-789-3979  
www.summitmedicalusa.com