- (EN) DISSOLVABLE VISIPLUG® LACRIMAL PLUG** Instructions for Use
(CS) ROZPUSTNÁ LAKRIMÁLNÍ ZÁTKA VISIPLUG® Návod k použití
(DA) OPLØSELIG VISIPLUG® TÅREPROP brugsanvisning
(NL) VISIPLUG® OPLOSBAAR TRAANBUISSSENT Gebruiksaanwijzing
(FI) SULAVA VISIPLUG®-KYYNELTIETULPPA - Käyttöohjeet
(FR) BOUCHON LACRYMAL RÉSORBABLE VISIPLUG® Instructions d'utilisation
(DE) LÖSLICHER VISIPLUG® TRÄNENWEGVERSCHLUSS Gebrauchsanweisung
(EL) ΔΙΑΛΥΟΜΕΝΟ ΔΑΚΡΥΪΚΟ ΒΥΣΜΑ VISIPLUG® Οδηγίες Χρήσης
(IT) TAPPO LACRIMALE SOLUBILE VISIPLUG® Istruzioni per l'uso
(LV) IZŠĶĪDINĀMS VISIPLUG® LAKRIMĀLAIS SPRAUDNIS Lietošanas instrukcija
(LT) TIRPSTANTIS „VISIPLUG®“ AŠARŲ TAKŲ KAMŠTUKAS Naudojimo instrukcija
(NO) OPPLØSELIG VISIPLUG® TÅREKJERTELPLUGG Bruksanvisning
(PL) ROZPUSZCZALNA ZATYCZKA DO KANALIKÓW ŁZOWYCH VISIPLUG®
Instrukcja obsługi
(PT) PLUG LACRIMAL DISSOLÚVEL VISIPLUG® Instruções de Utilização
(ES) TAPÓN LAGRIMAL SOLUBLE VISIPLUG® Instrucciones de uso
(SR) RASTVORLJIVI VISIPLUG® ČEPIĆ ZA SUZNI KANAL Uputstva za upotrebu
(SV) VISIPLUG® - UPPLÖSBAR TÅRKANALSPLUGG Bruksanvisning
(TR) ERİYEBİLEN VISIPLUG® LAKRİMAL TİKAÇ Kullanım Talimatları
(KO) 용해성 VISIPLUG® 누점 플러그 사용 지침







TABLE OF CONTENTS / OBSAH / INNHOLDSFORTEGNELSE / INHOUDSOPGAVE / SISÄLLYSLUETTELO / TABLE DES MATIÈRES / INHALT / ΠΙΝΑΚΑΣ ΠΕΡΙΕΧΟΜΕΝΩΝ / INDICE / SATURA RĀDĪTĀJS / TURINYS / INNHOLDSFORTEGNELSE / SPIS TREŚCI / ÍNDICE / TABLA DE CONTENIDOS / SADRŽAJ / INNEHÅLL / İÇİNDEKİLER / 목차

EN	English	3
CS	Česky	5
DA	Dansk	7
NL	Nederlands	9
FI	Suomi	11
FR	Français	13
DE	Deutsch	15
EL	Ελληνικά	17
IT	Italiano	19
LV	Latviešu	21
LT	Lietuvių k.	23
NO	Norsk	25
PL	Polski	27
PT	Português	29
ES	Español	31
SR	Srpski	33
SV	Svenska	35
TR	Türkçe	37
KO	한국어	39

Read Before Use

INTRODUCTION

The Dissolvable VisiPlug® Lacrimal Plug is an ophthalmic device commonly referred to as a lacrimal plug, placed by a practitioner into the lacrimal duct to restrict the natural lubricating tears from being pumped off the eye. Lacrimal Plugs are prescribed for temporary treatment of certain eye conditions collectively referred to as 'Lacrimal System Dysfunction' or 'Dry Eye Syndrome' and other conditions of tear insufficiency.

Dissolvable Plugs:

- Are intended to treat Dry Eye Syndrome
- Are effective in the horizontal canaliculi
- Never touch the eye
- Are comfortable after proper placement, and will not fall out of the punctum

Punctal dilation and topical anesthesia are not usually required for insertion.

INDICATIONS FOR USE

Dissolvable VisiPlug Lacrimal Plugs may be used:

- As a diagnostic aid to determine the potential effectiveness of Occlusion Therapy with non-dissolvable plugs.
- To temporarily enhance the efficacy of topical medications or ocular lubricants.
- After ocular surgery to prevent complications due to dry eyes.
- To evaluate treatment of ocular dryness secondary to contact lens use.
- In the treatment of Dry Eye Syndrome and the dry eye components of varying Ocular Surface Diseases.

INTENDED USER & PATIENT TARGET GROUP

The devices are intended to be inserted by suitably qualified ophthalmic surgeons/physicians. The patient population is determined by the surgeon/physician in accordance with the intended use.

CONTRAINDICATIONS

This device is not designed, sold, or intended for use except as indicated or prescribed by a physician. Below are examples where Dissolvable VisiPlug Lacrimal Plugs may not be the therapy of choice:

- Chronic tearing (epiphora) secondary to canalicular obstruction
- Dacryocystitis with or without mucopurulent discharge.

PACKAGING AND STORAGE CONDITIONS

Dissolvable VisiPlug Lacrimal Plugs are made of polydioxanone and are provided in a Tyvek pouch which acts as the sterile barrier, two plugs per pouch. After sterilization the devices are packaged inside a moisture barrier pouch with a desiccant as the plugs are moisture sensitive. To achieve the stated shelf life and performance claims, Visiplugs must be stored within the moisture barrier pouch until such time as they are ready to be used.

INTENDED LIFETIME OF THE DEVICE

Dissolvable VisiPlug Lacrimal Plugs have a medium term degradation rate of approximately 180 days.

BEFORE INSERTION

Dissolvable VisiPlug Lacrimal Plugs are available in two sizes: 0.4mm & 0.5mm. The 0.4mm plugs are appropriate for use with most adult patients; use 0.5mm when 0.4mm plugs will not stay in position.

A physician must use medical judgement and consider a patient's medical history prior to implanting a Dissolvable VisiPlug Lacrimal Plug. Patients with pre-existing intermittent tearing should receive pressure irrigation to rule out canalicular obstruction.

INSERTION

- Inspect the patient's punctum to determine which size plug to use.
- Remove a sterile packet from the box, peel open the pouch and remove the foam holder.
- While using magnifying loupes or a slit-lamp, have the patient look away from the point of insertion.
- Use a jeweler's forceps to remove one of the plugs from between the grooves in the foam holder.
- Apply traction with a cotton tip applicator to evert the eyelid and clearly view the punctum.
- Guide the plug partially into the punctum and release the forceps (Fig. 1).
- Apply lateral traction with the cotton tip applicator to straighten the angle between the vertical and horizontal canaliculus (Fig. 2).
- Use the tips of the forceps to push the plug out of sight into the horizontal canaliculus (Fig. 3).
- Blinking and normal tear flow cause the plug to migrate down into the horizontal canaliculus (adjacent to the common canaliculus).
- Repeat the same procedure for each remaining punctum. After insertion, inspect each punctum to ensure the plugs are not sticking out (Fig. 4).
Unused plugs should be discarded.

Plug insertion (Fig.1-4) is the same for inferior and superior puncta.

FOLLOWING PLACEMENT

Patients should feel little or no discomfort following plug placement.

For patients experiencing irritation or epiphora after plug insertion, use a sterile probing device to confirm proper placement within the horizontal canaliculus (4-6mm beneath the punctum).

Use a single drop of topical antibiotic to prevent complications related to plug insertion.

WARNINGS

During insertion be careful not to perforate the canaliculus with the jeweler's forceps. Perforation can increase the risk of infection. If perforation occurs, delay plug insertion until the wound heals.

 This device is designed, intended and distributed for single use only. **DO NOT RE-STERILIZE OR REUSE THE DEVICE.** There are no data to support the sterility and functionality of the device after reprocessing.

 Do not use if package is damaged. Please contact the manufacturer with the product code and lot number listed on the device and dispose of the package in accordance with the "Disposal" section of this instruction for use.

VisiPlugs are moisture sensitive devices and will begin degradation when exposed to ambient air. Premature degradation of the lacrimal plugs may decrease the length of treatment.

IMPLANT REMOVAL

Pressure irrigation (Fig. 5) is the preferred method for plug removal. In rare cases, when this method proves unsuccessful, surgical removal may be required.

DISPOSAL

VisiPlug devices leave the body through the digestive tract via natural waste. Any unused product should be disposed of via clinical waste in accordance with local procedures. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

- The Dissolvable VisiPlug Lacrimal Plugs provide temporary occlusion of the tear drainage system and have a medium term degradation rate of approximately 180 days.
- A summary of safety and clinical performance can be found at the following link, when available via EUDAMED: <https://ec.europa.eu/tools/eduamed/#/screen/home>.

POSSIBLE ADVERSE EFFECTS

- Infection
- Canaliculitis
- Tearing or Perforation of the canaliculus
- Eye watering or excessive tearing (epiphora)
- Irritation, itching, swelling, or pain
- Need to surgically remove plug that doesn't flush out of canaliculus

MATERIALS: POLYDIOXANONE AND D&C VIOLET #2 (<0.25%)

STERILITY

Lacrimal plugs are provided terminally sterilized by ethylene oxide (EO) gas. The product is sterile in an undamaged, unopened package.

Adhesive product labels are included on the pouch for use on patient and hospital records.

INCIDENT REPORTING


Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the FDA or the competent authority of the Member State in which the user and/or patient is established.

PATIENT IMPLANT CARD: INSTRUCTIONS FOR COMPLETION

All steps are to be completed by the healthcare provider.

1. Implantation date
2. Patient name or identification number
3. Healthcare institution/provider name and address
4. Remove device information sticker from pouch label and adhere to this card in the space provided
5. Provide patient with completed implant card, which provides a link to the patient information leaflet

**FOR PROFESSIONAL USE ONLY.
RX ONLY**

Summit Medical LLC
815 Vikings Parkway, Suite 100
St. Paul, MN 55221 | USA
www.summitmedicalusa.com

1 / /
2
3

summitmedicalusa.com/patientinfo

EN DISSOLVABLE LACRIMAL PLUG IMPLANT CARD
CS KARTA IMPLANTÁTU ROZPUSTNÝCH UZÁVERU SLZŇNÝCH CEST
DA IMPLANTÁTKORT TIL OPPLØSELIG TÅREKIRTELPROP
NL IMPLANTAATKAART OPLOSBARE LACRIMALE PLUG
FI SULAVAN KYNNELTIETULPAN IMPLANTTIKORTTI
FR CARTE D'IMPLANT DU BOUCHON LACRYMAL RESORBABLE
DE IMPLANTÁTKARTE FÜR LOSLICHE TRÄNENKANALSTÖPSEL
EL KARTA EMOYEMIAZOZ DĀNYOMENOY ĀKPYKOY BYZMATOZ
IT SCHEDA DI IMPIANTO TAPPO LACRIMALE SOLUBILE
LV IZSKIDINĀMA LAKRIMĀLĀ SPRĀIDNA IMPLANTĀTA KARTE
LT TIRPIOJO AŠARŲ LATAKO KAISCIO IMPLANTO KORTELE
NO IMPLANTÁTKORT FOR OPPLØSELIG TÅREKANALPLUGG
PL KARTA IMPLANTU ROZPUSTZALNEJ ZATYCKI DO KANAŁKÓW ŁZOWYCH
PT CARTÃO DE IMPLANTE DE PLUG LACRIMAL DISSOLÚVEL
ES TARJETA DE IMPLANTE DE TAPON LAGRIMAL SOLUBLE
SR KARTICA IMPLANTATA ZA RASTVORLJIVI ČEPIĆ ZA SUZNI KANAL
SV ERIYEBILEN LAKRIMAL TIKÅC IMPLANT KARTI
TR IMPLANTÁTKORT FÖR UPPLÖSBAR TÅRPROPP
KO 분해성 눈점 플러그 임플란트 카드

267381 REVA | JULY2022

TO ORDER, REQUEST INFORMATION, OR REPORT PRODUCT PERFORMANCE, CONTACT YOUR LOCAL DISTRIBUTOR, OR:

Summit Medical LLC.
ATTN: Customer Service
815 Vikings Parkway, Suite 100
St. Paul, MN 55121
Tel: 1.888.229.2875 | 651.789.3939
Fax: 1.888.229.1941 | 651.789.3979
E-mail: customerservice@innoviamedical.com
WWW: <https://summitmedicalusa.com/>

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