Summit Medical Internal Nasal Splints

i Before using, read the following information:

DESCRIPTION

Internal Nasal Splints are designed to splint the septum following surgery or trauma.

INTENDED USE

The Internal Nasal Splints are intended to splint the septum following surgery or trauma.

INDICATIONS FOR USE

The Internal Nasal Splints are indicated for the temporary, internal splinting of the nasal septum.

TARGET PATIENT GROUP

Any adult in the human population as determined by a suitably qualified healthcare professional.

INTENDED USER

The Internal Nasal Splints are intended to be used by appropriately trained healthcare professionals.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

The following clinical benefits and performance characteristics are anticipated with the internal nasal splints:

- Designed to reduce complications associated with nasal septal surgery or trauma (e.g., adhesion, obstruction, bleeding, crusting, and synechia) Anatomically designed to reduce patient discomfort on insertion, while in situ and on removal.
- Designed to be sutured in place to reduce the possibility of accidental displacement.

CONTRAINDICATIONS

The Breeze™ Nasal Splints with Airway and Breeze™ Nasal Splints without Airway are made of silicone which may cause an allergic reaction in patients with sensitivity to this material. User should risk assess the use of the device in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to silicone.

Breeze™ Nasal Splint with Airway

Airway geometry creates a pressure gradient force that increases airflow.

Bivalve Nasal Splint

A medial slit aids insertion and removal

Available in standard or oversized with thick or thin options

PRECAUTIONS

This device is not designed, sold or intended for use except as indicated or prescribed by a physician. Intended for single use only, do not re-sterilize. Do not use if pouch is damaged.

POSSIBLE ADVERSE EFFECTS

- Minor scar tissue
- Localized infection
- Granuloma

WARNINGS

- The device is SINGLE USE ONLY. DO NOT re-sterilize or re-use.

- DO NOT use if the package is opened or damaged or the date of expiration has passed.
- DO NOT use if the device is damaged

INSTRUCTIONS FOR USE:

- Insert one internal nasal splint into each of the patient’s nasal cavities.
- Suture the internal nasal splint to the nasal septal wall.
- During the wear time, nasal irrigation may be used to prevent blockages through the airway tube.
- Typically sutures and internal nasal splints are removed after two to ten days but can be left in-situ for up to 30 days.

INTENDED LIFETIME OF THE DEVICE

The Internal Nasal Splints are intended for short term use and can be left in situ for up to 30 days.

STERILITY **STERILIZED**

Internal Nasal Splints are provided terminally sterilized by ethylene oxide (ETO). The product is sterile in an undamaged, unopened package.

MATERIALS

Bi-Valve Nasal splints: Fluoroplastic (FEP)

Breeze™ Nasal Splints: Medical Grade Silicone

DISPOSAL

After use the Internal Nasal Splints should be disposed of in accordance with the healthcare settings disposal procedure. The Internal Nasal Splints come into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the Internal Nasal Splints after use to prevent contamination.

INCIDENT REPORTING

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the FDA or the competent authority of the Member State in which the user and/or patient is established.Member State in which the user and/or patient is established.

- X** Single-use, one-patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of the patient.

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