Summit Medical Internal Nasal Splints

i

Before using, read the following information:

DESCRIPTION

Nasal splints are made of silicone or fluoroplastic and designed to fill splinting needs following various nasal surgical procedures.

INTENDED USE

The internal nasal splints are intended to be used postoperatively on septal surgery patients to reduce complication rates.

INDICATIONS FOR USE

Where temporary, post-surgical splinting of the nasal septum is desired.

TARGET PATIENT GROUP

The patient population is determined by the surgeon/clinician in accordance with the intended use.

INTENDED USER

This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.

EXPECTED CLINICAL BENEFITS AND PERFORMANCE CHARACTERISTICS

The following clinical benefits and performance characteristics are anticipated with the internal nasal splints:

- Reduction in rate of postoperative complications associated with nasal septal surgery (e.g., adhesion, obstruction, bleeding, crusting, and synechia)
- Reduction in patient pain/discomfort compared to nasal packing

CONTRAINDICATIONS

None known. This device is not designed, sold, or intended for use except as indicated or prescribed by a physician.

PRECAUTIONS

This device is not designed, sold or intended for use except as indicated or prescribed by a physician. Intended for single use only, do not resterilize. Do not use if pouch is damaged.

POSSIBLE ADVERSE EFFECTS

- Minor scar tissue
- Localized infection
- Granuloma

WARNINGS

- Nasal splints are to be sutured into place to prevent aspiration or movement.
- The device is for SINGLE USE ONLY. DO NOT re-sterilize or re-use.

- Prior to use, check the product expiration date and verify the sterile packaging has not been compromised. DO NOT use if the package is opened or damaged or the date of expiration has passed. DO NOT use if the device is damaged.

- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.

INSTRUCTIONS FOR USE:

- The nasal splint is inserted into the nose and positioned as needed.
- The nasal splint is sutured to the nasal septal wall.

INTENDED LIFETIME OF THE DEVICE

The Nasal Splints are intended for short term use and should not be left in place for more than 30 days.

STERILITY ^{STERILIZED}

Internal Nasal Splints are provided terminally sterilized by ethylene oxide (ETO). The product is sterile in an undamaged, unopened package.

MATERIALS

Bi-Valve Nasal splints: Fluoroplastic (FEP)
Breeze™ Nasal Splints: Medical Grade Silicone

DISPOSAL

Used devices should be disposed of as clinical waste in accordance with local procedures. If unused, product in original packaging which has exceeded the declared shelf life should be removed from packaging and disposed of as clinical waste.

INCIDENT REPORTING

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the FDA or the competent authority of the Member State in which the user and/or patient is established.

- Single-use, one-patient device will degrade if reprocessed. No effective cleaning process has been developed to prevent cross contamination. Contamination of a reprocessed device may lead to injury, illness or death of the patient.

Summit Medical Attelle nasale interne

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