

# Face Shield

Our full Face Shield is 9" in length and constructed using a P.E.T. plastic shield, polyurethane foam strip, and non-latex elastic headband for comfort.



## PROTECT YOUR PATIENTS & EMPLOYEES

- Clinics
- Hospitals
- Veterinary Practices
- Ambulatory Surgery Centers



## PRODUCT INFORMATION

### ORDERING INFORMATION

PART NUMBER	QTY
CM-5001-25	25
CM-5001	100

*\*Custom labels available for large orders*

### SPECIAL FEATURES

- FDA compliant PPE
- Full face protection
- Proudly made in the USA
- Lightweight and comfortable for long shifts
- Helps protect against splashes, sprays, and bodily fluids
- Non-latex elastic band
- Single-Use

## SPECIFICATIONS

**FOAM STRIP**  
9 in (L) x 1 in (W) x 1 ½ in (H)

**ELASTIC BAND**  
14 in (L) x ¾ in (H)



**ANTI-FOG P.E.T. PLASTIC SHIELD**  
13 in (L) x 9 in (H)

## CONTACT US

 1.888.229.2875

 [customerservice@innoviamedical.com](mailto:customerservice@innoviamedical.com)