

# Face Shield

Face Shields are a great alternative to face masks. They allow for better comfort over long shifts and the ability to more clearly communicate with customers. Face shields can also be worn in addition to masks for added protection.



## PROTECT YOUR CUSTOMERS & EMPLOYEES

- Health Care
- Office
- Retail
- Industrial
- Everyday Use



## PRODUCT INFORMATION

### ORDERING INFORMATION

PART NUMBER	QTY
CM-5001-02	2
CM-5001-25	25
CM-5001	100

*\*Custom labels available for large orders*

### SPECIAL FEATURES

- FDA compliant PPE
- Proudly made in the USA
- Lightweight and comfortable for long shifts
- Unobstructed communication and hearing
- Breathable for active jobs
- Effective protection for people with beards
- Single-Use
- Full face protection
- Non-latex elastic band

## SPECIFICATIONS

### FOAM STRIP

9 in (L) x 1 in (W) x 1 ½ in (H)

### ELASTIC BAND

14 in (L) x ¾ in (H)



### ANTI-FOG P.E.T. PLASTIC SHIELD

13 in (L) x 9 in (H)

## CONTACT US

 1.888.229.2875

 [customerservice@innoviamedical.com](mailto:customerservice@innoviamedical.com)