



# INSTRUSAFE

CARE + MAINTENANCE PRODUCTS

## Instrument Identification Tape








The Instrument Identification Tapes are intended to identify surgical instruments through processing and use. The Tape material is vinyl with acrylic adhesive.

They have been evaluated for the following sterilization methods:

- Conventional Steam Autoclave Cycle for 121°C for 30 minutes
- Steam Autoclave Cycle at 132°C for 3 minutes
- Conventional Hospital Ethylene Oxide Cycle

The tape is gas and steam permeable.

### Symbol Reference Key

	Consult Instructions For Use
	Quantity
	Catalog Number
	Lot Number
	Non-Sterile
	For Us Audiences Only
	Date Of Manufacture

Distributed by:



815 Vikings Parkway, Suite 100  
St. Paul, MN 55121 | USA

P : 1-888-229-2875 | 651-789-3939  
F : 1-888-229-1941 | 651-789-3979

[www.summitmedicalusa.com](http://www.summitmedicalusa.com)