



INTRUSAFE

INSTRUMENT PROTECTION TRAYS

Your Trusted Leader for Instrument Longevity.



Innovia Medical is a **family of specialty surgical companies** that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.















We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com.

THE INSTRUSAFE ADVANTAGE

360° OF PROTECTION	4
INDUSTRY DISADVANTAGES	6
FEATURES + BENEFITS	8
ORDER PROCESS: THE INSTRUSAFE WAY	10
STERILIZATION BARRIERS + FOOT OPTIONS	12

TRAYS BY PROCEDURE

 ARTHROSCOPY	14
 BARIATRIC	16
 CARDIOVASCULAR	18
 DENTAL	20
 ENT	22
 GENERAL	34
 GYNECOLOGY	38
 NEUROLOGY/SPINE	40
 OPHTHALMOLOGY	46
 ORTHOPEDIC	50
 PLASTICS	52
 RETRACTORS	54
 ROBOTICS	58
 SCOPES	79

TRANSPORT CONTAINERS

XL TRANSPORT CONTAINERS	90
HEAVY DUTY TRANSPORT CONTAINERS	92
TRANSPORT CONTAINERS	93

ACCESSORIES

SILICONE SMALL PARTS BOXES	94
PLAQUE OPTIONS	95
FINGER MATS	96
GUARDIAN TRAY	96
INSTRU-ZYME	97
ROUNDERS POINT GUARDS	97

RESOURCES



SUBMIT PURCHASE ORDERS TO:

Summit Medical LLC
815 Vikings Parkway, Suite 100
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939
F: 1-888-229-1941 | +1 651-789-3979
www.summitmedicalusa.com



The InstruSafe® Advantage

360° of Protection

The use of complex and delicate instrument sets is a new reality in which sterile processing and the operating room live. Intricacy of surgical instrumentation continues to increase, and so does the cleaning and sterilization process. Now, healthcare facilities face the challenge of reprocessing instruments safely and efficiently while also meeting the demand for fast turnaround. Keeping these sets organized and protected upon arrival in the operating room creates a challenge for all involved. Relieve yourself of these daily obstacles with InstruSafe Trays—an everyday solution for the intricacies of your instrument sets.

We have spent over 30 years understanding the use cycle and improving outcomes. With the use cycle in mind, we have created a product that eliminates common frustrations associated with instrument sets. The products in this catalog are a compilation of our hard work. As you journey through the pages, let us help you solve your everyday instrument problems.

Summit Medical's most recent FDA 510(k) clearances and product advancements are a testament of our dedication to our customers' satisfaction. We wouldn't have continuing success without the healthcare community we serve. Focus, determination and your support will continue our impact on patient care for decades to come.




Designed to organize and protect from the point of operating room use to decontamination, our solutions benefit everyone who handles our products.



FDA 510(k) Clearance

Designed to ensure optimal sterilization, InstruSafe Trays have undergone rigorous testing before receiving multiple FDA 510(k) clearances. Cleared for various sterilization cycles, these reliable trays give hospitals and healthcare facilities even more options when it comes to organizing, sterilizing, and protecting their instruments.

For the most up-to-date list of cleared cycles, visit www.instrusafe.com/validations.



POINT CAMERA TO
SEE VALIDATIONS

**DECONTAMINATION**

Instruments are individually protected going into the wash cycle and the trays are designed with large perforations, allowing maximum detergent and water flow.

TRAY PREPARATION

Each tray is designed specifically for the instruments within it, providing individual locations for each one. This protects them from harming the technician, from damaging themselves and from puncturing or jeopardizing the integrity of the sterile barrier.

STERILIZATION

The large amount of perforations in our trays provide superior sterilant penetration. Many of our silicone parts are also designed with minimal contact points and additional ventilation in each slot that match our tray holes to improve sterilant penetration to each individual instrument.

TRANSPORTATION

Soft and durable silicone holders in the base and cover secure instruments to prevent shifting and provide 360 degrees of protection. InstruSafe Trays help to avoid unnecessary damage and to present an organized instrument set every time.

OPERATING ROOM USE

Each instrument is easily identifiable and placed in a singular slot for surgical presentation and tracking at the end of a procedure.


Look Familiar?

Unfortunately, this is what most sterile processing departments have come to know. The familiarity and acceptance of these messy, hazardous baskets makes us cringe! Why, you ask? Because the minimal organization and protection of surgical instruments practiced in many facilities presents numerous risks to the individuals that come in contact with them, as well as to the instruments themselves.



DISADVANTAGES OF MESSY BASKETS

- POTENTIAL TO COMPROMISE STERILE BARRIER BY POKING THROUGH WRAP □
- DANGEROUS TO TRANSPORT FROM THE O.R. TO THE SPD □
- MINIMAL INSTRUMENT PROTECTION AND ORGANIZATION □
- DIFFICULT TO COUNT AND TRACK INSTRUMENTS □
- BASKETS NEED FREQUENT REPLACEMENT □
- ADDS TIME TO SURGICAL SET-UP □
- BROKEN INSTRUMENTS IF DROPPED □
- HIGH RISK FOR FINGER POKES □
- POOR SURGICAL PRESENTATION □
- HARD TO FIND WHAT YOU NEED □
- MISSING/LOST INSTRUMENTS □
- UNNECESSARY REPROCESSING □
- DAMAGED INSTRUMENTS □
- WASTE OF SPACE □



EXPECT MORE.

For these reasons and many more, InstruSafe Trays offer custom solutions to meet each customer's specific needs. We think it's time for everyone to expect more from their trays.


Absolute Protection

Trusted Solutions + Superior Results




QUARTER TURN LATCH

Secures the tray cover to the tray base, providing visual indication of the open/close position. Our off-set latch hole aligns the cover in the correct position ensuring the silicone holders match up with the instruments below.



RETRACTABLE HANDLES

Our multi-functional space-saving design fits into rigid containers making it effortless to remove the tray.




TF FOOT

Our Tall Foot fits into rigid containers providing clearance for the filter in the base.


FOOT OPTIONS


Sterile barrier and storage needs determine the foot type built into a tray.



R FOOT

Rounders Corner Guards prevent tears in sterilization wrap during storage and transportation.





POINT CAMERA
TO READ IFU

PERFORATIONS (GRID PATTERN)

Each tray contains a large quantity of perforations that provide maximum sterilant penetration and ventilation.

VENTILATION

Our unique silicone design has channels that align with our perforations allowing detergent and sterilant to flow through the holder and access the entire instrument.

SILICONE HOLDERS

InstruSafe silicone holders are soft yet durable, with a history of longevity. Our variety of holders accommodate instruments large or small, allowing each instrument to have its own location within the tray.

Holders in the tray base and cover provide 360 degrees of protection reducing instrument repair and replacement while delivering instruments to the O.R. in their optimal working condition.

REDUCED CONTACT POINTS

Each silicone holder is individually designed with ridges and channels built to reduce contact with the instrument while securely holding the instrument in place.

5 Simple Steps

If one of our preconfigured tray designs does not fit your needs let our product experts design a custom tray with our effortless process below.




1 ASSESS

An InstruSafe Representative conducts a free, on-site assessment of your instrument set.




2 GATHER

Your instrument set specifics and inputs are gathered by our team of experts.



3 DESIGN

Using your inputs, our product specialists build a tray design with you and your instrument set in mind.





EVALUATE

You be the judge! Evaluate the new tray design while it's loaded with your instrument set. Let us know your thoughts and we will make any and all adjustments necessary!



COMPLETE

It's all yours! Once complete, each custom tray is assigned its own unique part number, making it simple and easy for you to reorder in future.

HOW TO USE THE CATALOG


The colored headers identify which procedure the instrument set you're looking at is used for. The trays on each page are differentiated by preconfigured or custom configured with the symbols shown in the key below. Some trays have several layers including a cover, shelves, baskets, and a base. The custom configured tray examples included in the catalog are a small representation of the trays our local representatives have built for facilities around the world.

CUSTOM + PRECONFIGURED SYMBOLS




THIS SYMBOL IDENTIFIES A PRECONFIGURED TRAY

Order by part number from your local InstruSafe representative or call us today.



THIS SYMBOL IDENTIFIES A CUSTOM CONFIGURED TRAY EXAMPLE

Call us today to schedule your free assessment.



Sterilization Barriers







InstruSafe Instrument Protection Trays are FDA cleared for several sterilization cycles using wrap and/or rigid containers.

For the most up-to-date list of FDA 510(k) cleared cycles, visit www.instrusafe.com.

Foot Options

R Foot: Rounders Corner Guards—to be used with wrap

TF Foot: Tall Foot fits into rigid containers

STERILE BARRIER		FOOT OPTION
<div>WRAP</div> <div></div>		<div>Rounders Corner Guards -R</div> <div></div>
<div>CONTAINER</div> <div></div>		<div>Tall Foot -TF</div> <div></div>

HOW TO ORDER:

EXAMPLE: TRAY IN-2890

Wrapped Foot IN-2890- R	Container Foot IN-2890- TF
-----------------------------------	--------------------------------------

FOOT OPTION IS ADDED TO THE END OF THE PART NUMBER

Rounders® Corner Guards

Rounded, medical grade corner guards affix to InstruSafe Tray corners to reduce potential damage to the sterile wrap barrier. Rounders are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to all InstruSafe Trays.



ROUNDERS CORNER GUARDS

Small Rounders

Kit of Four

IN-7334

Large Rounders

Kit of Four

IN-7337

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit all InstruSafe Trays
- Rounded, medical grade silicone corner guard
- Help avoid compromised sterility in the O.R.
- Environmentally and economically efficient
- Impact resistant
- US Patent #9119611 | #D783856

Tall Foot

Designed to affix to trays that are intended to be sterilized within a rigid container. Tall Feet are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to most InstruSafe Trays.



TF FOOT KIT

Kit of Four

IN-7338


- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit most InstruSafe Trays



Preconfigured Tray for
Arthroscopy

This three-level tray is equipped for both small and large arthroscopic instrument sets. Our tray is sized to fit most three-quarter size containers or can be wrapped for your convenience.


Cover




Shelf



Basket



Base



ARTHROSCOPY INSTRUMENT TRAY

Arthroscopy Instrument Tray


Holds:
A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)	387 L x 244 W x 146 H (mm)
Wrapped Foot IN-2846-R	Container Foot IN-2846-TF

Expand capacity with positionable Hold-Its and Hold-Downs™ silicone instrument holders.
Multi-level to accommodate a high quantity of instruments.

Basket and base can be customized to accommodate specific instrumentation.





Custom Configured Tray Example for Arthroscopy Example

ARTHROSCOPY TMJ TRAY | 17 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Preconfigured Tray for Bariatric

The bariatric tray is designed to accommodate even the longest instruments. Our three-level system allows you to get organized and protect all of your bariatric instrumentation in one tray.



BARIATRIC TRAY

Bariatric Instrument Tray

Holds:
A combination of up to 15 (32cm to 45cm+) of the following instruments:

- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars

28 L x 13 W x 8 H (in)

711 L x 330 W x 203 H (mm)

Wrapped Foot
IN-8864-R


Container Foot
IN-8864-TF

Expand capacity with positionable Hold-Its and Hold-Downs™ silicone instrument holders.


Multi-level to accommodate a high quantity of instruments.

CALL US TODAY FOR CUSTOMIZATION!

CUSTOM CONFIGURED BARIATRIC TRAY | 31 INSTRUMENTS




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





Custom Configured Tray Example for Bariatric

BARIATRIC LAPAROSCOPIC SET | 59 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Cardiovascular

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

MICROVASCULAR SET | 19 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for
Cardiovascular

EVANS ROTATABLE TRAY | 32 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Dental

ARTHO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

DENTAL SET | 34 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Dental

OPERATIVE TRAY SET | 11 INSTRUMENTS

CUSTOM
CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ENT SYSTEMS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

Since the mid-1980's, we have offered a complete solution for fine and delicate ear instrumentation through our Middle Ear System. We are proud to claim this as our very first InstruSafe product. Our modules are designed to specifically fit each instrument and can be used individually or in combination, to organize and protect an entire instrument set.

MIX AND MATCH OUR MODULES FROM PAGES 23-24 TO FIT YOUR INSTRUMENT SET, CREATING YOUR UNIQUE ENT SYSTEM!




IN-2004




IN-2001



IN-2003



IN-2008



Tray shown without cover for illustrative purposes.



Preconfigured Trays for ENT Modules



□ MIDDLE EAR TRAY

Middle Ear Tray - 2 Level

Holds:

- 2 Alligator Forcep Modules
- 2 Pick + Hook Modules
- 1 Micro-Suction Module
- 1 Speculum Module
- 9½ x 13" Finger Mat

19 ⅝ L x 9 ⅝ W x 6 H (in)

498 L x 244 W x 152 H (mm)

Wrapped Foot
IN-2890-R

Container Foot
IN-2890-TF



□ EAR MODULE TRAY

Alligator Module


Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

5 ⅞ L x 5 ⅞ W x 1 H (in)

149 L x 130 W x 25 H (mm)

IN-2001



□ EAR MODULE TRAY

Pick + Hook Module, Wide

Holds:

- 10 Small Handle Picks, Hooks, Dissectors or Curettes

7 ⅝ L x 5 ⅞ W x 1 H (in)

194 L x 130 W x 25 H (mm)

IN-2002



□ EAR MODULE TRAY

Speculum Module

Holds:

- 6-12 Speculas

3 ½ L x 5 ⅞ W x 1 H (in)

89 L x 130 W x 25 H (mm)

IN-2003

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for ENT Modules

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



EAR MODULE TRAY

Micro-Suction Module

Holds:

- 12 Micro-Suction Irrigators with Luer Lock Hubs

5 1/8 L x 3 1/2 W x 1 H (in)

130 L x 89 W x 25 H (mm)

IN-2004



EAR MODULE TRAY

Suction Module

Holds:

- 10 Ferguson Style Suctions

5 7/8 L x 5 1/8 W x 1 H (in)

149 L x 130 W x 25 H (mm)

IN-2005



EAR MODULE TRAY

Alligator Module With Cover


Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

7 5/8 L x 5 1/8 W x 3 H (in)

194 L x 130 W x 76 H (mm)

IN-2006



EAR MODULE TRAY

Pick + Hook Module

Holds:

- 9 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 4 5/8 W x 1 H (in)


194 L x 117 W x 25 H (mm)

IN-2008



Preconfigured Trays for Picks + Hooks

Similar to our original module designs, these standalone trays will protect and organize pick, hook and knife instrumentation for ENT procedures. Instruments can either be positioned at a 45 degree angle for procedural presentation or secured flat for sterilization and storage.



PICK + HOOK TRAY

10 Pick + Hook Tray


Holds:

- 10 Pick + Hook Style Instruments

5 ¼ L x 8 ¾ W x 1 ½ H (in)

133 L x 222 W x 38 H (mm)

IN-4010



PICK + HOOK TRAY

20 Pick + Hook Tray


Holds:

- 20 Pick + Hook Style Instruments

11 L x 8 ¾ W x 1 ½ H (in)

279 L x 222 W x 38 H (mm)

IN-4000



PICK + HOOK TRAY

30 Pick + Hook Tray

Holds:

- 30 Pick + Hook Style Instruments

15 L x 8 ¾ W x 1 ½ H (in)

381 L x 222 W x 38 H (mm)

IN-4003

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Custom Configured Tray Example for Ear

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS



ACCESSORIES

MYRINGOTOMY SET | 26 INSTRUMENTS

BEFORE




**CUSTOM
CONFIGURED**



AFTER


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Custom Configured Tray Example for Ear

EAR TRAY | 86 INSTRUMENTS

CUSTOM
CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Custom Configured Tray Example for Ear

MAJOR EAR TRAY | 66 INSTRUMENTS

CUSTOM
CONFIGURED



Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Tray for Nose/Sinus



Basket and base can be customized to accommodate specific instrumentation.

ENDOSCOPY SINUS TRAY

Endoscopy Sinus Instrument Tray

Holds:

A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)

387 L x 244 W x 146 H (mm)

Wrapped Foot
IN-2846-R

Container Foot
IN-2846-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES





Custom Configured Tray Example for Nose/Sinus

SINUS SET | 32 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED




AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Nose/Sinus

SINUS TRAY | 34 INSTRUMENTS



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Tray for Throat

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES


Cover




Shelf



Basket



Base



ENDOSCOPY TRAY

Full Size Endoscopy Instrument Tray


Holds:

- 14 Endoscopy Style Instruments

19 5/8 L x 9 5/8 W x 7 H (in)


498 L x 244 W x 178 H (mm)

 Wrapped Foot
 IN-2847-R

 Container Foot
 IN-2847-TF


Small Parts Boxes are perfect for your small delicate instruments or accessories.


Basket and base can be customized to accommodate specific instrumentation.





Custom Configured Tray Example for Throat

THYROID SET | 25 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Tray for General

Maximize the protection and organization of your laparoscopic instrument sets with this three-level tray that is excellently constructed to fit both small and large sets of instrumentation.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

Cover




Shelf



Basket



Base



GENERAL/GYNECOLOGY TRAY

Extra Long Laparoscopic Instrument Tray

Holds:
A combination of up to 11 of the following instruments:


- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars

24 3/4 L x 9 5/8 W x 6 H (in)	629 L x 244 W x 152 H (mm)
Wrapped Foot IN-8865-R	Container Foot IN-8865-TF

Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

Basket and base can be customized to accommodate specific instrumentation.





Custom Configured Tray Example for General

TENDON STRIPPER | 13 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Custom Configured Tray Example for General

LAPAROSCOPIC | 47 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for General

BASIC UROLOGY TRAY | 35 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS



ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

CUSTOM
CONFIGURED




Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Hysteroscopy

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS



ACCESSORIES

HYSTEROSCOPE SET | 10 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Hysteroscopy

HYSTERECTOMY RESECTOSCOPE TRAY | 15 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Shelf

Base



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Trays for Kerrisons

Our Neurology and Spine Trays provide maximum accessibility and flexibility when handling and sterilizing kerrisons and forceps.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS



SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes.



KERRISON INSTRUMENT TRAY + RACK

Ten Position Instrument Tray

Holds:

- 10 Kerrisons

15 ¼ L x 9 ⅝ W x 6 H (in)

387 L x 244 W x 152 H (mm)

Wrapped Foot
IN-5012-R

Container Foot
IN-5012-TF

KERRISON INSTRUMENT RACK

Ten Position Instrument Rack

Holds:

- 10 Kerrisons

14 ¼ L x 9 ½ W x 6 H (in)

362 L x 241 W x 152 H (mm)

Wrapped Foot
N/A

Container Foot
IN-5013-TF

KERRISON INSTRUMENT RACK

Five Position Instrument Rack

Holds:


- 5 Kerrisons

14 ¼ L x 5 W x 6 H (in)

362 L x 127 W x 152 H (mm)


Wrapped Foot
N/A

Container Foot
IN-5014-TF




Preconfigured Trays for

Bi-Polar Forceps




BAYONET FORCEPS MODULE TRAY

<i>Dual Forceps Module Tray with Cover</i>	
Holds: <ul style="list-style-type: none">• 2 Forceps	
9 ¾ L x 2 ⅛ W x 1 ⅞ H (in)	248 L x 54 W x 29 H (mm)
IN-5401-02	



BAYONET FORCEPS MODULE TRAY

<i>Dual Extra Long Forceps Module Tray with Cover</i>	
Holds: <ul style="list-style-type: none">• 2 Forceps	
12 L x 2 ⅛ W x 1 ½ H (in)	305 L x 54 W x 38 H (mm)
IN-5401-03	



BAYONET FORCEPS MODULE TRAY

<i>Six Forceps Module Tray with Cover</i>	
Holds: <ul style="list-style-type: none">• 6 Forceps	
9 ¾ L x 4 ⅝ W x 1 ⅞ H (in)	248 L x 248 W x 29 H (mm)
IN-5401	

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for

Bi-Polar Forceps

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



BAYONET FORCEPS MODULE TRAY

Eight Forceps Module Tray with Cover


Holds:

- 8 Forceps

12 L x 7 ¼ W x 1 ½ H (in)

305 L x 184 W x 38 H (mm)

IN-5401-08



BAYONET FORCEPS MODULE TRAY

Twelve Forceps Module Tray with Cover


Holds:

- 12 Forceps

12 L x 9 ⅝ W x 1 ½ H (in)

305 L x 244 W x 38 H (mm)


IN-5401-12





Custom Configured Tray Example for
Neurology/Spine

NEUROLOGY TRAY | 13 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES




Custom Configured Tray Example for
Neurology/Spine

HIGH SPEED PNEUMATIC DRILL TRAY | 9 INSTRUMENTS





Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for
Neurology/Spine

DURAL REPAIR KIT TRAY | 18 INSTRUMENTS



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Ophthalmology

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

CATARACT-PHACO SET | 28 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Custom Configured Tray Example for Ophthalmology

EYE TRAY | 57 INSTRUMENTS

CUSTOM
CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Custom Configured Tray Example for Ophthalmology

CORNEA TRAY | 61 INSTRUMENTS

CUSTOM
CONFIGURED




Shelf



Base



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Ophthalmology

PHACO HAND PIECE TRAY | 23 INSTRUMENTS

CUSTOM
CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Trays for Orthopedic

While not usually delicate, orthopedic instrumentation needs to be kept sharp and straight to get the job done. Our Ortho Trays are built to safely accommodate and protect instruments that are commonly found in almost all surgical facilities.



ORTHOPEDIC TRAY

Lambotte Osteotome Tray

Holds:

- 12 Instruments—straight or curved, up to 1 ½" width (38mm)

8 ⅜ L x 9 ⅝ W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot IN-5010-R	Container Foot IN-5010-TF



CURETTE TRAY

Curette Tray

Holds:

- 14 Instruments—7-13" in length

15 ¼ L x 9 ⅝ W x 2 H (in)	387 L x 244 W x 51 H (mm)
Wrapped Foot IN-5011-R	Container Foot IN-5011-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES





Custom Configured Tray Example for Orthopedic

SMALL FRAG SET | 23 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Plastics

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

MICRO SET | 32 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Plastics

PLASTIC TRAY | 38 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Trays for Bookwalter® Retractors

This tray system is engineered to have a balanced weight distribution that meets the AAMI 25 lb. weight limit standard. We offer a two-level tray design made to accommodate the ring and blade components with an additional tray to store the table post and horizontal bar. Our two-tray system is organized for easy counting and safe for transporting. **Note:** The Mediflex Bookler Retractor also fits in this tray.



BOOKWALTER RETRACTOR TRAY

Ring + Blade Tray

Holds:

- 1 (1") Post Coupling
- 1 Balfour Retractor
- 1 Horizontal Flex Bar
- 1 Harrington Retractor
- 1 Medium Oval Ring
- 1 Segmented Ring
- 1 Gelpi Retractor
- 2 (2x6"/51mm x 152mm) Malleable Retractors
- 1 (3x6"/76mm x 152mm) Malleable Retractor
- 2 (2x3"/51mm x 76mm) Kelly Retractors
- 1 (2x4"/51mm x 110mm) Kelly Retractor
- 1 (2x5"/51mm x 127mm) Kelly Retractor
- 1 (2x6"/51mm x 152mm) Kelly Retractor
- 6 Tilt Ratchet Mechanisms
- 2 Ratchet Mechanisms

22 ½ L x 13 W x 4 ½ H (in)	572 L x 330 W x 114 H (mm)
----------------------------	----------------------------

Wrapped Foot IN-8700-R	Container Foot IN-8700-TF
---------------------------	------------------------------



BOOKWALTER RETRACTOR TRAY

Post Tray


Holds:

- Table Post and Horizontal Bar

25 L x 5 ¾ W x 3 H (in)	635 L x 137 W x 76 H (mm)
-------------------------	---------------------------

Wrapped Foot IN-8701-R	Container Foot IN-8701-TF
---------------------------	------------------------------


The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.





Custom Configured Tray Example for Retractors

GREENBERG SET | 17 INSTRUMENTS

BEFORE




CUSTOM
CONFIGURED



AFTER


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Retractors

BOOKWALTER RETRACTOR TRAY | 38 INSTRUMENTS

CUSTOM
CONFIGURED




Shelf

Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.



Custom Configured Tray Example for Retractors

CARDIAC MIS RETRACTOR HOLDER TRAY | 2 INSTRUMENTS

CUSTOM
CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

Preconfigured Trays for da Vinci SP®

The storage and reprocessing of the da Vinci SP isn't limited to a single solution. That's why we've designed multiple trays specific to the system's instrumentation to accommodate a wide variety of customer preferences.




INSTRUMENT TRAY

Single Instrument Tray

Holds:

- 1 Instrument

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8921-R	Container Foot IN-8921-TF




INSTRUMENT TRAY

Dual Instrument Tray

Holds:

- 2 Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
Wrapped Foot IN-8922-R	Container Foot IN-8922-TF



INSTRUMENT TRAY

Four Instrument Tray

Holds:


- 4 Instruments

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8923-R	Container Foot IN-8923-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci SP




INSTRUMENT TRAY

Six Instrument Tray

Holds:

- 6 Instruments

28 L x 11 W x 5 $\frac{3}{8}$ H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8925-R	Container Foot IN-8925-TF




PROCEDURE TRAY

Four Instrument Procedure Tray

Holds:

- 4 Instruments + Associated Accessories

28 L x 8 W x 5 $\frac{3}{8}$ H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8924-R	Container Foot IN-8924-TF



PROCEDURE TRAY

Six Instrument Procedure Tray

Holds:

- 6 Instruments + Associated Accessories

28 L x 11 W x 5 $\frac{3}{8}$ H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8926-R	Container Foot IN-8926-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for
da Vinci SP




ACCESSORY TRAY

Six Accessory Tray

Holds:

- 2 Cannulas
- 2 Obturators
- 2 Cables

15 ¼ L x 9 ⅝ W x 2 ⅜ H (in)	387 L x 244 W x 60 H (mm)
Wrapped Foot IN-8927-R	Container Foot IN-8927-TF



CAMERA TRAY

Single Camera Tray

Holds:

- 1 Camera

28 L x 11 W x 3 H (in)	711 L x 279 W x 76 H (mm)
Wrapped Foot IN-8920-R	Container Foot IN-8920-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Xi® Single-Site®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Single-Site instruments—from single instruments to procedural sets and their accommodating accessories—our trays are the recommended solution for instrument reprocessing.



INSTRUMENT TRAY

Single Instrument Tray

Holds:


- 1 (5mm) Xi Single-Site Instrument

27 1/8 L x 3 3/8 W x 3 1/2 H (in)

689 L x 86 W x 89 H (mm)

Wrapped Foot
IN-8911-R

Container Foot
N/A



PROCEDURE TRAY

22 Instrument Procedure Tray

Holds:


- 6 (5mm) Xi Single-Site Instruments
- 3 Rigid Obturators
- 2 Flexible Obturators
- 4 Curved Cannulas
- 3 Rigid Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

26 3/4 L x 13 W x 3 5/8 H (in)

679 L x 330 W x 92 H (mm)

Wrapped Foot
IN-8912-R

Container Foot
N/A



ACCESSORY TRAY

16 Accessory Tray

Holds:

- 3 Rigid Obturators
- 2 Flexible Obturators
- 3 Rigid Cannulas
- 4 Curved Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

15 1/4 L x 9 5/8 W x 4 1/8 H (in)

387 L x 244 W x 105 H (mm)

Wrapped Foot
IN-8913-R

Container Foot
IN-8913-TF


All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

TRANSPORT CONTAINERS

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

Preconfigured Trays for
**da Vinci Xi
Single-Site**



INSTRUMENT TRAY

Six Instrument Tray

Holds:

- 6 (5mm) Xi Single-Site Instruments

26 3⁄4 L x 10 W x 3 5⁄8 H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8914-R	Container Foot N/A

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Xi®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Instruments—Single or Dual Staplers and their accommodating accessories—our trays are the recommended solutions for instrument reprocessing.



INSTRUMENT TRAY

Single Instrument Tray

Holds:


- 1 EndoWrist Instrument

27 1/8 L x 3 1/2 W x 2 1/8 H (in)

689 L x 89 W x 54 H (mm)

Wrapped Foot
IN-8931-R

Container Foot
N/A



INSTRUMENT TRAY

Two Instrument Tray

Holds:


- 2 EndoWrist instruments

27 1/8 L x 3 3/8 W x 3 1/2 H (in)

689 L x 86 W x 89 H (mm)

Wrapped Foot
IN-8947-R

Container Foot
N/A



INSTRUMENT TRAY

Three Instrument Tray

Holds:

- 3 EndoWrist instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)

679 L x 175 W x 94 H (mm)

Wrapped Foot
IN-8948-R

Container Foot
IN-8948-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




INSTRUMENT TRAY

Four Instrument Tray

Holds:
• 4 EndoWrist Instruments

26 $\frac{3}{4}$ L x 6 $\frac{7}{8}$ W x 3 $\frac{3}{4}$ H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8950-R	Container Foot IN-8950-TF



INSTRUMENT TRAY

Six Instrument Tray

Holds:
• 6 EndoWrist Instruments

26 $\frac{3}{4}$ L x 10 W x 3 $\frac{5}{8}$ H (in)	679 L x 254 W x 92 H (mm)
Wrapped Foot IN-8936-R	Container Foot N/A



INSTRUMENT TRAY

Eight Instrument Tray


Holds:
• 8 EndoWrist Instruments

26 $\frac{3}{4}$ L x 13 W x 3 $\frac{5}{8}$ H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8938-R	Container Foot N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Xi



INSTRUMENT TRAY

Three Instrument + Accessory Tray

Holds:

- 3 EndoWrist Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8949-R	Container Foot IN-8949-TF




PROCEDURE TRAY

Four Instrument Procedure Tray

Holds:

- 4 EndoWrist Instruments + Associated Accessories

28 L x 8 W x 5 7/16 H (in)	711 L x 203 W x 138 H (mm)
Wrapped Foot IN-8915-R	Container Foot N/A



PROCEDURE TRAY

Six Instrument Procedure Tray

Holds:

- 6 EndoWrist Instruments + Associated Accessories

26 3/4 L x 13 W x 3 3/4 H (in)	679 L x 330 W x 95 H (mm)
Wrapped Foot IN-8937-R	Container Foot N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



PROCEDURE TRAY

Eight Instrument Procedure Tray

Holds:

- 8 EndoWrist Instruments + Associated Accessories

26 $\frac{3}{4}$ L x 13 W x 3 $\frac{5}{8}$ H (in)

679 L x 330 W x 92 H (mm)

 Wrapped Foot
IN-8935-R

 Container Foot
N/A

STAPLER + ACCESSORY TRAY

Single Stapler + Accessory Tray

Holds:

- 1 Stapler Instrument + Associated Accessories

27 $\frac{1}{8}$ L x 3 $\frac{3}{8}$ W x 3 $\frac{1}{2}$ H (in)

689 L x 86 W x 89 H (mm)

 Wrapped Foot
IN-8933-R

 Container Foot
N/A

STAPLER + ACCESSORY TRAY

Dual Stapler + Accessory Tray

Holds:

- 2 Stapler Instruments + Associated Accessories

26 $\frac{3}{4}$ L x 6 $\frac{7}{8}$ W x 3 $\frac{11}{16}$ H (in)

679 L x 175 W x 94 H (mm)

 Wrapped Foot
IN-8932-R

 Container Foot
N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Xi



ACCESSORY TRAY

Ten Accessory Tray

Holds:

- 6 Cannulas or Obturators
- 1 Release Tool
- 3 Cables

19 5/8 L x 9 5/8 W x 2 H (in)	498 L x 244 W x 51 H (mm)
Wrapped Foot IN-8939-R	Container Foot IN-8939-TF



SCOPE TRAY

Single Scope Tray

Holds:

- 1 (8mm) 0° or 30° Scope

24 7/8 L x 13 1/2 W x 4 3/8 H (in)	632 L x 343 W x 111 H (mm)
Wrapped Foot INS-8930-R	Container Foot N/A



CAMERA TRAY

Single Camera Tray

Holds:

- 1 Hand Held Camera

15 1/4 L x 9 5/8 W x 3 H (in)	387 L x 244 W x 76 H (mm)
Wrapped Foot INS-8934-R	Container Foot INS-8934-TF



NEED A PLACE TO SOAK?
SEE PAGE 90 FOR DETAILS

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Si[®] Single-Site[®]

Coordinate your instrument sets with multi-level trays that accommodate a full Single-Site instrument and accessory procedure set all the way down to one Single-Site instrument.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




INSTRUMENT TRAY

Single Instrument Tray

Holds:

- 1 (5mm) Si Single-Site Instrument

25 5/8 L x 3 1/8 W x 2 3/8 H (in)	651 L x 79 W x 60 H (mm)
Wrapped Foot IN-8901-R	Container Foot N/A




PROCEDURE TRAY

Seven Instrument + Accessory Tray

Holds:

- 7 (5mm) Si Single-Site Instruments
- Flexible Obturators
- 1 Single-Site Suction Irrigator
- Curved Cannulas
- Alignment Target

22 1/8 L x 13 W x 4 1/8 H (in)	562 L x 330 H x 105 H (mm)
Wrapped Foot IN-8902-G2-R	Container Foot IN-8902-G2-TF



PROCEDURE TRAY

Ten Instrument + Accessory Tray

Holds:

- 10 (5mm) Si Single-Site Instruments
- Alignment Target
- Flexible Obturators
- Curved Cannulas
- 1 Single-Site Suction Irrigator

23 7/8 L x 9 5/8 W x 6 H (in)	606 L x 244 W x 152 H (mm)
Wrapped Foot IN-8904-R	Container Foot IN-8904-TF

All da Vinci Si Single-Site Trays have green covers.

Excludes IN-8901

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Si Single-Site



ACCESSORY TRAY

Accessory Tray

Holds:

- Alignment Target
- Curved Cannulas
- Flexible Obturators

22 1/8 L x 13 W x 2 3/8 H (in)	562 L x 330 W x 60 H (mm)
Wrapped Foot IN-8903-R	Container Foot IN-8903-TF



INSTRUMENT + SUCTION TRAY

Seven Instrument + Suction Tray

Holds:

- 7 (5mm) Si Single-Site Instruments
- 1 Single-Site Suction Irrigator

22 1/2 L x 13 W x 2 3/8 H (in)	512 L x 330 W x 60 H (mm)
Wrapped Foot IN-8907-R	Container Foot IN-8907-TF

All da Vinci Si Single-Site Trays have green covers.

Excludes IN-8901

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Si EndoWrist®

No two sets are alike. That is why we have created a wide range of trays to accommodate single or multiple EndoWrist instruments and their accompanying accessories for the Si Surgical System.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



INSTRUMENT TRAY

Single Instrument Tray

Holds:


- 1 EndoWrist Instrument

25 ¼ L x 3 ¼ W x 2 ¾ H (in)

641 L x 83 W x 60 H (mm)

Wrapped Foot
IN-8885-R

Container Foot
N/A



INSTRUMENT TRAY

Five Instrument Tray

Holds:

- 5 EndoWrist Instruments

23 ¾ L x 9 ¾ W x 2 ¾ H (in)

606 L x 244 W x 60 H (mm)

Wrapped Foot
IN-8894-R

Container Foot
IN-8894-TF



INSTRUMENT TRAY

Six Instrument Tray

Holds:

- 6 EndoWrist Instruments

22 ½ L x 13 W x 2 ¾ H (in)

562 L x 330 W x 60 H (mm)

Wrapped Foot
IN-9999-160-01-R


Container Foot
IN-9999-160-01-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for

da Vinci Si EndoWrist



INSTRUMENT TRAY

Ten Instrument Tray

Holds:

- 10 EndoWrist Instruments

23 7/8 L x 9 5/8 W x 4 1/8 H (in)	606 L x 244 W x 105 H (mm)
Wrapped Foot IN-8898-01-R	Container Foot IN-8898-01-TF




INSTRUMENT TRAY

Twelve Instrument Tray

Holds:

- 12 EndoWrist Instruments

22 1/2 L x 13 W x 4 1/2 H (in)	572 L x 330 W x 114 H (mm)
Wrapped Foot IN-8892-R	Container Foot IN-8892-TF



PROCEDURE TRAY

Eight Instrument Procedure Tray

Holds:

- 8 EndoWrist Instruments
- 1 Alignment Target
- 3 Obturators
- 3 Cannulas

23 7/8 L x 9 5/8 W x 4 1/8 H (in)	606 L x 244 W x 105 H (mm)
Wrapped Foot IN-8984-01-R	Container Foot IN-8984-01-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Si EndoWrist

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



PROCEDURE TRAY

Ten Instrument Procedure Tray


Holds:

- 10 EndoWrist Instruments
- 1 Alignment Target
- 3 Obturators
- 3 Cannulas

22 ½ L x 13 W x 4 ½ H (in)

572 L x 330 W x 114 H (mm)

 Wrapped Foot
IN-8982-R

 Container Foot
IN-8982-TF


ACCESSORY TRAY

Eight Accessory Tray

Holds:

- 1 Alignment Target
- 3 Obturators
- 1-2 Cables
- 3 Cannulas

15 ¼ L x 9 ⅝ W x 3 H (in)

387 L x 244 W x 76 H (mm)

 Wrapped Foot
IN-8899-R

 Container Foot
IN-8899-TF


ACCESSORY TRAY

Seven Accessory Tray

Holds:

- 1 Alignment Target
- 3 Obturators
- 3 Cannulas

8 ⅜ L x 9 ⅝ W x 3 H (in)

213 L x 244 W x 76 H (mm)

 Wrapped Foot
IN-8897-R

 Container Foot
IN-8897-TF


da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for

da Vinci Si[®]
Stapler

Get more flexibility to meet the demands of your procedure schedule with trays designed specifically for use with the Si Series da Vinci Surgical Stapler System. Options include various feet choices for use with wrap or rigid containers.



STAPLER TRAY

Single Stapler + Accessory Tray

Holds:

- 1 EndoWrist Stapler Instrument
- 1 EndoWrist 13mm Stapler Cannula
- 1 EndoWrist 13mm Stapler Obturator
- 1 EndoWrist 13-8mm Stapler Reducer
- 1 EndoWrist Stapler Cable
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

25 L x 5 $\frac{3}{8}$ W x 3 H (in)

635 L x 137 W x 76 H (mm)

 Wrapped Foot
IN-8940-R

 Container Foot
IN-8940-TF


STAPLER TRAY

Single Stapler + Accessory Tray


Holds:

- 1 EndoWrist Stapler Instrument
- 1 EndoWrist 13mm Stapler Cannula
- 1 EndoWrist 13mm Stapler Obturator
- 1 EndoWrist 13-8mm Stapler Reducer
- 1 EndoWrist Stapler Cable
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

23 $\frac{7}{8}$ L x 9 $\frac{5}{8}$ W x 2 $\frac{3}{8}$ H (in)

606 L x 244 W x 60 H (mm)

 Wrapped Foot
IN-8944-R

 Container Foot
IN-8944-TF


STAPLER TRAY

Dual Stapler Instrument Tray

Holds:

- 1-2 EndoWrist Stapler Instruments

25 L x 5 $\frac{5}{8}$ W x 2 $\frac{3}{8}$ H (in)

635 L x 137 W x 60 H (mm)

 Wrapped Foot
IN-8945-R

 Container Foot
N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Preconfigured Trays for da Vinci Si Stapler



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS




RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



STAPLER TRAY

Dual Stapler + Accessory Tray

Holds:

- 2 EndoWrist Stapler Instruments
- 2 EndoWrist 13mm Stapler Cannulas
- 2 EndoWrist 13mm Stapler Obturators
- 2 EndoWrist 13-8mm Stapler Reducers
- 2 EndoWrist Stapler Cables
- 1 EndoWrist Stapler Release Kit
- 1 Instruction Plaque

23 7/8 L x 9 5/8 W x 2 3/8 H (in)	606 L x 244 W x 60 H (mm)
Wrapped Foot IN-8942-R	Container Foot IN-8942-TF

ACCESSORY TRAY

Single Si Stapler Motor Pack Accessory Tray

Holds:

- 1 EndoWrist Stapler 45 Motor Pack

8 1/4 L x 9 5/8 W x 3 H (in)	210 L x 244 W x 76 H (mm)
Wrapped Foot IN-8943-R	Container Foot IN-8943-TF

ACCESSORY TRAY

Stapler Accessory Tray

Holds:

- 2 EndoWrist Stapler Cannulas
- 2 EndoWrist Stapler Cannula Reducers
- 2 EndoWrist Stapler Cannula Blunt Obturators
- 2 EndoWrist Stapler Cables
- 1 EndoWrist Stapler Release Kit


25 L x 5 3/8 W x 2 3/8 H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot IN-8946-R	Container Foot IN-8946-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for da Vinci Si Scope

One of our original da Vinci products, Single and Dual Scope Trays are designed to protect the full variety of available endoscopes for the Si Surgical System and are 510(k) cleared for use in the STERRAD® 100NX® Express cycle.




SCOPE TRAY

Single Scope Tray

Holds:

- 1 (12mm) 0° or 30° Scope

22 L x 5 5/8 W x 3 H (in)	559 L x 143 W x 76 H (mm)
Wrapped Foot IN-8891-SI-12-S-R	Container Foot N/A




SCOPE TRAY

Single Scope Tray

Holds:

- 1 (8.5mm) 0° or 30° Scope

22 L x 5 5/8 W x 3 H (in)	559 L x 143 W x 76 H (mm)
Wrapped Foot IN-8891-SI-85-S-R	Container Foot N/A



SCOPE TRAY

Dual Scope Tray

Holds:

- 2 (12mm) 0° or 30° Scopes

25 L x 5 5/8 W x 3 H (in)	635 L x 143 W x 76 H (mm)
Wrapped Foot IN-9999-172-S-R	Container Foot N/A

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for da Vinci Si Scope

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



SCOPE TRAY

Dual Scope Tray


Holds:

- 2 (8.5mm) 0° or 30° Scope

25 L x 5 5/8 W x 3 H (in)

635 L x 143 W x 76 H (mm)

 Wrapped Foot
IN-9999-178-S-R

 Container Foot
N/A


SCOPE TRAY

Dual Scope Tray


Holds:

- 2 (12mm) 0° or 30° Scopes

21 1/4 L x 9 W x 3 H (in)

540 L x 229 W x 76 H (mm)

 Wrapped Foot
IN-8989-S-R

 Container Foot
IN-8989-S-TF


SCOPE TRAY

Dual Scope Tray

Holds:

- 2 (8.5mm) 0° or 30° Scopes

21 1/4 L x 9 W x 3 H (in)

540 L x 229 W x 76 H (mm)

 Wrapped Foot
IN-8988-S-R

 Container Foot
IN-8988-S-TF

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



Preconfigured Trays for

Senhance™ Surgical Robotic Trays

A sophisticated robotic system deserves a solution just as sophisticated to hold its instruments. That's why TransEnterix collaborated with us, the instrument protection tray experts. This line of trays is not only designed and constructed to fit the Senhance Surgical Robotic System instruments and adapters like a glove, but it is also validated specifically for this instrumentation.

INSTRUMENT TRAY

Instrument Tray

Holds:

- 8 Instrument Adapters
- 12 Instrument Shafts
- Cables

19 5/8 L x 9 5/8 W x 4 1/8 H (in)	498 L x 244 W x 105 H (mm)
Wrapped Foot IN-8750-R	Container Foot IN-8750-TF

ULTRASONIC INSTRUMENT TRAY

Instrument Tray

Holds:

- 1 Ultrasonic Adapter
- 1 Ultrasonic Transducer

25 L x 3 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot IN-8753-R	Container Foot IN-8753-TF

ENDOSCOPE ADAPTER TRAY

Single Instrument Tray

Holds:

- 1 Endoscope Adapter

19 5/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8751-R	Container Foot IN-8751-TF




ENDOSCOPE ADAPTER TRAY

Dual Instrument Tray


Holds:

- 2 Endoscope Adapters

19 5/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8752-R	Container Foot IN-8752-TF



Senhance is a trademark of TransEnterix Surgical, Inc. Please refer to the Senhance™ Surgical Robotic System IFU for instrument reprocessing.



Custom Configured Tray Example for Robotic Instrumentation

ARTHRO

ROBOTIC LAPAROSCOPIC SET | 9 INSTRUMENTS

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS



ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

BEFORE

CUSTOM
CONFIGURED


AFTER

da Vinci, SP, Xi, Si, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.
Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Trays for Rigid Scopes

Our rigid scope trays are universal and designed to fit scopes of any size for any procedure.




RIGID SCOPE TRAY 4-7"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 4-7" in working length

11 L x 3 ¼ W x 2 ¾ H (in)	279 L x 83 W x 60 H (mm)
Wrapped Foot IN-6011-R	Container Foot N/A



RIGID SCOPE TRAY 5-13½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 5-13½" in working length

17 ½ L x 3 ¼ W x 2 ¾ H (in)	445 L x 83 W x 60 H (mm)
Wrapped Foot IN-6012-R	Container Foot N/A



RIGID SCOPE TRAY 7-16½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 7-16½" in working length

20 ¼ L x 3 ¼ W x 2 ¾ H (in)	515 L x 83 W x 60 H (mm)
Wrapped Foot IN-6013-R	Container Foot N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Trays for
Rigid Scopes

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



RIGID SCOPE TRAY 9 1/2-21 1/2"

Holds: <ul style="list-style-type: none">• Up to 3 (2-10mm) Scopes• Accommodates scopes from 9 1/2-21 1/2" in working length	
25 1/4 L x 3 1/4 W x 2 3/8 H (in)	641 L x 83 W x 60 H (mm)
Wrapped Foot IN-6014-R	Container Foot N/A




RIGID SCOPE TRAY 11-26"

Holds: <ul style="list-style-type: none">• Up to 3 (2-10mm) Scopes• Accommodates scopes from 11-26" in working length	
30 1/8 L x 3 1/4 W x 2 3/8 H (in)	765 L x 83 W x 60 H (mm)
Wrapped Foot IN-6015-R	Container Foot N/A



Preconfigured Trays for Scope + Cable


Our scope + cable trays allow for faster loading before sterilization while also protecting the delicate fiber optics within.



SCOPE + CABLE TRAY 5-6"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 5-6" in working length


11 L x 6 3/8 W x 2 3/8 H (in)	279 L x 162 W x 60 H (mm)
Wrapped Foot IN-6031-R	Container Foot N/A



SCOPE + CABLE TRAY 5-13½"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 5-13½" in working length

20 1/2 L x 5 3/8 W x 2 3/8 H (in)	521 L x 143 W x 60 H (mm)
Wrapped Foot IN-6032-R	Container Foot N/A



SCOPE + CABLE TRAY 6½-17"

- Holds:
- Up to 3 (2-10mm) Scopes
 - 1 Cable
 - Accommodates scopes from 6½-17" in working length

25 L x 5 3/8 W x 2 3/8 H (in)	635 L x 137 W x 60 H (mm)
Wrapped Foot IN-6033-R	Container Foot N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





Custom Configured Tray Example for Scope + Cable

RESECTOSCOPE | 11 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Preconfigured Trays for Flexible Scopes

Our durable flexible scope trays are sure to protect your delicate and expensive scopes. Our Small Parts Boxes are perfect to keep all your accessories together and protected.



FLEXIBLE SCOPE TRAY

Holds:


- 1 Flexible Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot
IN-6021-R

Container Foot
N/A



FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Scope with Attached Light Source
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot
IN-6022-R

Container Foot
N/A



FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Video Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot
IN-6023-R

Container Foot
N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Preconfigured Tray for
Flexible Scopes



FLEXIBLE SCOPE TRAY

Holds: <ul style="list-style-type: none">• 1 Flexible Bronchoscope• Accommodates scopes from 12-27" in working length	
19 5⁄8 L x 9 5⁄8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6024-R	Container Foot N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES




Custom Configured Tray Example for Flexible Scopes

FLEXIBLE URETEROSCOPE | 2 INSTRUMENTS

BEFORE



CUSTOM
CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Custom Configured Tray Example for Flexible Scopes

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

VIDEO LARYNGOSCOPE | 1 INSTRUMENT

BEFORE




CUSTOM
CONFIGURED



AFTER


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Flexible Scopes

COLONOSCOPY TRAY | 1 INSTRUMENT

CUSTOM
CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




Custom Configured Tray Example for Flexible Scopes

FLEXIBLE SCOPE TRAY | 3 INSTRUMENTS

CUSTOM
CONFIGURED




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



Custom Configured Tray Example for Flexible Scopes

FLEXIBLE SCOPE TRAY | 1 INSTRUMENT

CUSTOM
CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

XL Transport Containers

Transporting da Vinci SP, Xi scopes, and EndoWrist instruments appropriately and effectively is a difficult task. Use the InstruSafe XL Transport Container to get the job done right! This durable container is the perfect size, enabling you to soak and protect your instruments while they travel from the O.R. to sterile processing. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



IN-PIC-20-D



IN-PIC-20-D-LT



Single Endoscope Configuration

Multiple Instrument Configuration

XL TRANSPORT CONTAINER

XL Transport Container + Optional Drain

Holds:

- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 ¼ W x 6 H (in)
864 L x 438 W x 152 H (mm)

With Drain
IN-PIC-20-D

Without Drain
IN-PIC-20

XL TRANSPORT CONTAINER SYSTEM

XL Transport Container with Lifting Tray + Optional Drain

Holds:

- Lifting Tray
- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 ¼ W x 6 H (in)
864 L x 438 W x 152 H (mm)

With Drain
IN-PIC-20-D-LT

Without Drain
IN-PIC-20-LT

LIFTING TRAY

Single Endoscope Configuration

Holds:

- 1 da Vinci Xi or SP Endoscope + Associated Accessories

Multiple Instrument Configuration

Holds:

- Up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)

699 L x 343 W x 98 H (mm)

IN-PIC-34

XL Transport Container

Special Features

- Highly durable material designed to withstand everyday use
- Cart wash safe—maximum temperature of 230°F
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Made in the USA

Inner gasket and strong silicone latches reduce spills

Built-in biohazard symbol on cover



Dedicated dry erase "soak start time" space on cover




Water volume feature to accurately measure detergent solutions

Quick drain feature—no need to tip container to drain



LIFTING TRAY

Single Endoscope Configuration
Holds: 1 da Vinci Xi or SP Endoscope




Multiple Instrument Configuration
Holds: up to 8 da Vinci Xi or SP EndoWrist Instruments

Accessory Basket



Large Handles



Loading Identifiers

CONFIGURATION EXAMPLES

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Heavy Duty Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. They are similar to our other transport containers but have a more robust construction and feature stainless steel hinges and latches that can be easily undone with just one hand.



HEAVY DUTY TRANSPORT CONTAINER

<i>Short Container</i>	
20 ½ L x 12 ⅝ W x 3 ¾ H (in)	521 L x 321 W x 95 H (mm)
IN-PIC-15	



HEAVY DUTY TRANSPORT CONTAINER

<i>Tall Container</i>	
20 L x 12 ¼ W x 6 ⅞ H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	

SPECIAL FEATURES

- Autoclavable
- Made of PEI material (high-temperature plastic)
- High durability against automated washing
- Transparent amber color for easy instrument identification
- Easy locking, removable cover with removable silicone gasket that can be replaced
- Lightweight
- Combine with Instru-Zyme for optimal pre-cleaning

Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. Also, you can easily identify content through the transparent amber color. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



TRANSPORT CONTAINER

Small, Short Container

13 $\frac{9}{16}$ L x 15 $\frac{5}{8}$ W x 4 $\frac{5}{16}$ H (in)	345 L x 398 W x 110 H (mm)
--	----------------------------

IN-PIC-14



TRANSPORT CONTAINER

Large, Short Container

20 $\frac{5}{8}$ L x 15 $\frac{5}{8}$ W x 5 $\frac{5}{16}$ H (in)	524 L x 397 W x 135 H (mm)
---	----------------------------

IN-PIC-12



TRANSPORT CONTAINER

Small, Tall Container

13 $\frac{9}{16}$ L x 15 $\frac{5}{8}$ W x 7 $\frac{1}{2}$ H (in)	345 L x 398 W x 190 H (mm)
---	----------------------------

IN-PIC-13



TRANSPORT CONTAINER

Large, Tall Container

20 $\frac{5}{8}$ L x 15 $\frac{5}{8}$ W x 7 $\frac{1}{2}$ H (in)	524 L x 397 W x 191 H (mm)
--	----------------------------

IN-PIC-11

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Silicone Small Parts Boxes

Don't sweat the small stuff. Made of silicone, these boxes are perfect for keeping extra delicate small parts and accessories organized and protected.

ARTHO

BIARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



SILICONE SMALL PARTS BOX


<i>3x3" Small Parts Box</i>	
3 L x 3 W x ¾ H (in)	76 L x 76 W x 19 H (mm)
IN-2652	



SILICONE SMALL PARTS BOX

<i>3x6" Small Parts Box</i>	
6 L x 3 W x ¾ H (in)	152 L x 76 W x 19 H (mm)
IN-2653	

- Made entirely of medical grade silicone
- Maximum protection for small parts
- Tethered cover to prevent accidental loss
- Can be attached to trays or placed on posts for easy removal
- Textured internal surface and proper ventilation for efficient sterilant penetration




- Small Parts Box is shown within tray. Tray is shown without cover for illustrative purposes only.

Plaque Options

A variety of plaque colors are available to help identify specific sets by procedure or service type. Affix them to our trays and modules using our standard screws. Common identification uses are for bar code placement, tray name labels, or other custom information. Laser engraving of plaques is provided at an additional cost.

Our in-house laser engraving capabilities allow us to personalize trays with custom identification and/or specifications, such as interior instrument descriptions, illustrations, exterior procedure types, specialty, or facility information.



PLAQUES

Plaque Assembly for Trays

1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-01	SILVER
1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-05	BLACK
1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-15	BLACK
Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.		
4 3/8 L x 2 W (in) 111 L x 51 W (mm)	IN-2600-XX	●●●
2 L x 2 1/8 W (in) 51 L x 56 W (mm)	IN-2300-XX	●●●

COLOR OPTIONS

Steam Sterilization

These anodized aluminum plaques for use with steam only.

Low-Temp Sterilization


Plaques have a coating that specially withstands low-temperature sterilization.

01	SILVER		N/A	
02	RED		RED	12
03	GREEN		GREEN	13
04	BLUE		BLUE	14
05	BLACK		BLACK	15
06	GOLD		GOLD	16
07	PURPLE		PURPLE	17
08	ORANGE		N/A	
10	MINT GREEN		N/A	



Finger Mats

These textured silicone instrument cushion mats are available in a variety of sizes and textures to fit all of your instrument protection needs.




FINE FINGER MAT

Full Sheet

19 L x 13 W (in)

483 L x 330 W (mm)

IN-2616




FINE FINGER MAT

Half Sheet

9 ½ L x 13 W (in)

241 L x 330 W (mm)

IN-2617



FINE FINGER MAT

Quarter Sheet

4 ¾ L x 13 W (in)

121 L x 330 W (mm)

IN-2618

Guardian™ Tray

Simple, small instrument sets now have a new home. Our silicone finger mats with sidewalls provide a safe and pliable environment for them to reside in.



GUARDIAN TRAY

Anodized Aluminum Tray with Finger Mat

10 L x 6 W x 1 ¼ H (in)

254 L x 152 W x 32 H (mm)

IN-1310

- Accommodate a wide variety of small instrumentation
- Feature built-in and vented silicone finger mats
- Reduce instrument repair and replacement costs
- Allow maximum sterilant penetration
- Prevent moisture accumulation

Instru-Zyme®

This powerful enzyme and detergent formulation is designed to assist in the removal of blood, tissue, mucus and other proteinaceous materials. Use this product with our transport containers for a simple solution to water and detergent.



INSTRU-ZYME

Properties:

- Specific Gravity: 1.200
- pH: 7.8
- Tropical Fragrance
- Color: Light Blue
- Form: Liquid Gel

6 Bottles/Case

IN-1032

- Will not stain instruments or equipment made of: metal, plastic, glass or rubber
- Non-corrosive
- Biodegradable
- Sustains moisture during transport
- Provides a ready-to-use, neutral pH enzyme solution
- Reduces heavy instrument transport weight as compared to liquid-filled trays
- Minimizes potential exposure due to splashing
- Spray applicator for even coverage of soiled areas
- No rinsing necessary before automated processing

Rounders® Point Guards

Wrap your trays with confidence. Point Guards are engineered to sit securely on almost any tray corner. Prevent holes and tears in your sterilization wrap caused by surgical trays that are heavy, have sharp corners or are improperly handled.



ROUNDERS POINT GUARDS


80/Bag

IN-7330

- Air cushion tray corner protection
- Help avoid compromised sterility in the O.R.
- Can be used in steam, ethylene oxide and low temperature plasma sterilization
- Reusable and economical
- Patent #9119611
- Non-sterile

MADE IN THE USA





InstruSafe

Sterile Processing

Summit Medical LLC**An Innovia Medical Company**

815 Vikings Parkway

Suite 100

St. Paul, MN 55121 | USA

1.888.229.2875

www.innoviamedical.com/instrusafe

InstruSafe Instrument Protection Trays are proudly made in St. Paul, Minnesota, USA.

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and CMDCAS certified.

InstruSafe Trays are 510(k) cleared for multiple sterilization cycles to give our customers more options.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

da Vinci, SP, S, Si, Xi, Single-Site and EndoWrist are registered trademarks of Intuitive Surgical, Inc. InstruSafe, Instru-Zyme, Hold-Its, Hold-Downs, Rounders Point Guards, and Guardian are trademarks and registered trademarks of Summit Medical.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

All brands are trademarks or registered trademarks of their respective owners.

©2020 Innovia Medical. All rights reserved. M-5 | RevB | DEC2020