



By Summit Medical LLC



Elevating Care. Improving Outcomes.

Innovia Medical is a worldwide manufacturer of medical devices across the specialty surgery spectrum including:

- Aesthetic & Reconstructive
- Ear, Nose & Throat
- Dermatology
- General & Vascular
- Gynecology
- Neurosurgical
- Ophthalmic
- Orthopedic
- Personal Protective Equipment
- Sterile Processing

Innovia Medical exports to over 80 countries worldwide. Every Innovia Medical company is passionately committed to partnering with medical professionals to help elevate the delivery of care and improve patient outcomes.

Summit Medical

Summit Medical works closely with healthcare professionals to develop unique, patented products that improve surgical outcomes, make healthcare more efficient, and meet industry demands through all stages of design, engineering, and manufacturing.

Network Medical

Network Medical has 20+ years of experience designing and manufacturing medical products for the ophthalmic and ENT sectors. They work with industry leaders to develop quality products sold in over 72 countries.

Eagle Labs

Eagle Labs continually strives to offer high quality and cost-effective single-use medical devices to the Ophthalmic Medical Community. Their products are manufactured in the USA and delivered to customers through a worldwide network of distributors.

DTR Medical

DTR Medical specializes in the design and manufacture of sterile single-use surgical instruments that deliver both patient and clinical value. Exporting to over 30 countries, they built their award-winning reputation with high-quality, innovative products and exceptional customer service.

Learn more at www.innoviamedical.com.

4

TISSU-TRANS

Filtron
Filtron Stand
AquaVage

8

ASPIRATORS

Waste Canisters
Aspirator Filters

9

TUBING

Infiltration
Aspiration

11

FAT GRAFTING GUIDE

12

CANNULAS

Cannula Cleaning
Cannula
Accessories
Infusion Handle
Infusion Cannula
Disposable Cannulas

21

WOUND DRAINAGE

Kits

22

COMPRESSION GARMENTS

Breasts
Face
Body

28

RHINOPLASTY

Nasal Packing
Epistaxis Balloons
Internal Nasal Splints
External Nasal Splints
Cross Reference Guide

SUBMIT PURCHASE ORDERS TO:



Summit Medical LLC
815 Vikings Parkway, Suite 100
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939
F: 1-888-229-1941 | +1 651-789-3979
www.summitmedicalusa.com

Our products are available through distributors and Group Purchasing Organization contracts. Medline, Henry Schein, McKesson, Cardinal Health and Owens & Minor are some of the participating companies. Please call customer service with any questions.

FILTRON®


The FILTRON is a closed, sterile filtering system that filters and collects viable fat, removing waste, non-viable cells and excess fluid during fat harvest. The FILTRON does not change the normal flow of the physician’s technique and does not require decant time. The Filtron unit is placed in-line between the harvest cannula and the waste canister/suction source. Suction pressure is reduced to 15-18’ Hg and fat is drawn into the mesh filter within the canister. Waste (blood, epinephrine, xylocaine) and non-viable cells are pulled through the filter and are suctioned out to the waste canister. When the desired volume is collected, the transfer tube is attached and the Filtron unit is inverted. The quality fat is then drawn directly into a luer lock syringe and is ready to reinject.

Each unit includes six feet of harvest tubing and a transfer tube with luer lock tip.

PART NO.	VOLUME	QTY/BOX
3-TT-FILTRON 100	100cc	1
3-TT-FILTRON 250	250cc	1
3-TT-FILTRON 500	500cc	1
3-TT-FILTRON 1000	1000cc	1
3-TT-FILTRON 2000	2000cc	1



AVAILABLE IN FIVE VOLUMES
VOLUME IS TOTAL FAT READY FOR RE-INJECTION WITH THE WASTE FRACTION REMOVED



Filtron Stand

Canister stand has four universal brackets and has five legs for stability. FILTRON devices can be hung upright for harvest or upside down for fat removal. The extended height FILTRON stand raises the sterile field.

Tip: Hang a sterile drape or mayo stand cover over the stand prior to hanging a FILTRON to maintain sterility.

PART NO.	DESCRIPTION	HEIGHT	QTY/BOX
3-T-STAND-DLX	Standard Height	34"	1
3-T-STAND-TALL	Extended Height	52"	1



Filtron IV Pole Bracket

This reusable, autoclavable universal bracket with clamp fits any IV pole. Devices can be hung upright for harvest or upside down for fat removal.

PART NO.	QTY/BOX
3-TT-BRACK-HK	1



AQUAVAGE®

The AQUAVAGE is a closed, sterile decanting system that separates viable fat from waste, non-viable cells and excess fluid during fat harvest. The AQUAVAGE does not change the normal flow of the physician’s technique. The Aquavage unit is placed in-line between the harvest cannula and the waste canister/suction source. Suction pressure is reduced to 15-18’ Hg. Lactated Ringers solution is placed in the bottom of the Aquavage canister to enable separation of the fat as it flows into the canister. When the fat hits the Lactated Ringers, it separates from the waste portion and floats to the top. Once gravity decantation process is complete, fat is removed via the extraction tube with a Luer Lock syringe and is ready to reinject.

Each unit includes six feet of harvest tubing and a transfer tube with luer lock tip.

PART NO.	VOLUME	QTY/BOX
AV1200D	1200cc	1
AV2000D	2000cc	1

COMING
2022



- Closed, Single-Use System
- Sterile, Disposable
- Separates Tissue and Fluids by Decantation
- Compatible with Most Aspirators and Instruments
- Eliminates Centrifuge
- Fat Immediately Ready for Reinjection

AVAILABLE
IN TWO
VOLUMES



1200cc

2000cc

AquaVage Stand

The AquaVage Stand is autoclavable and made completely of stainless steel.

PART NO.	DIMENSIONS (IN)	QTY/BOX
AVSTAND	7 ¾ W x 8 ⅝ L x ½ H	1



WASTE CANISTERS

Implosion resistant waste collection canisters are disposable, liner-free and have sealable ports. Canisters offer easy view of aspirate and have a high pathogen filter.

PART NO.	VOLUME	QTY/BOX
LS7020*	2000cc	10
LS7030	3000cc	10

**most popular*



Waste Canister Holder for IV Pole

- Autoclavable holder connects the LS7020 and LS7030 Waste Canisters to any IV Pole.
- Allows for varying placement of the waste containers to maximize surgical space
 - Reusable

PART NO.	QTY/BOX
TT-HOLDER NR	1

ASPIRATOR FILTERS



Replacement Absolute Biofilter

Replacement Absolute Biofilter prevents aerosolized pathogens from your exhaust from seeping into the air. Biofilters have a 15 hour replacement interval and each is 99.999% DOP rated for 0.02 micron particles.

PART NO.	QTY/BOX
FL3000-05	5
FL3000M-05*	5

**works with Microaire aspirator*

INFILTRATION

Infiltration Tubing with Pump Segment

Sterile auto-fuse infiltration tubing has a luer connector and a spike on opposite ends. Disposable tubing is compatible with all peristaltic head infiltration pumps.




PART NO.	DESCRIPTION	OUTER DIAMETER	QTY/BOX
AFTD	Single Spike	3/8" (.375")	10
AFTDY	Double Spike	3/8" (.375")	10

Soft Tube

Sterile tubing is soft and pliable for use with BPI and standard infiltration pumps. Disposable tubing is spiked on one end.

PART NO.	LENGTH	OUTER DIAMETER	QTY/BOX
SOFTTUBE	13 ft	.315"	10



ASPIRATION

High Pressure Tubing for Liposuction

Sterile PVC tubing is latex-free and has thick walls for high pressure suction. Disposable tubing is compatible with all available aspirators and has excellent drapability.

PART NO.	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY/BOX
LS7010	3/8" (.375")	5/8" (.625")	8 ft	10
LS7012	3/8" (.375")	3/8" (.625")	10 ft	10



PAL Microaire® Tapered Tubing

Sterile tapered tubing is individually packaged, 12 feet in length, and disposable. Pair tubing with the luer adapter to use with LL cannula and PAL hand piece.

PART NO.	DESCRIPTION	LENGTH	QTY/BOX
PAL TUBING	Tubing	12 ft	10
3-PAL-CPLR	Luer Adapter	N/A	1

CANNULA CLEANING BRUSHES

Non-sterile stainless steel rod with medical grade nylon bristles.

PART NO.	SIZE	QTY/BOX
3-C-BR2	2mm x 42cm (0.07" x 16.5")	1
3-C-BR3	3mm x 48cm (0.11" x 19")	1
3-C-BR4	4mm x 48cm (0.15" x 19")	1
3-C-BR5	5mm x 48cm (0.2" x 19")	1
3-C-BR6	6mm x 48cm (0.23" x 19")	1
3-C-BR8	8mm x 48cm (0.31" x 19")	1
3-C-BR10	10mm x 48cm (0.39" x 19")	1



FAT GRAFTING GUIDE

			<100cc	Eyelids	Lips	Face	Body (Small)	Body (Large)	Breast	Buttocks
Infusion Cannulas	3-C-BAN-16x15			•	•	•				
	3-C-GAN-16x15									
	3-C-BAN-14x30						•	•	•	•
	3-C-GAN-14x30									
Harvest Cannulas	3-C-24H-3x15			•	•	•				
	3-H-BA-24H-3x15 (handle)									
	3-H-BA-4H-3x30 (handle)						•	•	•	•
Collection Units	3-TT-FILTRON 100		100cc		•	•	•			
	3-TT-FILTRON 250		250cc				•			
	3-TT-FILTRON 500		500cc				•	•	•	
	3-TT-FILTRON 1000		1000cc					•	•	
	AV1200D		1200cc					•	•	
	AV2000D		2000cc					•	•	•
	3-TT-FILTRON 2000		2000cc					•	•	•
Injection Cannulas	3-C-1H2-18x8			•	•	•				
	3-C-1H2-16x8			•			•			
	3-C-1H2-14x8			•			•		•	
	3-C-1H2-3x15							•	•	•

This is a guide only. Treatment plan is determined by healthcare professional.

HOW TO ORDER

A multitude of hole patterns, diameters and length combinations are available. Use the example below to help build your cannula. Replace **XX** with desired hole pattern, diameter, and length.

Diameter

18gauge16gauge14gauge2mm3mm4mm5mm6mm

Length

8cm10cm12cm15cm20cm25cm30cm35cm40cm45cm

LUER LOCK

3 - C - XX - XX - XX
Part NumberHole PatternDiameterLength

HANDLE

3 - H - BA - XX - XX - XX
Part NumberHole PatternDiameterLength

Blue Aluminum Handle with Cannula

3-H-BA-*Hole Pattern-DiameterxLength*
1" Autoclavable Blue Aluminum Handle with Cannula of choice.

Luer Lock Handle

3-H-BAL
1" Autoclavable Blue Aluminum Handle with Luer Lock Hub.


Luer Lock Cannulas

3-C-*Hole Pattern-DiameterxLength*
Autoclavable stainless steel Luer Lock Cannula.



Toomey 60cc Cannulas

3-C-T60-*Hole Pattern-DiameterxLength*
Autoclavable Stainless Steel 60cc hub fits 60cc Toomey Syringes and Handles.



Toomey Handle

3-H-TH60
Autoclavable Black Delrin 60cc Toomey Handle.

CANNULA HOLE PATTERNS

INFILTRATION



BAN 6 Hole Body Spray



GAN 6 Hole Garden Spray



MH 6 Multi Hole Spray

HARVEST



4H Large Volume 4 Hole Harvester



24H Small Volume 24 Hole Harvester


REINJECTION



1H2 One Hole Type II (Coleman)



1H3 One Hole Type III (Spoon Tip)



1HV One Hole V Type

TRANSFER



OP Open-Ended

STANDARD



MC Mercedes



TP Triport



SP Spatula



2SP 2 Port Spatula



LV Las Vegas



MAN Bucket




MS Modern Spiral



TV Toledo V



DS Disrupter



BB Becker Basket



BTD Becker Tear Drop

Contact customerservice@innoviamedical.com for custom cannulas. Please allow 5-7 business days for cannula orders. No returns on cannula orders.

ACCESSORIES



Disposable Luer to Luer Transfer

Sterile female Luer Lock to Luer Lock transfer system is made of Polypropylene and has large finger grips to protect against contamination during tissue transfer or emulsification.

PART NO.	QTY/BOX
3-S-13901	50

Reusable Luer to Luer Transfer

Reusable stainless steal female Luer Lock to Luer Lock transfer system. Used to transfer tissue or emulsify between syringes.

PART NO.	SIZE	QTY/BOX
3-T-TLL2	2mm	1
3-T-TLL3	3mm	1
3-T-TLL4	4mm	1

Reusable Toomey to Luer Lock Syringe Adapter

Reusable stainless steel 60cc Toomey Monoject syringe hub to Luer Lock Syringe transfers tissue between syringes and is autoclavable.


PART NO.	QTY/BOX
3-S-T60-LL-SS	1

ACCESSORIES

Monoject™ Syringes

Monoject™ Syringes are sterile, individually packaged, single-use syringes without a needle. They have a polypropylene barrel and plunger rod with graduated barrel markings.


PART NO.	DESCRIPTION	VOLUME	QTY/BOX
3-S-T60-20	Toomey Tip	60cc	20



Syringe Racks

Syringe Racks are made of a unique polypropylene that is autoclavable up to 121°C (250°F) and maintain stability while submerged in water. Racks are easy to assemble and stackable.

Dimension of assembled:
24.6 L x 10.4 W x 6.4 H (mm)
9.7 L x 4.09 W x 2.52 H (in)




PART NO.	SIZE	HOLE	HOLE DIAMETER	FOR SYRINGE	QTY/BOX
3-S-R5	5cc	90	13mm	5cc, 3cc & 1cc BD	1
3-S-R10	10cc	60	16mm	10cc BD	1
3-S-R30	30cc	24	25mm	30cc BD, 20cc Monoject	1
3-S-R60	60cc	21	30mm	60cc BD, 35cc & 60cc Monoject Syringe	1

ACCESSORIES

Luer Lock Syringe Tip Caps

Sterile red polypropylene luer lock caps for standard luer lock syringes are individually packaged and disposable

PART NO.	QTY/BOX
3-S-65001	100



Horseshoe Locks

Locks are used to maintain negative pressure in syringes during manual tissue aspiration. Horseshoe locks are aluminum and autoclavable.

PART NO.	DESCRIPTION	FOR SYRINGE	QTY/BOX
3-S-HL10	10cc	10cc BD Syringe	1
3-S-HL10M	10cc	10cc Monoject Syringe	1
3-S-HL20	20cc	20cc BD and Terumo	1
3-S-HL30	30cc	30cc BD Syringe	1
3-S-HL35	35cc	35cc Monoject Syringe	1
3-S-HL60	60cc	60cc BD or Monoject Syringe	1

Clicker Locks

Locks are used to maintain negative pressure in syringes during manual tissue aspiration. Clicker locks are aluminum with stainless steel springs and are autoclavable.


PART NO.	DESCRIPTION	FOR SYRINGE	QTY/BOX
3-S-CL10	10cc	10cc BD Syringe	1
3-S-CL10-M	10cc	10cc Monoject Syringe	1
3-S-CL20	20cc	20cc BD	1
3-S-CL30	30cc	300cc BD	1
3-S-CL35	35cc	35cc BD and Monoject Syringe	1
3-S-CL60	60cc	60cc BD and Monoject Syringe	1

ACCESSORIES

MicroAire Pal Luer Adapter

The MicroAire Pal Luer Adapter connects a PAL hand piece to a luer lock aspiration cannula. It is superiorly constructed with both aluminum and stainless steel and is autoclavable. The luer lock connection has a larger 3.8mm orifice which is the same as a 4mm diameter cannula. From the luer connection, through the body of the adapter to the aspiration tubing connection, is 5mm.

PART NO.	LENGTH	QTY/BOX
3-PAL-CPLR	5mm	1



Syringe Power Adaptor 60cc

Syringe Power adaptor connects a 60cc syringe to tubing.

PART NO.	QTY/BOX
SPA60CC	1



INFUSION HANDLES


Infiltration Handles

Autoclavable blue aluminum stainless steel handles with Luer Lock fittings on both ends.

PART NO.	SIZE
3-H-IHSL	Large 1" diameter
3-H-IHSS	Small ¾" diameter



3-H-IHSL



3-H-IHSS

INFUSION CANNULAS



GAN 6 Hole Garden Spray



BAN 6 Hole Body Spray



MH 6 Multi Hole Spray

Tumescent Infiltration Cannulas


Autoclavable stainless steel Luer Lock cannulas for infusion and infiltration.

PART NO.	DESCRIPTION
3-C-GAN	6 Hole Garden Spray
3-C-BAN	6 Hole Body Spiral
3-C-MH	6 Multi Hole Spray

Infiltration IV Spike

Connects an IV bag to infiltration tubing through spike connection.

PART NO.	QTY/BOX
IVSPIKE	1




DISPOSABLE CANNULA

Luer Lock Cannulas

Disposable sterile hub with standard Luer Lock connection and Lipocyte friendly hole patterns used to remove fat from collection devices. 10/box.

One Hole Type II	Description
CF301590	1 Hole II Injector, 3mmx15cm
CF201590	1 Hole II Injector, 2mmx15mm
CF10790	1 Hole II Injector, 1mmx7mm
CF16890	1 Hole II Injector, 16gaugex8cm
CF12790	1 Hole II Injector, 18gaugex7cm
Open-Ended	
3COE3300	OP, 3mmx30cm
Multi-Hole	
CP301524H	24 Hole, 3mmx15cm
CP141508H	8 Hole, 2.1mmx15cm
CP141506HS	6 Hole, Spiral, 2.1mmx15cm
CP30154H	4 Hole, 3mmx15cm



Handle and Cannula

Lightweight sterile white handle with thumb groove. 6/box.

3 Hole	Description
LC32303	3 Hole, 3mmx23cm
LC42803	3 Hole, 4mmx28cm
LC52803	3 Hole, 5mmx28cm
24 Hole	Description
LC328024H	24 Hole, 3mmx28cm
LC428024H	24 Hole, 4mmx28cm



3 Hole







24 Hole

DISPOSABLE CANNULA

Micro Reinjection Cannula

Disposable sterile Luer Lock hub used for fat or filler reinjection with blunt atraumatic tip and one side port hole pattern.


PART NO.	DESCRIPTION	
MI18110	18gaugex11cm	
MI2290	22gaugex9cm	
MI2270	22gaugex7cm	
MI2250	22gaugex5cm	
MI2370	23gaugex7cm	
MI2350	23gaugex5cm	
MI2550	25gaugex5cm	
MI2638	26gaugex3.8cm	
MI2625	26gaugex2.5cm	
MI2738	27gaugex3.8cm	



WOUND DRAINAGE

SeromaCath® Wound Drainage Kits

Sterile wound drainage systems for post-operative seromas following mastectomy, axillary dissection, groin dissection, latissimus dorsi flap, TRAM flap, liposuction, abdominoplasty, breast implants and other procedures. Kit comes with connecting drainage tube for a bulb suction reservoir; bulb sold separately.



PART NO.	KIT CONTENTS	QTY/BOX
GR-1275	2.75" (7cm) Blunt for Seromas surrounding Breast implants 2" (5cm) 14gauge Teflon catheter/stylet needle/plug with two side holes 2.75" (7cm) 17gauge blunt metal catheter with two side holes Connecting drainage tube for a bulb suction reservoir	5 kits
GR-1035	3.5" (8.9cm) for Deep Seromas common with abdominoplasty 3.5" (8.9cm) 14gauge Teflon catheter/stylet needle/plug with six side holes Connecting drainage tube for a bulb suction reservoir	5 kits
GR-1002	2" (5cm) for Common Seromas 2" (5cm) 14gauge Teflon catheter/stylet needle/plug with six side holes Connecting drainage tube for a bulb suction reservoir	5 kits

WOUND DRAINAGE
BULBS ARE NOT
INCLUDED IN THE
DRAINAGE KITS

SeromaCath® Sapphire Wound Drainage Bulb

Sterile bulbs with anti-reflux valve and volume markers for ideal fluid management used in general, orthopedic, urologic and neurosurgery. High silicone content gives the bulb the highest suction of any commercially available 100cc silicone reservoir offering 55% more suction than the popular "JP Bulb." Stainless steel clip allows for convenient and discreet wear.

PART NO.	SIZE	QTY/BOX
GR-2001	100ml	1 bulb



VERONIQUE SIZE CHART

Female Size Chart

Our compression garments are latex-free and come in a variety of sizes and colors for both males and females. Garment stock is limited to availability. We recommend placing your order at least one week prior to the date the garment is needed. We do not offer returns on garment orders.

SIZE	WAIST	HIP	BUST	BREAST	UNDER BICEP CIRCUMFERENCE
XS	21-31"	31-33"	27-29"	25-26"	9.50-10.25"
Small	24-26"	34-36"	30-32"	27-30"	10.50-11.25"
Medium	27-29"	37-39"	33-36"	31-34"	11.50-12.25"
Large	30-32"	40-42"	37-40"	35-38"	12.50-13.25"
XL	33-35"	43-45"	41-44"	39-42"	13.50-14.25"
2XL	36-39"	45-48"	45-48"	43-46"	14.50-15.25"
3XL	40-43"	49-52"	49-52"	47-49"	15.75-16.75"
4XL	44-47"	53-56"	53-56"	50-52"	17.00-18.00"

NO RETURNS ON GARMENT ORDERS

Garments are not guaranteed in stock. Please order one week prior to procedure.
All Veronique Garments are LATEX-FREE.



Male Size Chart

SIZE	WAIST	HIP	CHEST
Small	31-33"	36-38"	35-37"
Medium	34-36"	39-41"	38-41"
Large	37-39"	42-44"	42-45"
XL	40-43"	45-47"	46-49"
2XL	44-48"	48-51"	50-53"
3XL	49-53"	52-55"	54-58"
4XL	54-56"	56-59"	59-62"

NO RETURNS ON GARMENT ORDERS

Garments are not guaranteed in stock. Please order one week prior to procedure.
All Veronique Garments are LATEX-FREE.



BREASTS

Breast/Body Wraps

All of our breast wraps are made of nylon and spandex, offer adjustable downward compression, and have adjustable front hooks and loop clasps. They are one size fits all.



V-510-3
3" Breast Wrap



V-510-3U
3" Breast Wrap with Cotton Lining



V-510-6
6" Breast Wrap



V-510-9
9" Breast Wrap

Specialty Breast Wraps



V-510-J
3" Breast Wrap with Lace

Breast wrap with lace is made of nylon and spandex, offers adjustable downward compression, and has adjustable front hooks and loop clasps. They are one size fits all.



V-510-X
Versatile Stabilizing Breast Wrap

Versatile stabilizing breast wrap is made of nylon and spandex, offers customized band positioning with, and has adjustable front hooks and loop clasps. They are one size fits all.

BREASTS

Breast/Body Wraps

Panel binders have adjustable front hook and loop clasps and reinforced back stays for a strong and durable stretch. They are available in XS-4XL.



V-350-9
9" Double Panel Binder



V-350-6
6" Single Panel Binder



Support Bra

V-458
Vianna Lightweight Support Bra

Vianna lightweight support bra is made of nylon and spandex. It has seamless molded contouring cups, a smooth, seamless back, and a front hook and eye clasp. It is available in XS-3XL.

NO RETURNS ON GARMENT ORDERS

FACE

Universal Facial Bands

Our universal facial bands have adjustable contouring face and neck bands and hook and loop clasps for a soft and durable stretch. They come in white and champagne and are available in S-XL.



V-210-2
2" Universal Facial Band



V-2100B
Occipital Universal Facial Band



V-210
4" Universal Facial Band



V-210S
Universal Split Ear Facial Band

Universal split ear facial band has adjustable and detachable anchor straps, split ear design, and hook and loop clasps for a soft and durable stretch. It is one size fits all.



V-210U
Universal Facial Band with Cotton Lining

Universal facial band has adjustable and detachable anchor straps and hook and loop clasps for a soft and durable stretch. It has a soft cotton lining and is one size fits all.



V-210K
Universal Facial Band with Cold/Hot Compress

Universal facial band with adjustable and detachable anchor straps as well as detachable cotton knit pouches for cold or hot compresses. Hook and loop clasps create a soft and durable stretch. It is one size fits all.



V-210N
Universal Facial Band

Universal facial band has adjustable contouring face and neck bands and hook and loop clasps for a soft and durable stretch. It is available in S-XL.

NO RETURNS ON GARMENT ORDERS

BODY

Zippered Body Girdles

Zippered high back girdles with bra have front zipper with cotton guard tape, soft molded cotton knit bra without underwire and targeted compression panels. They have adjustable hook & eye shoulder straps, reinforced open crotch, and gripper lace on leg. Breathable powernet fabric retains compressive memory. They are available in XS-4XL.



V-B853
Zippered High Back Girdle with Bra - Above Knee



V-B855
Zippered High Back Girdle with Bra - Below Knee



V-850
Zippered Abdominal Girdle



V-853
Zippered Body Girdle

Zippered Body Girdles

Zippered girdles with side zippers and cotton guard tape, abdominal compression panels, and adjustable and detachable suspenders. Girdles have gripper lace on leg, reinforced open crotch, and breathable powernet fabric that retains compressive memory. They are available in XS-4XL.

NO RETURNS ON GARMENT ORDERS

BODY



V-753
Zippered Above Knee Girdle

Zippered Above Knee Girdle has side separating zippers with cotton guard tape, abdominal compression panels, and 3”-wide waistband. It has gripper lace on leg, a reinforced open crotch, and breathable powernet fabric that retains compressive memory. It is available in XS-4XL.



V-B85
Zippered High Back Torso Brief with Bra

Zippered high back torso brief with bra has a front zipper with cotton guard tape, soft molded cotton knit bra without underwire and targeted compression panels. It has adjustable hook & eye shoulder straps and crotch closure and breathable powernet fabric that retains compressive memory. It is available in XS-4XL.



V-640
Adjustable Compression Vest

Adjustable Compression Vest has adjustable front hook & eye clasp, wide, adjustable hook & loop shoulder straps, and 2” wide waist band. It has a breathable powernet fabric that retains compressive memory and is available in S-4XL.



V-840
Zippered Abdominal Garment

Zippered Abdominal Garment has side zippers with cotton guard tape, targeted compression panels in key surgical sites, and side, adjustable hook & loop shoulder straps. It has a reinforced open crotch, breathable powernet fabric that retains compressive memory and is available in S-4XL.

NO RETURNS ON GARMENT ORDERS


NASAL AND SINUS PACKING

Our nasal and sinus packing products are the perfect product for facial aesthetic procedures like rhinoplasties and septoplasties. They are sterile, expand to various sizes, and are specifically designed for comfortable insertion and fast removal.

AbsorbENT™ Nasal and Sinus Packs

Single-use nasal and sinus packs are made of sterile Polyvinyl Alcohol (PVA) foam that provides gentle but firm compression to nasal mucosa. They are lint and fiber free and have interconnected open cells that are soft and pliable for epistaxis management.

Packs expand quickly and evenly and have a string attached for retention and easy removal. Dimensions listed are for the hydrated state while inserted.




PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
RH-7530-10	Ultra Thin	3.5 x 1.2 x .6	10
RH-7532-10	Ultra Thin	3.5 x 1.2 x .9	10
RH-2315-20	Small	5.5 x 1.5 x 2	20
RH-7401-10	Small	4.5 x 1.5 x 2	10
RH-7408-10	Small with Airway	4.5 x 1.5 x 2	10
RH-7404-10	Large	8 x 1.5 x 2	10
RH-7405-10	Large with Airway	8 x 1.5 x 2	10
RH-7409-10	Large Anatomical	8 x 1.5 x 3	10
RH-7410-10	Large Anatomical with Airway	8 x 1.5 x 3	10
RH-7411-10	Large Anatomical - Thin	8 x 1 x 3	10
RH-7406-10	Large - Thin	8 x 1 x 2	10
RH-2312-10	Extra Large without String	10 x 1.5 x 2.5	10

Series 5000 Non-Stick


PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
15-105	5000 Non-Stick, Small	6 x 1.5 x 2	10
15-110	5000 Non-Stick, Large	8 x 1.5 x 2	10



Compressed



Hydrated




Anatomical with Airway

NASAL AND SINUS PACKING

Rhino Rocket® Nasal Packs with Applicator

Rhino Rocket Nasal Packs are made of sterile PVA Expandacell® foam that provides epistaxis management after trauma, nasal or sinus surgery. Rhino Rockets provide gentle but firm mucosal compression with a tamponade of anterior or middle segment epistaxis. They take approximately 30 seconds to position and insert and swell to up to six times their compressed diameter when wet. Attached string for retention and easy removal. Dimensions listed are for hydrated state while inserted.



- Aids in Insertion
- Patented Applicator
- Protects Foam from Moisture
- Minimizes Blood Containment Concerns

Slimline

PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
11S-S0300-08AS	Small	3 x 1 x 2	8
11S-S0500-08AS	Medium	5 x 1 x 2	8
11S-S0800-08AS	Large	8 x 1 x 2	8
11S-S0800-4PAS	Large Pairs	8 x 1 x 2	8
11S-S1000-08AS	Extra Large	10 x 1 x 2	8

Wideline

PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
11W-W0300-08AS	Small	3 x 1.5 x 2	8
11W-W0500-08AS	Medium	5 x 1.5 x 2	8
11W-W0800-08AS	Large	8 x 1.5 x 2	8
11W-W1000-08AS	Extra Large	10 x 1.5 x 2	8

Starter Pack


PART NO.	SIZE	DIMENSIONS	QTY/BOX
11S-STARTERRHINO	Variety pack for new users. One of each size.		8

EPISTAXIS BALLOONS

Epistaxis® Catheters

Double balloon catheter made of medical-grade silicone that controls intranasal bleeding. Single-use catheter is completely sterile and is available with or without a syringe.

PART NO.	SIZE	DESCRIPTION	QTY/BOX
RH-4200-01	97mm	No Syringe	1
RH-4200-10	97mm	No Syringe	10
RH-4220-01	97mm	With 20cc Syringe	1
RH-4400-01	102mm	No Syringe	1
RH-4400-10	102mm	No Syringe	10
RH-4420-01	102mm	With 20cc Syringe	1




INTERNAL NASAL SPLINTS

Breeze™ Nasal Airway Splint

The Breeze Nasal Airway Splint is a single-use product made of medical grade silicone that is commonly referred to as a “Doyle Splint.” Its unique airway geometry creates a pressure gradient force that increases airflow. It has suture indicators located at the anterior tip to aid in suture placement. The splint has a lubricious coating straight out of the package that aids insertion and removal.

PART NO.	DESCRIPTION	QTY/BOX
SP-79300	With Airway, Sterile	1 Pair
SP-79300-05	With Airway, Sterile	5 Pairs
SP-79300-50	With Airway, Sterile	50 Pairs
SP-79200-B	With Airway, Non-Sterile	50 Pairs
SP-79400	Without Airway, Sterile	1 Pair



Nasal Septal Button

Sterile, single-use Nasal Septal Buttons are made of medical grade silicone and can be easily trimmed to ensure proper fit using the 1 cm measurement rings. The Nasal Septal Button’s one-piece design with 2 flange provides maximum security and cleanliness for nasal septal perforations. 6mm posts on all sizes.

PART NO.	DESCRIPTION	DIAMETER AND THICKNESS	QTY/BOX
SP-78100	Small	3cm x 4mm	1
SP-78105	Medium	5cm x 4mm	1
SP-78107	Large	7cm x 4mm	1

Bi-Valve Nasal Splint

Sterile, single-use Bi-Valve Nasal Splints are made of clear medical grade fluoroplastic. Their non-stick surface reduces adhesion and/or hematoma, has a medial slit that facilitates insertion and removal and holes for easy suturing.

Standard Size: 26.1 mm x 42.6 mm

Oversize: 39.11 mm x 52 mm

PART NO.	DESCRIPTION	QTY/BOX
SP-27000	Standard Thin (.25mm)	6 Pairs
SP-27005	Standard Thick (.5mm)	6 Pairs
SP-27010	Oversized Thin (.25mm)	6 Pairs
SP-27015	Oversized Thick (.5mm)	6 Pairs





DENVER SPLINTS®

The Denver Nasal Splint stabilizes the nose after surgery reducing edema after nasal surgery or trauma. The moldable aluminum retains nasal bones in position and multiple layers provide cushioned pressure on subcutaneous tissue. The dorsal pad reduces dead space and the possibility of polybeak. The Denver Splints come in four sizes: mini, petite, small/medium, large and are sold standalone or in kits.

Denver Splint® Series 1500

The Denver Splint Series 1500 is made of soft moldable aluminum and Velcro™ and can be worn for 7-10 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-1500-05KM	Kit, Mini	50mm x 28mm	5
10-1500-05KP	Kit, Petite	50mm x 30mm	5
10-1500-05KS	Kit, Small/Medium	65mm x 35mm	5
10-1500-05KL	Kit, Large	69mm x 43mm	5
10-1500-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-1500-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-1500-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-1500-10NKL	Non-Kit, Large	69mm x 43mm	10



Adhesive

Moldable Aluminum

No Water Mess


Inconspicuous

DENVER SPLINTS®

Denver Splint® Series 2000 - One Piece


The Denver Splint Series 2000 is made of soft moldable aluminum and 1/8” of ventilating foam and can be worn for 5-7 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-2000-05KM	Kit, Mini	50mm x 28mm	5
10-2000-05KP	Kit, Petite	50mm x 30mm	5
10-2000-05KS	Kit, Small/Medium	65mm x 35mm	5
10-2000-05KL	Kit, Large	69mm x 43mm	5
10-2000-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-2000-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-2000-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-2000-10NKL	Non-Kit, Large	69mm x 43mm	10



KIT INCLUDES:

- 1 Denver Splint®
- 2 Pre-cut Protecto Strips, sterile
- 1 Alcohol wipe, sterile
- 1 Skin prep, sterile
- 1 Dorsal pad



Denver Splint® Series 5000

The Denver Splint Series 5000 is made of soft moldable aluminum and Velcro™ and has an attached dorsal pad. It can be worn for 7-10 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-5000-05KM	Kit, Mini	50mm x 28mm	5
10-5000-05KP	Kit, Petite	50mm x 30mm	5
10-5000-05KS	Kit, Small/Medium	65mm x 35mm	5
10-5000-05KL	Kit, Large	69mm x 43mm	5
10-5000-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-5000-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-5000-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-5000-10NKL	Non-Kit, Large	69mm x 43mm	10

Denver Splint® Series 1800

The Denver Splint Series 1800 is made of soft moldable aluminum and Velfoam. It can be worn for 7-10 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-1800-05KM	Kit, Mini	50mm x 28mm	5
10-1800-05KP	Kit, Petite	50mm x 30mm	5
10-1800-05KS	Kit, Small/Medium	65mm x 35mm	5
10-1800-05KL	Kit, Large	69mm x 43mm	5
10-1800-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-1800-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-1800-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-1800-10NKL	Non-Kit, Large	69mm x 43mm	10

THERMOPLASTIC NASAL SPLINTS




Pre-Cut Denver Splint® Series 4000


- Rigid, opaque, lightweight material
- Soften with hot water and mold over taped nose

KIT INCLUDES:

- 1 Denver Splint®
- 2 Pre-cut Protecto Strips, sterile
- 1 Alcohol wipe, sterile
- 1 Skin prep, sterile
- 1 Dorsal pad



INDICATE SIZE



PART NO.	SIZE	DESCRIPTION	QTY/BOX
10-4021-05KL	Large	Kit, Ivory 1/16”	5
10-4021-05KM	Mini	Kit, Ivory 1/16”	5
10-4021-05KP	Petite	Kit, Ivory 1/16”	5
10-4021-05KS	Small	Kit, Ivory 1/16”	5
10-4021-10NKL	Large	Non-Kit, Ivory 1/16”	10
10-4021-10NKM	Mini	Non-Kit, Ivory 1/16”	10
10-4021-10NKP	Petite	Non-Kit, Ivory 1/16”	10
10-4021-10NKS	Small	Non-Kit, Ivory 1/16”	10
10-4022-05KL	Large	Kit, Blush 1/16”	5
10-4022-05KM	Mini	Kit, Blush 1/16”	5
10-4022-05KP	Petite	Kit, Blush 1/16”	5
10-4022-05KS	Small	Kit, Blush 1/16”	5
10-4022-10NKL	Large	Non-Kit, Blush 1/16”	10
10-4022-10NKM	Mini	Non-Kit, Blush 1/16”	10
10-4022-10NKP	Petite	Non-Kit, Blush 1/16”	10
10-4022-10NKS	Small	Non-Kit, Blush 1/16”	10
10-4031-05KL	Large	Kit, Ivory 3/32”	5
10-4031-05KM	Mini	Kit, Ivory 3/32”	5
10-4031-05KP	Petite	Kit, Ivory 3/32”	5
10-4031-05KS	Small	Kit, Ivory 3/32”	5
10-4031-10NKL	Large	Non-Kit, Ivory 3/32”	10
10-4031-10NKM	Mini	Non-Kit, Ivory 3/32”	10
10-4031-10NKP	Petite	Non-Kit, Ivory 3/32”	10
10-4031-10NKS	Small	Non-Kit, Ivory 3/32”	10

THERMOPLASTIC NASAL SPLINTS

Thermoplast Splints


Thermoplast splints are made of rigid but lightweight material that softens with hot water allowing them to be affixed to the patient. Splints are available in pre-cut sizes (3"x3") or in sheets to cut your own shape (9"x12").

3x3" Squares

PART NO.	DESCRIPTION	QTY/BOX
AQ-F-05-21	Ivory 1⁄16"	5
AQ-F-10-21	Ivory 1⁄16"	10
AQ-F-05-31	Ivory 3⁄32"	5
AQ-F-10-31	Ivory 3⁄32"	10
AQ-F-05-22	Blush 1⁄16"	5
AQ-F-10-22	Blush 1⁄16"	10

9x12" Sheets

PART NO.	DESCRIPTION	QTY/BOX
AQ-9x12-21	Ivory 1⁄16"	1
AQ-9x12-31	Ivory 3⁄32"	1
AQ-9x12-22	Blush 1⁄16"	1



SINUS PACK AND SPLINTING CROSS REFERENCE

All nasal and sinus packs are made from sterile, polyvinyl alcohol (PVA) foam. This chart offers easy cross referencing and ordering.

OUR PART NO.	COMPETITORS PART NO.	DESCRIPTION
NASAL M-TYPE		
RH-7401-10	440400/400400	M-Type Small
RH-7408-10	440409/400409	M-Type Small with Airway
RH-7404-10	440402/400402	M-Type Large
RH-7405-10	440410/400410	M-Type Large with Airway
NASAL P-TYPE		
RH-2310-10	440406 (Pope)/400406	P-Type XL
RH-2315-20	400421 (Weimert)	P-Type Small
NASAL D-TYPE		
RH-7409-10	440424 (Doyle)/400424	D-Type Wideline Anatomical
RH-7411-10	440413 (Doyle)/400413	D-Type Slimline Anatomical
SINUS NASAL K & T TYPES		
RH-7530-10	400530 (Kennedy)	K-Type, Ultra-Slim
RH-7532-10	400428 (Toffel)	T-Type, Slim
EAR PACKS		
OT-5201-50	400141	Ear Speedi-Wick 15mm
OT-5202-50	400146	Ear Speedi-Wick 25mm
30-305	425008	Ear Fenestrated 15mm
OT-5204-10	400355	Ear Post-Op Small
OT-5203-10	400162	Ear Post-Op Large
DENVER SPLINTS		
10-2000-05KL	1528136	Denver Splint, 2000 Kit Large
10-2000-05KS	1528126	Denver Splint, 2000 Kit Sm/Med
10-2000-05KP	1528116	Denver Splint, 2000 Kit Petite
10-2000-05KM	1528116	Denver Splint, 2000 Kit Mini
10-2000-10NKL	1528137/1528138	Denver Splint, 2000 Non-Kit Large
10-2000-10NKS	1528120/1528121	Denver Splint, 2000 Non-Kit Sm/Med
10-2000-05NKP	1528110/1528111	Denver Splint, 2000 Non-Kit Petite
10-2000-05NKM	1528110/1528111	Denver Splint, 2000 Non-Kit Mini
THERMOPLAST SPLINTS		
10-4021-05KL	1529120	Denver Splint, 4000 Kit Large
10-4021-05KS	1529110	Denver Splint, 4000 Kit Sm/Med
10-4021-05KP	1529100	Denver Splint, 4000 Kit Petite
10-4021-10NKS	1529010	Denver Splint, 4000 Non-Kit Sm/Med
10-4021-10NKP	1529000	Denver Splint, 4000 Non-Kit Petite
10-4022-10NKL	1529020	Denver Splint, 4000 Non-Kit Large

*All of our packs have an attached string



Shippert Aesthetics
Aesthetic & Reconstructive

Summit Medical LLC
An Innovia Medical Company
815 Vikings Parkway
Suite 100
St. Paul, MN 55121 | USA
1.888.229.2875
www.innoviamedical.com/aesthetic-reconstructive

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and CMDCAS certified.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

Denver Splint, Tissu-Trans and Rhino Rocket are all registered trademarks of Shippert Medical. AbsorbENT and Breeze are all trademarks of Summit Medical.

All brands are trademarks or registered trademarks of their respective owners.

©2021 Innovia Medical. All rights reserved. M-81 | RevA | NOV2021

