



# INTRUSAFE

INSTRUMENT PROTECTION TRAYS

*Your Trusted Leader for Instrument Longevity.*



Innovia Medical is a **family of specialty surgical companies** that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

## **ELEVATING CARE. IMPROVING OUTCOMES.**

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

## **INNOVATION LEADERSHIP.**

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

## **UNSURPASSED PRODUCT QUALITY.**

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

## **RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.**

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.


We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at [www.innoviamedical.com](http://www.innoviamedical.com).


THE INSTRUSAFE ADVANTAGE	
360° OF PROTECTION	4
INDUSTRY DISADVANTAGES	6
FEATURES + BENEFITS	8
ORDER PROCESS: THE INSTRUSAFE WAY	10
STERILIZATION BARRIERS + FOOT OPTIONS	12
TRAYS BY PROCEDURE	
ARTHROSCOPY	14
BARIATRIC	16
CARDIOVASCULAR	18
DENTAL	20
ENT	22
GENERAL	34
GYNECOLOGY	38
NEUROLOGY/SPINE	40
OPHTHALMOLOGY	46
ORTHOPEDIC	50
PLASTICS	52
RETRACTORS	54
ROBOTICS	58
SCOPES	70

TRANSPORT CONTAINERS	
XL TRANSPORT CONTAINERS	81
HEAVY DUTY TRANSPORT CONTAINERS	83
TRANSPORT CONTAINERS	84
ACCESSORIES	
SILICONE SMALL PARTS BOXES	85
PLAQUE OPTIONS	86
FINGER MATS	87
GUARDIAN TRAY	87
INSTRU-ZYME	88
ROUNDERS POINT GUARDS	88

**RESOURCES**



POINT CAMERA TO  
SEE VALIDATIONS



POINT CAMERA  
TO READ IFU

**SUBMIT PURCHASE ORDERS TO:**

Summit Medical LLC  
 815 Vikings Parkway, Suite 100  
 St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939  
 F: 1-888-229-1941 | +1 651-789-3979  
[www.summitmedicalusa.com](http://www.summitmedicalusa.com)



# The InstruSafe® Advantage

## 360° of Protection

The use of complex and delicate instrument sets is a new reality in which sterile processing and the operating room live. Intricacy of surgical instrumentation continues to increase, and so does the cleaning and sterilization process. Now, healthcare facilities face the challenge of reprocessing instruments safely and efficiently while also meeting the demand for fast turnaround. Keeping these sets organized and protected upon arrival in the operating room creates a challenge for all involved. Relieve yourself of these daily obstacles with InstruSafe Trays—an everyday solution for the intricacies of your instrument sets.

We have spent over 30 years understanding the use cycle and improving outcomes. With the use cycle in mind, we have created a product that eliminates common frustrations associated with instrument sets. The products in this catalog are a compilation of our hard work. As you journey through the pages, let us help you solve your everyday instrument problems.

Summit Medical's most recent FDA 510(k) clearances and product advancements are a testament of our dedication to our customers' satisfaction. We wouldn't have continuing success without the healthcare community we serve. Focus, determination and your support will continue our impact on patient care for decades to come.



*Designed to organize and protect from the point of operating room use to decontamination, our solutions benefit everyone who handles our products.*






**FDA 510(k) Clearance**

Designed to ensure optimal sterilization, InstruSafe Trays have undergone rigorous testing before receiving multiple FDA 510(k) clearances. Cleared for various sterilization cycles, these reliable trays give hospitals and healthcare facilities even more options when it comes to organizing, sterilizing, and protecting their instruments.

For the most up-to-date list of cleared cycles, visit [www.instrusafe.com/validations](http://www.instrusafe.com/validations).



POINT CAMERA TO SEE VALIDATIONS



**DECONTAMINATION**

Instruments are individually protected going into the wash cycle and the trays are designed with large perforations, allowing maximum detergent and water flow.

**TRAY PREPARATION**

Each tray is designed specifically for the instruments within it, providing individual locations for each one. This protects them from harming the technician, from damaging themselves and from puncturing or jeopardizing the integrity of the sterile barrier.

**STERILIZATION**

The large amount of perforations in our trays provide superior sterilant penetration. Many of our silicone parts are also designed with minimal contact points and additional ventilation in each slot that match our tray holes to improve sterilant penetration to each individual instrument.

**TRANSPORTATION**

Soft and durable silicone holders in the base and cover secure instruments to prevent shifting and provide 360 degrees of protection. InstruSafe Trays help to avoid unnecessary damage and to present an organized instrument set every time.

**OPERATING ROOM USE**

Each instrument is easily identifiable and placed in a singular slot for surgical presentation and tracking at the end of a procedure.

# Look Familiar?

Unfortunately, this is what most sterile processing departments have come to know. The familiarity and acceptance of these messy, hazardous baskets makes us cringe! Why, you ask? Because the minimal organization and protection of surgical instruments practiced in many facilities presents numerous risks to the individuals that come in contact with them, as well as to the instruments themselves.





# DISADVANTAGES OF MESSY BASKETS

- POTENTIAL TO COMPROMISE STERILE BARRIER BY POKING THROUGH WRAP □
- DANGEROUS TO TRANSPORT FROM THE O.R. TO THE SPD □
- MINIMAL INSTRUMENT PROTECTION AND ORGANIZATION □
- DIFFICULT TO COUNT AND TRACK INSTRUMENTS □
- BASKETS NEED FREQUENT REPLACEMENT □
- ADDS TIME TO SURGICAL SET-UP □
- BROKEN INSTRUMENTS IF DROPPED □
- HIGH RISK FOR FINGER POKES □
- POOR SURGICAL PRESENTATION □
- HARD TO FIND WHAT YOU NEED □
- MISSING/LOST INSTRUMENTS □
- UNNECESSARY REPROCESSING □
- DAMAGED INSTRUMENTS □
- WASTE OF SPACE □



## EXPECT MORE.

For these reasons and many more, InstruSafe Trays offer custom solutions to meet each customer's specific needs. We think it's time for everyone to expect more from their trays.

# Absolute Protection

## Trusted Solutions + Superior Results



### QUARTER TURN LATCH

Secures the tray cover to the tray base, providing visual indication of the open/close position. Our off-set latch hole aligns the cover in the correct position ensuring the silicone holders match up with the instruments below.




### RETRACTABLE HANDLES

Our multi-functional space-saving design fits into rigid containers making it effortless to remove the tray.



### TF FOOT

Our Tall Foot fits into rigid containers providing clearance for the filter in the base.




### FOOT OPTIONS


Sterile barrier and storage needs determine the foot type built into a tray.

### R FOOT

Rounders Corner Guards prevent tears in sterilization wrap during storage and transportation.







POINT CAMERA  
TO READ IFU


**PERFORATIONS (GRID PATTERN)**

Each tray contains a large quantity of perforations that provide maximum sterilant penetration and ventilation.



**VENTILATION**




Our unique silicone design has channels that align with our perforations allowing detergent and sterilant to flow through the holder and access the entire instrument.



**SILICONE HOLDERS**


InstruSafe silicone holders are soft yet durable, with a history of longevity. Our variety of holders accommodate instruments large or small, allowing each instrument to have its own location within the tray.

Holders in the tray base and cover provide 360 degrees of protection reducing instrument repair and replacement while delivering instruments to the O.R. in their optimal working condition.



**REDUCED CONTACT POINTS**

Each silicone holder is individually designed with ridges and channels built to reduce contact with the instrument while securely holding the instrument in place.



# 5 Simple Steps

If one of our preconfigured tray designs does not fit your needs let our product experts design a custom tray with our effortless process below.




## 1 ASSESS

An InstruSafe Representative conducts a free, on-site assessment of your instrument set.




## 2 GATHER

Your instrument set specifics and inputs are gathered by our team of experts.



## 3 DESIGN

Using your inputs, our product specialists build a tray design with you and your instrument set in mind.





**EVALUATE**

You be the judge! Evaluate the new tray design while it's loaded with your instrument set. Let us know your thoughts and we will make any and all adjustments necessary!



**COMPLETE**

It's all yours! Once complete, each custom tray is assigned its own unique part number, making it simple and easy for you to reorder in future.

## HOW TO USE THE CATALOG


The colored headers identify which procedure the instrument set you're looking at is used for. The trays on each page are differentiated by preconfigured or custom configured with the symbols shown in the key below. Some trays have several layers including a cover, shelves, baskets, and a base. The custom configured tray examples included in the catalog are a small representation of the trays our local representatives have built for facilities around the world.

### CUSTOM + PRECONFIGURED SYMBOLS



**THIS SYMBOL IDENTIFIES A PRECONFIGURED TRAY**

Order by part number from your local InstruSafe representative or call us today.



**THIS SYMBOL IDENTIFIES A CUSTOM CONFIGURED TRAY EXAMPLE**

Call us today to schedule your free assessment.



# Sterilization Barriers





InstruSafe Instrument Protection Trays are FDA cleared for several sterilization cycles using wrap and/or rigid containers.

For the most up-to-date list of FDA 510(k) cleared cycles, visit [www.instrusafe.com](http://www.instrusafe.com).

# Foot Options

**R Foot:** Rounders Corner Guards—to be used with wrap

**TF Foot:** Tall Foot fits into rigid containers

STERILE BARRIER	FOOT OPTION
<p data-bbox="522 851 604 874">WRAP</p> 	<p data-bbox="1086 851 1420 874">Rounders Corner Guards -R</p> 
<p data-bbox="487 1184 638 1207">CONTAINER</p> 	<p data-bbox="1177 1184 1328 1207">Tall Foot -TF</p> 

## HOW TO ORDER:

### EXAMPLE: TRAY IN-2890

<p data-bbox="435 1919 652 2006">Wrapped Foot IN-2890-<b>R</b></p>	<p data-bbox="1043 1919 1260 2006">Container Foot IN-2890-<b>TF</b></p>
--	---

FOOT OPTION IS ADDED TO THE END OF THE PART NUMBER



# Rounders® Corner Guards

Rounded, medical grade corner guards affix to InstruSafe Tray corners to reduce potential damage to the sterile wrap barrier. Rounders are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to all InstruSafe Trays.



# Tall Foot

Designed to affix to trays that are intended to be sterilized within a rigid container. Tall Feet are included with the purchase of designated tray orders, but can also be ordered at any time and easily added to most InstruSafe Trays.



## ROUNDERS CORNER GUARDS

<i>Small Rounders</i>	
Kit of Four	IN-7334
<i>Large Rounders</i>	
Kit of Four	IN-7337

## TF FOOT KIT

Kit of Four	IN-7338
-------------	---------

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit all InstruSafe Trays
- Rounded, medical grade silicone corner guard
- Help avoid compromised sterility in the O.R.
- Environmentally and economically efficient
- Impact resistant
- US Patent #9119611 | #D783856

- Kits can be purchased to upgrade existing trays
- Affix to tray with screws (included in kit)
- Designed to fit most InstruSafe Trays



## Preconfigured Tray for Arthroscopy

This three-level tray is equipped for both small and large arthroscopic instrument sets. Our tray is sized to fit most three-quarter size containers or can be wrapped for your convenience.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

Cover



Shelf



Basket



Base



### ARTHROSCOPY INSTRUMENT TRAY

#### *Arthroscopy Instrument Tray*

Holds:

A combination of up to 12 of the following instruments:


- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)	387 L x 244 W x 146 H (mm)
Wrapped Foot IN-2846-R	Container Foot IN-2846-TF

Expand capacity with positionable Hold-Its and Hold-Downs™ silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

Basket and base can be customized to accommodate specific instrumentation.




# Custom Configured Tray Example for Arthroscopy Example

## ARTHROSCOPY TMJ TRAY | 17 INSTRUMENTS

### BEFORE



CUSTOM CONFIGURED



### AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





## Preconfigured Tray for Bariatric

The bariatric tray is designed to accommodate even the longest instruments. Our three-level system allows you to get organized and protect all of your bariatric instrumentation in one tray.



### BARIATRIC TRAY

#### *Bariatric Instrument Tray*

Holds:  
A combination of up to 15 (32cm to 45cm+) of the following instruments:

- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars


28 L x 13 W x 8 H (in)	711 L x 330 W x 203 H (mm)
Wrapped Foot IN-8864-R	Container Foot IN-8864-TF

Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.


**CALL US TODAY FOR CUSTOMIZATION!**

### CUSTOM CONFIGURED BARIATRIC TRAY | 31 INSTRUMENTS



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHRO
- BARIATRIC**
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



# Custom Configured Tray Example for Bariatric

BARIATRIC LAPAROSCOPIC SET | 59 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

## BEFORE




CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Cardiovascular

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS



ACCESSORIES

## MICROVASCULAR SET | 19 INSTRUMENTS

### BEFORE




### CUSTOM CONFIGURED






### AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Cardiovascular

EVANS ROTATABLE TRAY | 32 INSTRUMENTS



Shelf

Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Dental

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

## DENTAL SET | 34 INSTRUMENTS

# BEFORE




CUSTOM CONFIGURED



# AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.






# Custom Configured Tray Example for Dental

## OPERATIVE TRAY SET | 11 INSTRUMENTS

CUSTOM CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

# ENT SYSTEMS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

Since the mid-1980's, we have offered a complete solution for fine and delicate ear instrumentation through our Middle Ear System. We are proud to claim this as our very first InstruSafe product. Our modules are designed to specifically fit each instrument and can be used individually or in combination, to organize and protect an entire instrument set.

**MIX AND MATCH OUR MODULES FROM PAGES 23-24 TO FIT YOUR INSTRUMENT SET, CREATING YOUR UNIQUE ENT SYSTEM!**




IN-2004




IN-2001



IN-2003



IN-2008



Tray shown without cover for illustrative purposes.



## Preconfigured Trays for ENT Modules




### MIDDLE EAR TRAY

#### *Middle Ear Tray - 2 Level*

Holds:

- 2 Alligator Forcep Modules
- 2 Pick + Hook Modules
- 1 Micro-Suction Module
- 1 Speculum Module
- 9½ x 13" Finger Mat

19 5/8 L x 9 5/8 W x 6 H (in)	498 L x 244 W x 152 H (mm)
Wrapped Foot IN-2890-R	Container Foot IN-2890-TF



### EAR MODULE TRAY

#### *Alligator Module*

Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

5 7/8 L x 5 1/8 W x 1 H (in)	149 L x 130 W x 25 H (mm)
IN-2001	



### EAR MODULE TRAY

#### *Pick + Hook Module, Wide*

Holds:

- 10 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 5 1/8 W x 1 H (in)	194 L x 130 W x 25 H (mm)
IN-2002	



### EAR MODULE TRAY

#### *Speculum Module*

Holds:

- 6-12 Speculas

3 1/2 L x 5 1/8 W x 1 H (in)	89 L x 130 W x 25 H (mm)
IN-2003	

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





# Preconfigured Trays for ENT Modules

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## EAR MODULE TRAY

### *Micro-Suction Module*

Holds:

- 12 Micro-Suction Irrigators with Luer Lock Hubs

5 1/8 L x 3 1/2 W x 1 H (in)

130 L x 89 W x 25 H (mm)

IN-2004



## EAR MODULE TRAY

### *Suction Module*

Holds:

- 10 Ferguson Style Suctions

5 7/8 L x 5 1/8 W x 1 H (in)

149 L x 130 W x 25 H (mm)

IN-2005



## EAR MODULE TRAY

### *Alligator Module With Cover*


Holds:

- 8 Alligator Style Instruments (Cup, Scissors, etc.)

7 5/8 L x 5 1/8 W x 3 H (in)

194 L x 130 W x 76 H (mm)

IN-2006



## EAR MODULE TRAY

### *Pick + Hook Module*

Holds:

- 9 Small Handle Picks, Hooks, Dissectors or Curettes

7 5/8 L x 4 5/8 W x 1 H (in)


194 L x 117 W x 25 H (mm)

IN-2008



## Preconfigured Trays for Picks + Hooks

Similar to our original module designs, these standalone trays will protect and organize pick, hook and knife instrumentation for ENT procedures. Instruments can either be positioned at a 45 degree angle for procedural presentation or secured flat for sterilization and storage.



### PICK + HOOK TRAY

#### 10 Pick + Hook Tray

Holds:

- 10 Pick + Hook Style Instruments

5 3/4 L x 8 3/4 W x 1 1/2 H (in)

133 L x 222 W x 38 H (mm)

IN-4010



### PICK + HOOK TRAY

#### 20 Pick + Hook Tray

Holds:

- 20 Pick + Hook Style Instruments

11 L x 8 3/4 W x 1 1/2 H (in)

279 L x 222 W x 38 H (mm)

IN-4000



### PICK + HOOK TRAY

#### 30 Pick + Hook Tray

Holds:

- 30 Pick + Hook Style Instruments

15 L x 8 3/4 W x 1 1/2 H (in)

381 L x 222 W x 38 H (mm)

IN-4003

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# Custom Configured Tray Example for Ear

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS



ACCESSORIES

## MYRINGOTOMY SET | 26 INSTRUMENTS

# BEFORE




**CUSTOM CONFIGURED**



# AFTER

*Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.*








# Custom Configured Tray Example for Ear

EAR TRAY | 86 INSTRUMENTS

CUSTOM CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Ear

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS



ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

## MAJOR EAR TRAY | 66 INSTRUMENTS



Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





## Preconfigured Tray for Nose/Sinus



Basket and base can be customized to accommodate specific instrumentation.

### ENDOSCOPY SINUS TRAY

#### *Endoscopy Sinus Instrument Tray*

Holds:

A combination of up to 12 of the following instruments:

- Punches
- Ring Handle Instruments
- Other Associated Endoscopic Instruments
- Rongeurs
- Graspers

15 ¼ L x 9 ⅝ W x 5 ¾ H (in)	387 L x 244 W x 146 H (mm)
Wrapped Foot IN-2846-R	Container Foot IN-2846-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES






# Custom Configured Tray Example for Nose/Sinus

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT**
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES


## SINUS SET | 32 INSTRUMENTS

# BEFORE



# AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Nose/Sinus

SINUS TRAY | 34 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS


ACCESSORIES



Shelf



Basket



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Preconfigured Tray for Throat

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES


Cover




Shelf



Basket



Base



## ENDOSCOPY TRAY

### Full Size Endoscopy Instrument Tray

Holds:

- 14 Endoscopy Style Instruments

19 5/8 L x 9 5/8 W x 7 H (in)

498 L x 244 W x 178 H (mm)


Wrapped Foot  
IN-2847-R

Container Foot  
IN-2847-TF



Small Parts Boxes are perfect for your small delicate instruments or accessories.


Basket and base can be customized to accommodate specific instrumentation.





# Custom Configured Tray Example for Throat

THYROID SET | 25 INSTRUMENTS

## BEFORE



**CUSTOM CONFIGURED**



## AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

*Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.*




# Preconfigured Tray for General

Maximize the protection and organization of your laparoscopic instrument sets with this three-level tray that is excellently constructed to fit both small and large sets of instrumentation.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL**
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES

Cover




Shelf



Basket



Base



## GENERAL/GYNECOLOGY TRAY

### Extra Long Laparoscopic Instrument Tray

Holds:  
A combination of up to 11 of the following instruments:

- Laparoscopic Instruments
- General Dissecting Instruments
- Miscellaneous Parts
- Tubing/Cords
- Cannulas
- Sheaths
- Trocars


24 3/4 L x 9 5/8 W x 6 H (in)	629 L x 244 W x 152 H (mm)
Wrapped Foot IN-8865-R	Container Foot IN-8865-TF

Expand capacity with positionable Hold-Its and Hold-Downs silicone instrument holders.

Multi-level to accommodate a high quantity of instruments.

Basket and base can be customized to accommodate specific instrumentation.







# Custom Configured Tray Example for General

TENDON STRIPPER | 13 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# Custom Configured Tray Example for General

LAPAROSCOPIC | 47 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS



ACCESSORIES



Shelf




Basket



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.







# Custom Configured Tray Example for General

## BASIC UROLOGY TRAY | 35 INSTRUMENTS

CUSTOM CONFIGURED



Shelf



Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES


Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Custom Configured Tray Example for Hysteroscopy

## HYSTEROSCOPE SET | 10 INSTRUMENTS

### BEFORE



CUSTOM CONFIGURED




### AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN**
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES







## Custom Configured Tray Example for Hysteroscopy

### HYSTERECTOMY RESECTOSCOPE TRAY | 15 INSTRUMENTS

CUSTOM  
CONFIGURED



Shelf

Base

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





## Preconfigured Trays for Kerrisons

Our Neurology and Spine Trays provide maximum accessibility and flexibility when handling and sterilizing kerrisons and forceps.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes.



### KERRISON INSTRUMENT TRAY + RACK

#### Ten Position Instrument Tray

Holds:

- 10 Kerrisons

15 ¼ L x 9 ⅝ W x 6 H (in)

387 L x 244 W x 152 H (mm)

Wrapped Foot  
IN-5012-R

Container Foot  
IN-5012-TF

### KERRISON INSTRUMENT RACK

#### Ten Position Instrument Rack

Holds:

- 10 Kerrisons

14 ¼ L x 9 ½ W x 6 H (in)

362 L x 241 W x 152 H (mm)

Wrapped Foot  
N/A

Container Foot  
IN-5013-TF

### KERRISON INSTRUMENT RACK

#### Five Position Instrument Rack

Holds:

- 5 Kerrisons

14 ¼ L x 5 W x 6 H (in)

362 L x 127 W x 152 H (mm)

Wrapped Foot  
N/A

Container Foot  
IN-5014-TF



## Preconfigured Trays for Bi-Polar Forceps



### BAYONET FORCEPS MODULE TRAY

#### *Dual Forceps Module Tray with Cover*


Holds:

- 2 Forceps

9 ¾ L x 2 ⅛ W x 1 ⅞ H (in)

248 L x 54 W x 29 H (mm)

IN-5401-02



### BAYONET FORCEPS MODULE TRAY

#### *Dual Extra Long Forceps Module Tray with Cover*

Holds:

- 2 Forceps

12 L x 2 ⅛ W x 1 ½ H (in)

305 L x 54 W x 38 H (mm)

IN-5401-03



### BAYONET FORCEPS MODULE TRAY

#### *Six Forceps Module Tray with Cover*

Holds:

- 6 Forceps

9 ¾ L x 4 ⅝ W x 1 ⅞ H (in)

248 L x 248 W x 29 H (mm)

IN-5401

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# Preconfigured Trays for Bi-Polar Forceps



## □ BAYONET FORCEPS MODULE TRAY

### *Eight Forceps Module Tray with Cover*

Holds:

- 8 Forceps

12 L x 7 ¼ W x 1 ½ H (in)	305 L x 184 W x 38 H (mm)
---------------------------	---------------------------

IN-5401-08



## □ BAYONET FORCEPS MODULE TRAY

### *Twelve Forceps Module Tray with Cover*

Holds:

- 12 Forceps

12 L x 9 ⅝ W x 1 ½ H (in)	305 L x 244 W x 38 H (mm)
---------------------------	---------------------------

IN-5401-12

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES







# Custom Configured Tray Example for Neurology/Spine

NEUROLOGY TRAY | 13 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Neurology/Spine

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS



ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

## HIGH SPEED PNEUMATIC DRILL TRAY | 9 INSTRUMENTS




Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Neurology/Spine

## DURAL REPAIR KIT TRAY | 18 INSTRUMENTS

CUSTOM CONFIGURED



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.








# Custom Configured Tray Example for Ophthalmology

## CATARACT-PHACO SET | 28 INSTRUMENTS


### BEFORE



### AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE**
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



Custom Configured Tray Example for  
**Ophthalmology**

EYE TRAY | 57 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

**EYE**

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

**CUSTOM CONFIGURED**




Shelf



Base



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.






## Custom Configured Tray Example for Ophthalmology

### CORNEA TRAY | 61 INSTRUMENTS



Shelf




Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE**
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES





Custom Configured Tray Example for  
**Ophthalmology**

PHACO HAND PIECE TRAY | 23 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

**EYE**

ORTHO

PLASTICS



RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



## Preconfigured Trays for Orthopedic

While not usually delicate, orthopedic instrumentation needs to be kept sharp and straight to get the job done. Our Ortho Trays are built to safely accommodate and protect instruments that are commonly found in almost all surgical facilities.



### ORTHOPEDIC TRAY

#### *Lambotte Osteotome Tray*

Holds:

- 12 Instruments—straight or curved, up to 1 1/2" width (38mm)

8 3/8 L x 9 5/8 W x 3 H (in)	213 L x 244 W x 76 H (mm)
Wrapped Foot IN-5010-R	Container Foot IN-5010-TF



### CURETTE TRAY

#### *Curette Tray*

Holds:

- 14 Instruments—7-13" in length

15 1/4 L x 9 5/8 W x 2 H (in)	387 L x 244 W x 51 H (mm)
Wrapped Foot IN-5011-R	Container Foot IN-5011-TF

ARTHO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





## Custom Configured Tray Example for Orthopedic

SMALL FRAG SET | 23 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.



# Custom Configured Tray Example for Plastics





MICRO SET | 32 INSTRUMENTS

## BEFORE




CUSTOM  
CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

- ARTHO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS**
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



# Custom Configured Tray Example for Plastics

PLASTIC TRAY | 38 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Preconfigured Trays for Bookwalter® Retractors

This tray system is engineered to have a balanced weight distribution that meets the AAMI 25 lb. weight limit standard. We offer a two-level tray design made to accommodate the ring and blade components with an additional tray to store the table post and horizontal bar. Our two-tray system is organized for easy counting and safe for transporting. **Note:** *The Mediflex Bookler Retractor also fits in this tray.*



## BOOKWALTER RETRACTOR TRAY

### Ring + Blade Tray

Holds:

- 1 (1") Post Coupling
- 1 Balfour Retractor
- 1 Horizontal Flex Bar
- 1 Harrington Retractor
- 1 Medium Oval Ring
- 1 Segmented Ring
- 1 Gelpi Retractor
- 2 (2x6"/51mm x 152mm) Malleable Retractors
- 1 (3x6"/76mm x 152mm) Malleable Retractor
- 2 (2x3"/51mm x 76mm) Kelly Retractors
- 1 (2x4"/51mm x 110mm) Kelly Retractor
- 1 (2x5"/51mm x 127mm) Kelly Retractor
- 1 (2x6"/51mm x 152mm) Kelly Retractor
- 6 Tilt Ratchet Mechanisms
- 2 Ratchet Mechanisms

22 ½ L x 13 W x 4 ½ H (in)

572 L x 330 W x 114 H (mm)

Wrapped Foot  
IN-8700-R

Container Foot  
IN-8700-TF



## BOOKWALTER RETRACTOR TRAY

### Post Tray

Holds:

- Table Post and Horizontal Bar

25 L x 5 ¾ W x 3 H (in)


635 L x 137 W x 76 H (mm)

Wrapped Foot  
IN-8701-R

Container Foot  
IN-8701-TF

*The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.*








## Custom Configured Tray Example for Retractors

GREENBERG SET | 17 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





# Custom Configured Tray Example for Retractors

## BOOKWALTER RETRACTOR TRAY | 38 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES


TRANSPORT CONTAINERS

ACCESSORIES

CUSTOM CONFIGURED




Shelf



Base

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.



## Custom Configured Tray Example for Retractors

### CARDIAC MIS RETRACTOR HOLDER TRAY | 2 INSTRUMENTS

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.





# Preconfigured Trays for da Vinci SP<sup>®</sup>

The storage and reprocessing of the da Vinci SP isn't limited to a single solution. That's why we've designed multiple trays specific to the system's instrumentation to accommodate a wide variety of customer preferences.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## INSTRUMENT TRAY

### Single Instrument Tray

Holds:

- 1 Instrument

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
-----------------------------------	---------------------------

Wrapped Foot IN-8921-R	Container Foot IN-8921-TF
---------------------------	------------------------------



## INSTRUMENT TRAY


### Dual Instrument Tray

Holds:

- 2 Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 95 H (mm)
-----------------------------------	---------------------------

Wrapped Foot IN-8922-R	Container Foot IN-8922-TF
---------------------------	------------------------------



## INSTRUMENT TRAY

### Four Instrument Tray

Holds:

- 4 Instruments


28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
---------------------------	----------------------------

Wrapped Foot IN-8923-R	Container Foot IN-8923-TF
---------------------------	------------------------------

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



## Preconfigured Trays for da Vinci SP




### INSTRUMENT TRAY

#### *Six Instrument Tray*

Holds:

- 6 Instruments

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8925-R	Container Foot IN-8925-TF




### PROCEDURE TRAY

#### *Four Instrument Procedure Tray*

Holds:

- 4 Instruments + Associated Accessories

28 L x 8 W x 5 3/8 H (in)	711 L x 203 W x 137 H (mm)
Wrapped Foot IN-8924-R	Container Foot IN-8924-TF



### PROCEDURE TRAY

#### *Six Instrument Procedure Tray*

Holds:

- 6 Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8926-R	Container Foot IN-8926-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



# Preconfigured Trays for da Vinci SP

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS**
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



## ACCESSORY TRAY

### *Six Accessory Tray*

Holds:

- 2 Cannulas
- 2 Obturators
- 2 Cables

15 ¼ L x 9 ⅝ W x 2 ⅜ H (in)	387 L x 244 W x 60 H (mm)
Wrapped Foot IN-8927-R	Container Foot IN-8927-TF



## CAMERA TRAY

### *Single Camera Tray*

Holds:

- 1 Camera

28 L x 11 W x 3 H (in)	711 L x 279 W x 76 H (mm)
Wrapped Foot IN-8920-R	Container Foot IN-8920-TF

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





# Preconfigured Trays for da Vinci Xi® Single-Site®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Single-Site instruments—from single instruments to procedural sets and their accommodating accessories—our trays are the recommended solution for instrument reprocessing.




## INSTRUMENT TRAY

### Single Instrument Tray

Holds:

- 1 (5mm) Xi Single-Site Instrument

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
Wrapped Foot IN-8911-R	Container Foot N/A




## PROCEDURE TRAY

### 22 Instrument Procedure Tray

Holds:

- 6 (5mm) Xi Single-Site Instruments
- 3 Rigid Obturators
- 2 Flexible Obturators
- 4 Curved Cannulas
- 3 Rigid Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

26 3/4 L x 13 W x 3 5/8 H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8912-R	Container Foot N/A



## ACCESSORY TRAY

### 16 Accessory Tray

Holds:

- 3 Rigid Obturators
- 2 Flexible Obturators
- 3 Rigid Cannulas
- 4 Curved Cannulas
- 3 Cables
- 1 Suction Irrigator Tube and Connector

15 1/4 L x 9 5/8 W x 4 1/8 H (in)	387 L x 244 W x 105 H (mm)
Wrapped Foot IN-8913-R	Container Foot IN-8913-TF

All da Vinci Xi Single-Site Trays have green covers.

Excludes IN-8911

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# Preconfigured Trays for da Vinci Xi Single-Site



- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS**
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES

## INSTRUMENT TRAY

### *Six Instrument Tray*

Holds:

- 6 (5mm) Xi Single-Site Instruments

26 <sup>3</sup> / <sub>4</sub> L x 10 W x 3 <sup>5</sup> / <sub>8</sub> H (in)		679 L x 254 W x 92 H (mm)	
Wrapped Foot IN-8914-R		Container Foot IN-8914-TF	

**PAIRS WITH OUR XL CONTAINER**

**CALL YOUR LOCAL REP TO LEARN MORE!**

All da Vinci Xi Single-Site Trays have green covers.

*Excludes IN-8911*

*da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.*



## Preconfigured Trays for da Vinci Xi®

As robotic surgery evolves, so do our trays. Constructed exclusively for da Vinci Xi Instruments—Single or Dual Staplers and their accommodating accessories—our trays are the recommended solutions for instrument reprocessing.



### INSTRUMENT TRAY

#### Single Instrument Tray

Holds:  
• 1 EndoWrist Instrument

27 1/8 L x 3 1/2 W x 2 1/8 H (in)	689 L x 89 W x 54 H (mm)
Wrapped Foot IN-8931-R	Container Foot N/A




### INSTRUMENT TRAY

#### Two Instrument Tray

Holds:  
• 2 EndoWrist instruments

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
Wrapped Foot IN-8947-R	Container Foot N/A



### INSTRUMENT TRAY

#### Three Instrument Tray

Holds:  
• 3 EndoWrist instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8948-R	Container Foot IN-8948-TF

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





# Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## INSTRUMENT TRAY

### Four Instrument Tray

Holds:  
• 4 EndoWrist Instruments

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
-----------------------------------	---------------------------

Wrapped Foot IN-8950-R	Container Foot IN-8950-TF
---------------------------	------------------------------



## INSTRUMENT TRAY

### Six Instrument Tray

Holds:  
• 6 EndoWrist Instruments

26 3/4 L x 10 W x 3 5/8 H (in)	679 L x 254 W x 92 H (mm)
--------------------------------	---------------------------

Wrapped Foot IN-8936-R	Container Foot IN-8936-TF
---------------------------	------------------------------

**PAIRS WITH OUR XL CONTAINER  
CALL YOUR LOCAL REP TO LEARN MORE!**



## INSTRUMENT TRAY

### Eight Instrument Tray

Holds:  
• 8 EndoWrist Instruments


26 3/4 L x 13 W x 3 5/8 H (in)	679 L x 330 W x 92 H (mm)
--------------------------------	---------------------------

Wrapped Foot IN-8938-R	Container Foot N/A
---------------------------	-----------------------

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



## Preconfigured Trays for da Vinci Xi



### INSTRUMENT TRAY

#### Three Instrument + Accessory Tray

Holds:  
 • 3 EndoWrist Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 3/4 H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8949-R	Container Foot IN-8949-TF




### PROCEDURE TRAY

#### Four Instrument Procedure Tray

Holds:  
 • 4 EndoWrist Instruments + Associated Accessories

28 L x 8 W x 5 7/16 H (in)	711 L x 203 W x 138 H (mm)
Wrapped Foot IN-8915-R	Container Foot IN-8915-TF

**PAIRS WITH OUR XL CONTAINER  
 CALL YOUR LOCAL REP TO LEARN MORE!**



### PROCEDURE TRAY

#### Six Instrument Procedure Tray

Holds:  
 • 6 EndoWrist Instruments + Associated Accessories

28 L x 11 W x 5 3/8 H (in)	711 L x 279 W x 137 H (mm)
Wrapped Foot IN-8916-R	Tall Foot IN-8916-TF

**PAIRS WITH OUR XL CONTAINER  
 CALL YOUR LOCAL REP TO LEARN MORE!**

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.



# Preconfigured Trays for da Vinci Xi

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## PROCEDURE TRAY

### Six Instrument Procedure Tray

Holds:

- 6 EndoWrist Instruments + Associated Accessories

26 3/4 L x 13 W x 3 3/4 H (in)	679 L x 330 W x 95 H (mm)
Wrapped Foot IN-8937-R	Container Foot N/A

## PROCEDURE TRAY

### Eight Instrument Procedure Tray

Holds:

- 8 EndoWrist Instruments + Associated Accessories

26 3/4 L x 13 W x 3 5/8 H (in)	679 L x 330 W x 92 H (mm)
Wrapped Foot IN-8935-R	Container Foot N/A

## STAPLER + ACCESSORY TRAY

### Single Stapler + Accessory Tray

Holds:

- 1 Stapler Instrument + Associated Accessories

27 1/8 L x 3 3/8 W x 3 1/2 H (in)	689 L x 86 W x 89 H (mm)
Wrapped Foot IN-8933-R	Container Foot N/A

## STAPLER + ACCESSORY TRAY

### Dual Stapler + Accessory Tray

Holds:

- 2 Stapler Instruments + Associated Accessories

26 3/4 L x 6 7/8 W x 3 11/16 H (in)	679 L x 175 W x 94 H (mm)
Wrapped Foot IN-8932-R	Container Foot IN-8932-TF

**PAIRS WITH OUR XL CONTAINER  
CALL YOUR LOCAL REP TO LEARN MORE!**

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.





## Preconfigured Trays for da Vinci Xi



### ACCESSORY TRAY

#### Ten Accessory Tray

Holds:

- 6 Cannulas or Obturators
- 1 Release Tool
- 3 Cables

19 5/8 L x 9 5/8 W x 2 H (in)

498 L x 244 W x 51 H (mm)

Wrapped Foot  
IN-8939-R

Container Foot  
IN-8939-TF



### SCOPE TRAY

#### Single Scope Tray

Holds:

- 1 (8mm) 0° or 30° Scope

24 7/8 L x 13 1/2 W x 4 3/8 H (in)

632 L x 343 W x 111 H (mm)

Wrapped Foot  
INS-8930-R

Container Foot  
N/A



### CAMERA TRAY

#### Single Camera Tray

Holds:

- 1 Hand Held Camera

15 1/4 L x 9 5/8 W x 3 H (in)

387 L x 244 W x 76 H (mm)

Wrapped Foot  
INS-8934-R

Container Foot  
INS-8934-TF



NEED A PLACE TO SOAK?  
SEE PAGE 90 FOR DETAILS

da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# Preconfigured Trays for *Senhance*<sup>™</sup> Surgical Robotic Trays

A sophisticated robotic system deserves a solution just as sophisticated to hold its instruments. That's why TransEnterix collaborated with us, the instrument protection tray experts. This line of trays is not only designed and constructed to fit the Senhance Surgical Robotic System instruments and adapters like a glove, but it is also validated specifically for this instrumentation.

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS**
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



## INSTRUMENT TRAY

### *Instrument Tray*

Holds:

- 8 Instrument Adapters
- 12 Instrument Shafts
- Cables

19 5/8 L x 9 5/8 W x 4 3/8 H (in)	498 L x 244 W x 105 H (mm)
Wrapped Foot IN-8750-R	Container Foot IN-8750-TF




## ULTRASONIC INSTRUMENT TRAY

### *Instrument Tray*

Holds:

- 1 Ultrasonic Adapter
- 1 Ultrasonic Transducer

25 L x 3 3/8 W x 3 H (in)	635 L x 137 W x 76 H (mm)
Wrapped Foot IN-8753-R	Container Foot IN-8753-TF




## ENDOSCOPE ADAPTER TRAY

### *Single Instrument Tray*

Holds:

- 1 Endoscope Adapter

19 5/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8751-R	Container Foot IN-8751-TF



## ENDOSCOPE ADAPTER TRAY

### *Dual Instrument Tray*


Holds:

- 2 Endoscope Adapters

19 5/8 L x 9 5/8 W x 2 3/8 H (in)	498 L x 244 W x 60 H (mm)
Wrapped Foot IN-8752-R	Container Foot IN-8752-TF

*Senhance is a trademark of TransEnterix Surgical, Inc. Please refer to the Senhance<sup>™</sup> Surgical Robotic System IFU for instrument reprocessing.*






# Custom Configured Tray Example for Robotic Instrumentation

## ROBOTIC LAPAROSCOPIC SET | 9 INSTRUMENTS

BEFORE



CUSTOM CONFIGURED



AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection. da Vinci, SP, Xi, EndoWrist and Single-Site are all registered trademarks of Intuitive Surgical, Inc.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





## Preconfigured Trays for Rigid Scopes

Our rigid scope trays are universal and designed to fit scopes of any size for any procedure.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS


RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



### RIGID SCOPE TRAY 4-7"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 4-7" in working length

11 L x 3 ¼ W x 2 ⅜ H (in)

279 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6011-R

Container Foot  
N/A



### RIGID SCOPE TRAY 5-13½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 5-13½" in working length

17 ½ L x 3 ¼ W x 2 ⅜ H (in)

445 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6012-R

Container Foot  
N/A



### RIGID SCOPE TRAY 7-16½"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 7-16½" in working length

20 ¼ L x 3 ¼ W x 2 ⅜ H (in)

515 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6013-R

Container Foot  
N/A



## Preconfigured Trays for Rigid Scopes



### RIGID SCOPE TRAY 9 1/2-21 1/2"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 9 1/2-21 1/2" in working length

25 1/4 L x 3 1/4 W x 2 3/8 H (in)

641 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6014-R

Container Foot  
N/A



### RIGID SCOPE TRAY 11-26"

Holds:

- Up to 3 (2-10mm) Scopes
- Accommodates scopes from 11-26" in working length

30 1/8 L x 3 1/4 W x 2 3/8 H (in)

765 L x 83 W x 60 H (mm)

Wrapped Foot  
IN-6015-R

Container Foot  
N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## Preconfigured Trays for Scope + Cable

Our scope + cable trays allow for faster loading before sterilization while also protecting the delicate fiber optics within.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



### SCOPE + CABLE TRAY 5-6"

Holds:


- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 5-6" in working length

11 L x 6 3/8 W x 2 3/8 H (in)

279 L x 162 W x 60 H (mm)

Wrapped Foot  
IN-6031-R

Container Foot  
N/A



### SCOPE + CABLE TRAY 5-13 1/2"

Holds:

- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 5-13 1/2" in working length

20 1/2 L x 5 5/8 W x 2 3/8 H (in)

521 L x 143 W x 60 H (mm)

Wrapped Foot  
IN-6032-R

Container Foot  
N/A



### SCOPE + CABLE TRAY 6 1/2-17"

Holds:

- Up to 3 (2-10mm) Scopes
- 1 Cable
- Accommodates scopes from 6 1/2-17" in working length


25 L x 5 3/8 W x 2 3/8 H (in)

635 L x 137 W x 60 H (mm)

Wrapped Foot  
IN-6033-R

Container Foot  
N/A







# Custom Configured Tray Example for Scope + Cable

RESECTOSCOPE | 11 INSTRUMENTS

## BEFORE



CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



## Preconfigured Trays for Flexible Scopes

Our durable flexible scope trays are sure to protect your delicate and expensive scopes. Our Small Parts Boxes are perfect to keep all your accessories together and protected.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



### FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6021-R

Container Foot  
N/A

### FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Scope with Attached Light Source
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6022-R

Container Foot  
N/A

### FLEXIBLE SCOPE TRAY

Holds:

- 1 Flexible Video Scope
- Accommodates scopes from 12-27" in working length

19 5/8 L x 9 5/8 W x 3 H (in)

498 L x 244 W x 76 H (mm)

Wrapped Foot  
IN-6023-R

Container Foot  
N/A



## Preconfigured Tray for Flexible Scopes



### FLEXIBLE SCOPE TRAY

Holds: <ul style="list-style-type: none"> <li>• 1 Flexible Bronchoscope</li> <li>• Accommodates scopes from 12-27" in working length</li> </ul>	
19 5/8 L x 9 5/8 W x 3 H (in)	498 L x 244 W x 76 H (mm)
Wrapped Foot IN-6024-R	Container Foot N/A

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES





## Custom Configured Tray Example for Flexible Scopes

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES

### FLEXIBLE URETEROSCOPE | 2 INSTRUMENTS

## BEFORE




### CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.




# Custom Configured Tray Example for Flexible Scopes

VIDEO LARYNGOSCOPE | 1 INSTRUMENT

## BEFORE



CUSTOM CONFIGURED



## AFTER

Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS


ACCESSORIES






## Custom Configured Tray Example for Flexible Scopes

### COLONOSCOPY TRAY | 1 INSTRUMENT



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.






# Custom Configured Tray Example for Flexible Scopes

## FLEXIBLE SCOPE TRAY | 3 INSTRUMENTS

**CUSTOM CONFIGURED**



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS


ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

*Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.*



# Custom Configured Tray Example for Flexible Scopes

## FLEXIBLE SCOPE TRAY | 1 INSTRUMENT

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS



RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



Tray shown without cover for illustrative purposes. All custom InstruSafe Trays come with a cover to provide 360 degrees of protection.

# XL Transport Containers

Transporting da Vinci SP, Xi scopes, and EndoWrist instruments appropriately and effectively is a difficult task. Use the InstruSafe XL Transport Container to get the job done right! This durable container is the perfect size, enabling you to soak and protect your instruments while they travel from the O.R. to sterile processing. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.



IN-PIC-20-D

## XL TRANSPORT CONTAINER

### *XL Transport Container + Optional Drain*

Holds:

- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 ¼ W x 6 H (in)  
864 L x 438 W x 152 H (mm)

With Drain  
IN-PIC-20-D

Without Drain  
IN-PIC-20



IN-PIC-20-D-LT

## XL TRANSPORT CONTAINER SYSTEM

### *XL Transport Container with Lifting Tray + Optional Drain*

Holds:

- Lifting Tray
- 1 da Vinci Xi or SP Endoscope OR up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

34 L x 17 ¼ W x 6 H (in)  
864 L x 438 W x 152 H (mm)

With Drain  
IN-PIC-20-D-LT

Without Drain  
IN-PIC-20-LT

Single Endoscope Configuration



Multiple Instrument Configuration

## LIFTING TRAY

### *Single Endoscope Configuration*

Holds:

- 1 da Vinci Xi or SP Endoscope + Associated Accessories

### *Multiple Instrument Configuration*

Holds:

- Up to 8 da Vinci Xi or SP EndoWrist Instruments + Associated Accessories

27 ½ L x 13 ½ W x 3 ⅞ H (in)

699 L x 343 W x 98 H (mm)

IN-PIC-34

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



# XL Transport Container Special Features

- Highly durable material designed to withstand everyday use
- Cart wash safe—maximum temperature of 230°F
- Economical and creates a dedicated space for horizontally soaking robotic instruments
- Made in the USA

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES

**Inner gasket and strong silicone latches reduce spills**

**Built-in biohazard symbol on cover**

**Dedicated dry erase “soak start time” space on cover**

**Water volume feature to accurately measure detergent solutions**

VOLUME	
20 L	5 gal.
15 L	4 gal.
10 L	3 gal.
5 L	1.5 gal.

**Quick drain feature—no need to tip container to drain**

**Accessory Basket**

**Large Handles**

**Loading Identifiers**

## LIFTING TRAY

CONFIGURATION EXAMPLES

**Single Endoscope Configuration**  
Holds: 1 da Vinci Xi or SP Endoscope



**Multiple Instrument Configuration**  
Holds: up to 8 da Vinci Xi or SP EndoWrist Instruments

# Heavy Duty Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. They are similar to our other transport containers but have a more robust construction and feature stainless steel hinges and latches that can be easily undone with just one hand.



## HEAVY DUTY TRANSPORT CONTAINER

<i>Short Container</i>	
20 ½ L x 12 ⅝ W x 3 ¾ H (in)	521 L x 321 W x 95 H (mm)
IN-PIC-15	



## HEAVY DUTY TRANSPORT CONTAINER

<i>Tall Container</i>	
20 L x 12 ¼ W x 6 ⅞ H (in)	510 L x 310 W x 156 H (mm)
IN-PIC-16	

### SPECIAL FEATURES

- *Autoclavable*
- *Made of PEI material (high-temperature plastic)*
- *High durability against automated washing*
- *Transparent amber color for easy instrument identification*
- *Easy locking, removable cover with removable silicone gasket that can be replaced*
- *Lightweight*
- *Combine with Instru-Zyme for optimal pre-cleaning*

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

# Transport Containers

These containers are designed with a silicone gasket and locking cover to reduce spills and leakage while transporting contaminated instruments. Also, you can easily identify content through the transparent amber color. Pair this product with our Instru-Zyme Pre-Cleaning Enzyme Gel for best results.

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES



TRANSPORT CONTAINER

*Small, Short Container*

13 5/16 L x 15 5/8 W x 4 5/16 H (in)	345 L x 398 W x 110 H (mm)
--------------------------------------	----------------------------

IN-PIC-14



TRANSPORT CONTAINER

*Large, Short Container*

20 5/16 L x 15 5/8 W x 5 5/16 H (in)	524 L x 397 W x 135 H (mm)
--------------------------------------	----------------------------

IN-PIC-12



TRANSPORT CONTAINER

*Small, Tall Container*

13 5/16 L x 15 5/8 W x 7 1/2 H (in)	345 L x 398 W x 190 H (mm)
-------------------------------------	----------------------------

IN-PIC-13



TRANSPORT CONTAINER

*Large, Tall Container*

20 5/16 L x 15 5/8 W x 7 1/2 H (in)	524 L x 397 W x 191 H (mm)
-------------------------------------	----------------------------

IN-PIC-11



# Silicone Small Parts Boxes

Don't sweat the small stuff. Made of silicone, these boxes are perfect for keeping extra delicate small parts and accessories organized and protected.



## SILICONE SMALL PARTS BOX


<i>3x3" Small Parts Box</i>	
3 L x 3 W x ¾ H (in)	76 L x 76 W x 19 H (mm)
IN-2652	



## SILICONE SMALL PARTS BOX

<i>3x6" Small Parts Box</i>	
6 L x 3 W x ¾ H (in)	152 L x 76 W x 19 H (mm)
IN-2653	

- Made entirely of medical grade silicone
- Maximum protection for small parts
- Tethered cover to prevent accidental loss
- Can be attached to trays or placed on posts for easy removal
- Textured internal surface and proper ventilation for efficient sterilant penetration



*Small Parts Box is shown within tray. Tray is shown without cover for illustrative purposes only.*

ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS

SCOPES

TRANSPORT CONTAINERS

ACCESSORIES

# Plaque Options

A variety of plaque colors are available to help identify specific sets by procedure or service type. Affix them to our trays and modules using our standard screws. Common identification uses are for bar code placement, tray name labels, or other custom information. Laser engraving of plaques is provided at an additional cost.

*Our in-house laser engraving capabilities allow us to personalize trays with custom identification and/or specifications, such as interior instrument descriptions, illustrations, exterior procedure types, specialty, or facility information.*



## PLAQUES

### Plaque Assembly for Trays

1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-01	SILVER
1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-05	BLACK
1 L x 4 3/8 W (in) 25 L x 111 (mm)	IN-2550-15	BLACK
Please indicate desired color of plaque using corresponding number below in place of XX using the color options table.		
4 3/8 L x 2 W (in) 111 L x 51 W (mm)	IN-2600-XX	●●●
2 L x 2 1/8 W (in) 51 L x 56 W (mm)	IN-2300-XX	●●●

## COLOR OPTIONS ●●●


### Steam Sterilization

These anodized aluminum plaques for use with steam only.

### Low-Temp Sterilization

Plaques have a coating that specially withstands low-temperature sterilization.

01	SILVER		N/A	
02	RED	■	RED	12
03	GREEN	■	GREEN	13
04	BLUE	■	BLUE	14
05	BLACK	■	BLACK	15
07	PURPLE	■	PURPLE	17
08	ORANGE	■	N/A	
10	MINT GREEN	■	N/A	



ARTHRO

BARIATRIC

CARDIO

DENTAL

ENT

GENERAL

GYN

NEURO

EYE

ORTHO

PLASTICS

RETRACTORS

ROBOTICS


SCOPES

TRANSPORT CONTAINERS

ACCESSORIES


# Finger Mats

These textured silicone instrument cushion mats are available in a variety of sizes and textures to fit all of your instrument protection needs.




## FINE FINGER MAT

<i>Full Sheet</i>	
19 L x 13 W (in)	483 L x 330 W (mm)
IN-2616	



## FINE FINGER MAT

<i>Half Sheet</i>	
9 ½ L x 13 W (in)	241 L x 330 W (mm)
IN-2617	



## FINE FINGER MAT

<i>Quarter Sheet</i>	
4 ¾ L x 13 W (in)	121 L x 330 W (mm)
IN-2618	

# Guardian™ Tray

Simple, small instrument sets now have a new home. Our silicone finger mats with sidewalls provide a safe and pliable environment for them to reside in.



## GUARDIAN TRAY

<i>Anodized Aluminum Tray with Finger Mat</i>	
10 L x 6 W x 1 ¼ H (in)	254 L x 152 W x 32 H (mm)
IN-1310	

- Accommodate a wide variety of small instrumentation
- Feature built-in and vented silicone finger mats
- Reduce instrument repair and replacement costs
- Allow maximum sterilant penetration
- Prevent moisture accumulation

- ARTHRO
- BARIATRIC
- CARDIO
- DENTAL
- ENT
- GENERAL
- GYN
- NEURO
- EYE
- ORTHO
- PLASTICS
- RETRACTORS
- ROBOTICS
- SCOPES
- TRANSPORT CONTAINERS
- ACCESSORIES



# Instru-Zyme®

This powerful enzyme and detergent formulation is designed to assist in the removal of blood, tissue, mucus and other proteinaceous materials. Use this product with our transport containers for a simple solution to water and detergent.



## INSTRU-ZYME

Properties:

- Specific Gravity: 1.200
- pH: 7.8
- Tropical Fragrance
- Color: Light Blue
- Form: Liquid Gel


6 Bottles/Case

IN-1032

- Will not stain instruments or equipment made of: metal, plastic, glass or rubber
- Non-corrosive
- Biodegradable
- Sustains moisture during transport
- Provides a ready-to-use, neutral pH enzyme solution
- Reduces heavy instrument transport weight as compared to liquid-filled trays
- Minimizes potential exposure due to splashing
- Spray applicator for even coverage of soiled areas
- No rinsing necessary before automated processing

# Rounders® Point Guards

Wrap your trays with confidence. Point Guards are engineered to sit securely on almost any tray corner. Prevent holes and tears in your sterilization wrap caused by surgical trays that are heavy, have sharp corners or are improperly handled.



## ROUNDERS POINT GUARDS


80/Bag

IN-7330

- Air cushion tray corner protection
- Help avoid compromised sterility in the O.R.
- Can be used in steam, ethylene oxide and low temperature plasma sterilization
- Reusable and economical
- Patent #9119611
- Non-sterile

**MADE IN THE USA**





**0**  
INCH

**1**

**2**

**3**

**4**

**5**

**6**

**7**


**8**

**9**

**10**

**11**





**InstruSafe**

Sterile Processing

**Summit Medical LLC**

**An Innovia Medical Company**

815 Vikings Parkway

Suite 100

St. Paul, MN 55121 | USA

1.888.229.2875

[www.innoviamedical.com/instrusafe](http://www.innoviamedical.com/instrusafe)

---

InstruSafe Instrument Protection Trays are proudly made in St. Paul, Minnesota, USA.

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and CMDCAS certified.

InstruSafe Trays are 510(k) cleared for multiple sterilization cycles to give our customers more options.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

da Vinci, SP, S, Xi, Single-Site and EndoWrist are registered trademarks of Intuitive Surgical, Inc. InstruSafe, Instru-Zyme, Hold-Its, Hold-Downs, Rounders Point Guards, and Guardian are trademarks and registered trademarks of Summit Medical.

The Bookwalter trademark is a registered trademark of Symmetry Surgical, Inc., or its affiliates in the United States and other countries around the world. Original Bookwalter retractors are sold exclusively by Symmetry Surgical, Inc. The Mediflex Bookler Retractor is registered trademark of Mediflex.

All brands are trademarks or registered trademarks of their respective owners.

©2022 Innovia Medical. All rights reserved. M-5 | RevC | JULY2022