

Innovia Medical

OPHTHALMOLOGY

Product Guide 2022/2023.

innoviamedical.com/

OPHTHALMOLOGY



Innovia Medical

OPHTHALMOLOGY

A close-up photograph of a person's hands wearing white surgical gloves, holding a small, cylindrical surgical instrument. The background is a blurred operating room environment with blue drapes and white equipment.

ELEVATING CARE. IMPROVING OUTCOMES.

As part of Innovia Medical's effort to be a more sustainable and environmentally friendly company, we have begun our journey to becoming net zero by 2040. We believe in protecting the planet by reducing our carbon footprint.

innoviamedical.com/net-zero

**NET
ZERO
ROADMAP**

About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

Contents


 Click on any topic below to go directly to that page. The [BACK TO CONTENTS](#) button on the bottom of each page will return you to this page.

0
0
0

CANNULA

EAGLE LABS[®]

The Eagle Lab[®] range of ophthalmic cannula offers a portfolio of products for a variety of procedures including:





CANNULA

EAGLE LABS®

A full portfolio of products for every stage of Cataract and Glaucoma surgical procedures. Trusted by hundreds of surgeons around the world for their 'Made in the USA' reliability and performance.

- Colour coded hub for easy identification
- Continuous radius/lumen to reduce risk of blockages
- Smooth precise tip
- Supplied sterile or non-sterile



Anaesthesia Administration

Tri-Port Sub Tenon Cannulae

10 per box



Code	Hub	Gauge	Size	Description
100-19		19G x 1"	1.0 x 25 mm	2 x .4mm Side Ports, Open Ended, Angled
100-19C		19G x 1"	1.0 x 25 mm	2 x .4mm Side Ports, Open Ended, Curved
100-21		21G x 1"	0.8 x 25 mm	2 x .4mm Side Ports, Open Ended, Angled
100-21C		21G x 1"	0.8 x 25 mm	2 x .4mm Side Ports, Open Ended, Curved

Atkinson Tip Retrobulbar


10 per box



Code	Hub	Gauge	Size	Description
101-23		23G x 1 1/2"	0.6 x 38 mm	
101-25		25G x 1 1/2"	0.5 x 38 mm	
101-27		27G x 1 1/2"	0.4 x 38 mm	

Atkinson Tip Peribulbar

10 per box



Code	Hub	Gauge	Size	Description
102-23-32		23G x 1 1/4"	0.6 x 32 mm	
102-25-32		25G x 1 1/4"	0.5 x 32 mm	
102-27-32		27G x 1 1/4"	0.4 x 32 mm	
102-25-22		25G x 7/8"	0.5 x 22 mm	
102-27-22		27G x 7/8"	0.4 x 22 mm	

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Sharp Point Retrobulbar **10 per box**



103-25

Code	Hub	Gauge	Size	Description
103-25		25G x 1 1/2"	0.5 x 38 mm	Sharp Tip Retrobulbar

Anterior Chamber Maintenance

Anterior Chamber Maintainer **3 per box**



170-20

Code	Hub	Gauge	Size	Description
170-20		20G x 3/16"	0.9 x 4.76 mm	Infusion Cannulae with Self Retaining Serrated Tip, 8" Flexible Tubing

Irrigation, Air Injection & Viscoelastic

Anterior Chamber Washout **10 per box**



111-16




111-19

Code	Hub	Gauge	Size	Description
111-16		16G x 1"	1.4 x 25 mm	25mm Long, Flattened Tip, Angled
111-19		19G x 1"	1.1 x 25 mm	25mm Long, Flattened Tip, Angled

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Bishop Harmon Cannula

10 per box



111-16BH



111-19BH

Code	Hub	Gauge	Size	Description
111-16BH		16G x 1"	1.4 x 25 mm	6mm Bend to Tip
111-19BH		19G x 1"	1.1 x 25 mm	8mm Bend to Tip

Anterior Chamber Cannula (Rycroft) (Straight, Round Tip)

10 per box



112-25



112-30

Code	Hub	Gauge	Size	Description
112-25		25G x 1"	0.5 x 25 mm	25mm Straight, Round Tip
112-27		27G x 1"	0.4 x 25 mm	25mm Straight, Round Tip
112-30		30G x 1"	0.3 x 25 mm	25mm Straight, Round Tip

Anterior Chamber Cannula (Rycroft) (4mm Bend to Tip)

10 per box



110-25-4




110-30-4

Code	Hub	Gauge	Size	Description
110-25-4		25G x 1"	0.5 x 25 mm	4mm Bend to Tip, Angled 30°
110-27-4		27G x 1"	0.4 x 25 mm	4mm Bend to Tip, Angled 30°
110-30-4		30G x 1"	0.3 x 25 mm	4mm Bend to Tip, Angled 30°


Hub Gauge

Anterior Chamber Cannula (Rycroft) (7mm Bend to Tip) 10 per box





Code	Hub	Gauge	Size	Description
110-20		20G x 1"	0.9 x 25 mm	7mm Bend to Tip, Angled 30°
110-21		21G x 1"	0.8 x 25 mm	7mm Bend to Tip, Angled 30°
110-22		22G x 1"	0.7 x 25 mm	7mm Bend to Tip, Angled 30°
110-23		23G x 1"	0.6 x 25 mm	7mm Bend to Tip, Angled 30°
110-25		25G x 1"	0.5 x 25 mm	7mm Bend to Tip, Angled 30°
110-26		26G x 1"	0.45 x 25 mm	7mm Bend to Tip, Angled 30°
110-27		27G x 1"	0.4 x 25 mm	7mm Bend to Tip, Angled 30°
110-30		30G x 1"	0.3 x 25 mm	7mm Bend to Tip, Angled 30°

Anterior Chamber Cannula (Rycroft) (9mm Bend to Tip) 10 per box



Code	Hub	Gauge	Size	Description
110-27-9		27G x 1"	0.4 x 25 mm	9mm Bend to Tip, Angled 30°
110-30-9		30G x 1"	0.3 mm	9mm Bend to Tip, Angled 30°

Anterior Chamber Cannula (Rycroft) (11mm Bend to Tip) 10 per box



Code	Hub	Gauge	Size	Description
110-25-11		25G x 1"	0.5 x 25 mm	11mm Bend to Tip, Angled 30°
110-27-11		27G x 1"	0.4 x 25 mm	11mm Bend to Tip, Angled 30°

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Anterior Chamber Cannula (Mendez) 10 per box



110-27TM

Code	Hub	Gauge	Size	Description
110-27TM		27G x 3/8"	0.4 x 10mm	3mm Bend to Tip, Angled 45°


Anterior Chamber Cannula (Rycroft) (8mm Bend to Tip) 10 per box



110-26 1/2

Code	Hub	Gauge	Size	Description
110-26 1/2		26G x 1/2"	0.45 x 12mm	8mm Bend to Tip, Angled 25°

Anterior Chamber Cannula (Rycroft) (12mm Bend to Tip) 10 per box



110-30SP

Code	Hub	Gauge	Size	Description
110-30SP		30G x 9/16"	0.3 x 14 mm	12mm Bend to Tip, Angled 45°, Beveled

FLOvisc Viscoelastic Cannula 10 per box



109-27-25



109-22

Code	Hub	Gauge	Size	Description
109-27-25		20G x 1"	0.9/0.5 mm x 25mm	20G Stepped Down to 25G, 10mm Angled Tip
109-27		20G x 1"	0.9/0.4 mm x 25mm	20G Stepped Down to 27G, 10mm Angled Tip
109-22		27G x 1"	0.7/0.4 mm x 25mm	27G, 10mm Angled Tip

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

Infusion (Wright) 10 per box




Code	Hub	Gauge	Size	Description
170-20W		20G x 1"	0.9 x 25 mm	7mm Bend to Beveled Tip, Angled 45°

I and A Cannulae (Simcoe) 6 per box



Code	Hub	Gauge	Size	Description
150-23		23G	0.6 x 330 mm	Irrigation through Cannulae, 0.4mm Side Port, Aspiration through 13" (33cm) Flexible Tubing

Reverse I and A Cannulae (Simcoe) 6 per box



Code	Hub	Gauge	Size	Description
150-23R		21G	0.8 x 330 mm	Aspiration through Cannulae, 0.4mm Side Port, Irrigation through 13" (33cm) Flexible Tubing

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

Capsulotomy

Irrigating Cystotome (Straight Tip)

10 per box




120-23S

Code	Hub	Gauge	Size	Description
120-23S		23G x 1"	0.6 x 25 mm	Straight Tip
120-25S		25G x 5/8"	0.5 x 16 mm	Straight Tip

Irrigating Cystotome (Straight, Reverse Tip)

10 per box



122-25S

Code	Hub	Gauge	Size	Description
122-23S		23G x 3/4"	0.6 x 19 mm	Straight, Reverse Tip
122-25S		25G x 5/8"	0.5 x 16 mm	Straight, Reverse Tip

Irrigating Cystotome (8mm Formed Tip)

10 per box




120-27F-8

Code	Hub	Gauge	Size	Description
120-25F-8		25G x 5/8"	0.5 x 16 mm	8mm Formed Tip
120-27F-8		27G x 1/2"	0.4 x 13 mm	8mm Formed Tip

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Irrigating Cystotome (8mm Formed Tip, Push/Pull) 10 per box



121-25F-8

Code	Hub	Gauge	Size	Description
121-25F-8		25G x 5/8"	0.5 x 16 mm	8mm Formed Tip, Push/Pull
121-27F-8		27G x 1/2"	0.4 x 13 mm	8mm Formed Tip, Push/Pull


Irrigating Cystotome (10mm Angled Tip) 10 per box



123-27

Code	Hub	Gauge	Size	Description
123-25		25G x 5/8"	0.5 x 16 mm	10mm Angled Tip
123-27		27G x 1/2"	0.4 x 13 mm	10mm Angled Tip


Irrigating Cystotome (12mm Formed Tip) 10 per box



120-25F-12

Code	Hub	Gauge	Size	Description
120-25F-12		25G x 5/8"	0.5 x 16 mm	12mm Formed Tip
120-27F-12		27G x 1/2"	0.4 x 13 mm	12mm Formed Tip

Irrigating Cystotome (12mm Formed Tip, Push/Pull) 10 per box



121-25F-12

Code	Hub	Gauge	Size	Description
121-25F-12		25G x 5/8"	0.5 x 16 mm	12mm Formed Tip, Push/Pull

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----


Irrigating Cystotome (12mm Formed Tip, Reverse) 10 per box



122-27RF

Code	Hub	Gauge	Size	Description
122-25RF		25G x 5/8"	0.5 x 16 mm	12mm Formed Tip, Reverse
122-27RF		27G x 1/2"	0.4 x 13 mm	12mm Formed Tip, Reverse

Irrigating Cystotome (13mm Formed Tip) 10 per box




120-01

Code	Hub	Gauge	Size	Description
120-01		25G x 3/4"	0.5 x 19 mm	13mm Formed Tip

Hydrodissection

Hydrodissection (Chang) 10 per box



113-01

Code	Hub	Gauge	Size	Description
113-01		27G x 1"	0.4 x 25 mm	1.5mm Flat Beveled Tip, 90° Angled

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G


Nucleus Hydrodissection (7mm Flattened Tip) 10 per box



113-23

Code	Hub	Gauge	Size	Description
113-23		23G x 1"	0.6 x 25 mm	7mm Flattened Tip
113-25		25G x 1"	0.5 x 25 mm	7mm Flattened Tip
113-27		27G x 1"	0.4 x 25 mm	7mm Flattened Tip
113-30		30G x 1"	0.3 x 25 mm	7mm Flattened Tip

Nucleus Hydrodissection (7mm Flattened Tip, Vertical Flow) 10 per box



213-25

Code	Hub	Gauge	Size	Description
213-25		25G x 1"	0.5 x 25 mm	7mm Flattened Tip, Vertical Flow
213-27		27G x 1"	0.4 x 25 mm	7mm Flattened Tip, Vertical Flow

Nucleus Hydrodissection (8 & 10mm Flattened Tip) 10 per box




113-27-10

Code	Hub	Gauge	Size	Description
113-25-8		25G x 1"	0.5 x 25mm	8mm Flattened Tip
113-25-10		25G x 1"	0.5 x 25mm	10mm Flattened Tip
113-27-10		27G x 1"	0.4 x 25mm	10mm Flattened Tip

Hub Gauge


	19G	20G	21G	22G	23G	24	25G	26G	27G	
--	-----	-----	-----	-----	-----	----	-----	-----	-----	--

Nucleus Hydrodissection (12mm Formed, Flattened Tip) 10 per box




Code	Hub	Gauge	Size	Description
113-25F-12		25G x 1"	0.5 x 25mm	12mm Formed, Flattened Tip
113-27F-12		27G x 1"	0.4 x 25mm	12mm Formed, Flattened Tip

Nucleus Hydrodissection (Curved, Horizontal / Vertical Flattened Tip) 10 per box



Code	Hub	Gauge	Size	Description
118-25		25G x 1"	0.5 x 25 mm	Curved, Horizontal Flattened Tip
118-25S		25G x 1"	0.5 x 25 mm	Curved, Vertical Flattened Tip
118-27		27G x 1"	0.4 x 25 mm	Curved, Horizontal Flattened Tip
118-27S		27G x 1"	0.4 x 25 mm	Curved, Vertical Flattened Tip

Cortical Cleaving Hydrodissection (Grabow) (Round Tip) 10 per box



Code	Hub	Gauge	Size	Description
217-27		27G x 1"	0.4 x 25 mm	11mm Bend to Round Tip, Angled 45°


Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Cortical Cleaving Hydrodissection (Flattened Tip) 10 per box



Code	Hub	Gauge	Size	Description
218-27		27G x 1"	0.4 x 25 mm	11mm Flattened Tip, Angled 45°

Hydrodissection (Tuebingen) 10 per box



Code	Hub	Gauge	Size	Description
313-30		30G x 1"	0.3 x 25 mm	75° Curved, Round Tip

Hydrodelineation

Nucleus Hydrolineation 10 per box



Code	Hub	Gauge	Size	Description
114-30		30G x 1"	0.3 x 25 mm	Curved, Beveled Tip

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

Nucleus Removal

Irrigating Vectis

10 per box



140-25-1

Code	Hub	Gauge	Size	Description
140-25-1		25G	0.5 mm	5mm Wide (1 ea.), 0.3mm Open Port
140-25-2		25G	0.5 mm	5mm Wide (2 ea.), 0.3mm Open Port
140-25-3		25G	0.5 mm	5mm Wide (3 ea.), 0.3mm Open Port
140-27-1		27G	0.4 mm	3mm Wide (1 ea.), 0.3mm Open Port

Cortex Removal

Simcoe Cortex Extractor

10 per box



108-23

Code	Hub	Gauge	Size	Description
108-23		23G x 1"	0.6 x 25 mm	Curved, 0.3mm Inside Port

Aspirating J-Hook

10 per box



105-27S

Code	Hub	Gauge	Size	Description
105-27S		27G x 1"	0.4 x 25 mm	2.0mm x 2.0mm Straight Hook

Hub Gauge



Aspirating Cortex Extractor **10 per box**



106-27S

Code	Hub	Gauge	Size	Description
106-27S		27G x 1"	0.4 x 25 mm	1.5mm x 1.5mm Straight Hook


Cortex Aspiration (Male Adapter, Curved) **10 per box**



191-23

Code	Hub	Gauge	Size	Description
191-23		23G x 1"	0.6 x 25 mm	0.4mm Side Port, Male Adapter, Curved

Bimanual I/A Cortex Removal System **10 per box**



192-23

Code	Hub	Gauge	Size	Description
192-23		21G/23G x 1"	0.8/0.6 mm	0.4mm Side Port, Male/Female Adapters, Curved

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

Capsule Polishing

Capsule Polishers (Lu)

10 per box




124-23

Code	Hub	Gauge	Size	Description
124-23		23G x 1"	0.6 x 25 mm	30mm Top Port, Rounded Tip, Curved, Roughened Olive Tip, Open Port

Capsule Polishers (Terry)

10 per box




125-27

Code	Hub	Gauge	Size	Description
125-23		23G x 1"	0.6 x 25 mm	7mm Angled, Silicone Tip
125-27		27G x 1"	0.4 x 25 mm	7mm Angled, Silicone Tip

Capsule Polishers (Kratz)

10 per box




125-24

Code	Hub	Gauge	Size	Description
125-24		24G x 1"	0.55 x 25 mm	8mm Angled, 0.3mm Port

Hub Gauge



Capsule Polishers (10mm Angled, Olive Tip) 10 per box



Code	Hub	Gauge	Size	Description
126-23		23G x 1"	0.6 x 25 mm	10mm Angled, Olive Tip
126-25		25G x 1"	0.5 x 25 mm	10mm Angled, Olive Tip


Silicone Tip Capsule Polishers (Engels) 10 per box



Code	Hub	Gauge	Size	Description
128-27		27G x 1"	0.4 x 25 mm	40° Curve, Extended Silicone Sleeve

Corneal Abrasion

Stromal Puncture Needle 10 per box



Code	Hub	Gauge	Size	Description
165-25		25G x 1"	0.5 x 25 mm	19mm Bend to Tip, 30° Angled Sharp Tip

Hub Gauge 19G 20G 21G 22G 23G 24 25G 26G 27G 30G

Fluid Administration

Ocular Irrigation Tube

20 per box



169-30L-6

Code	Hub	Gauge	Size	Description
169-30L-6			6"	Extension Tube with Male/Female Luer Connections
169-30L-10			10"	Extension Tube with Male/Female Luer Connections

Glaucoma

Visco Canalostomy

5 per box



177-33

Code	Hub	Gauge	Size	Description
177-33		32G		5mm Beveled Tip, Angled

Glaucoma Revision Pick (Lederer)

10 per box



178-27

Code	Hub	Gauge	Size	Description
178-27		27G x 1/2"	0.4 x 13 mm	10mm Bend to Tip, Angled 30°

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

Oculoplastic

Lacrimal Cannula (Curved, Open End)


10 per box



Code	Hub	Gauge	Size	Description
160-21		21G x 1 1/2"	0.8 x 38 mm	Curved, Open End
160-25		25G x 1 1/2"	0.5 x 38 mm	Curved, Open End
160-26		26G x 1 1/2"	0.45 x 38 mm	Curved, Open End

Lacrimal Cannula (Curved, 0.3mm Side Port)



10 per box



Code	Hub	Gauge	Size	Description
161-23		23G x 1 1/2"	0.6 x 38 mm	Curved, 0.3mm Side Port
161-25		25G x 1 1/2"	0.5 x 38 mm	Curved, 0.3mm Side Port

Lacrimal Cannula (Straight, 0.3mm Side Port)

10 per box



Code	Hub	Gauge	Size	Description
163-21		21G x 1 1/2"	0.8 x 38 mm	Straight, 0.3mm Side Port
163-25		25G x 1 1/2"	0.5 x 38 mm	Straight, 0.3mm Side Port

Lacrimal Cannula (25mm Straight with 13mm 23G Tip)

10 per box



Code	Hub	Gauge	Size	Description
162-23		19G/23G x 1 1/2"	1.0 x 38 mm	25mm Straight with 13mm 23G Tip


Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	------------	------------	------------	------------	------------	-----------	------------	------------	------------	------------

Refractive (LASIK)

LASIK Irrigating Cannula (Guell)

10 per box




605-25

Code	Hub	Gauge	Size	Description
605-25		25G x 1"	0.5 x 25 mm	9mm Formed Bend to Tip, 6 Ports Open Ended

LASIK Irrigating Cannula (Banaji)

10 per box




604-25

Code	Hub	Gauge	Size	Description
604-25		25G x 1"	0.5 x 25 mm	12mm Formed, 6x ports, Closed End

LASIK Irrigating Cannula (Burrato)

10 per box




606-25

Code	Hub	Gauge	Size	Description
606-25		25G x 1"	0.5 x 25 mm	9mm Formed Bend, 2 Ports Open Ended

LASIK Irrigating Cannula (Lindstrom)

10 per box



603-25

Code	Hub	Gauge	Size	Description
603-25		25G x 1"	0.5 x 25 mm	7mm Angled, 3x 0.3mm Ports

Hub Gauge

	19G	20G	21G	22G	23G	24	25G	26G	27G	30G
--	-----	-----	-----	-----	-----	----	-----	-----	-----	-----

LASIK Flap Irrigator (Vidaaurri) 10 per box

Code	Hub	Gauge	Size	Description
600-16		25G x 1"	0.5 mm	9mm Angled Bend to Tip, Dual Multi-Port Irrigation Tips, 2 Ports (ea.)
600-06		25G x 1"	0.5 mm	13mm Angled Bend to Tip, Dual Multi-Port Irrigation Tips, 4 Ports (ea.)

LASIK Irrigating Cannula (Nichamin) 10 per box

Code	Hub	Gauge	Size	Description
600-25		25G x 1"	0.5 x 25 mm	12mm Formed, Flattened Tip, Open Irrigation
600-27		27G x 1"	0.4 x 25 mm	12mm Formed, Flattened Tip, Open Irrigation


LASIK Flap Spatula/Irrigator (Hersh) 10 per box

Code	Hub	Gauge	Size	Description
610-23		23G x 1"	0.6 x 25 mm	9mm Angled Bend to Tip
610-25		25G x 1"	0.5 x 25 mm	9mm Angled Bend to Tip

LASIK Irrigating Spatula (Slade) 10 per box


Code	Hub	Gauge	Size	Description
613-26		26G x 1"	0.45 x 25 mm	9mm Formed Bend to Tip


LASIK Irrigating Cannula (Soloman) 10 per box



Code	Hub	Gauge	Size	Description
614-25		25G x 1"	0.5 x 25 mm	Formed Spatula Tip, 3 Ports

LASIK Irrigating Cannula 10 per box



Code	Hub	Gauge	Size	Description
617-27		27G x 1"	0.4 x 25 mm	10mm Formed Bend to Tip, Flattened, 4 Ports

Vitreoretinal

Subretinal Fluid Cannula (Gizzard) 10 per box




Code	Hub	Gauge	Size	Description
166-25		20G	0.9 x 37 mm	Silicone Sleeve Extended 6mm Beyond Tip, Straight

KNIVES & BLADES

EAGLE LABS[®]

The Eagle Lab[®] range of ophthalmic knives & blades offers a portfolio of products for a variety of procedures including:





KNIVES & BLADES

EAGLE LABS®

A full portfolio of products for every stage of Cataract and Glaucoma surgical procedures. Trusted by hundreds of surgeons around the world for their 'Made in the USA' reliability and performance.

- Lightweight and ergonomic handle
- Supplied sterile, single use only, declared 5 year shelf life
- Precise micro surgical design
- Manufactured in the USA



Surgical Knives

Double Bevel Sideport Knife

6 per box



EG-5010

Code	Description
EG-5010	1.0mm, Angled, Full Handle, Double Bevel
EG-5012	1.2mm, Angled, Full Handle, Double Bevel

Spoon Knife

6 per box



EG-6821

Code	Description
EG-6821	3.0mm Angled, Full Handle, Double Bevel

MVR Knife

6 per box



EG-5555

Code	Description
EG-5555	19GA, Full Handle
EG-5560	20GA, Full Handle
EG-5563	23GA, Full Handle

Clear Corneal Knife **6 per box**



Code	Description
EG-1122	2.2mm, Angled, Full Handle, Double Bevel
EG-1124	2.4mm, Angled, Full Handle, Double Bevel
EG-1127	2.75mm, Angled, Full Handle, Double Bevel
EG-1130	3.0mm, Angled, Full Handle, Double Bevel
EG-1132	3.2mm, Angled, Full Handle, Double Bevel

Clear Corneal Knife, Parallel Sides **6 per box**



Code	Description
EG-1327	2.75mm, Angled, Full Handle, Double Bevel

Trapezoid Clear Corneal/ Implant Knife **6 per box**



Code	Description
EG-2022	2.0mm - 2.2mm, Angled, Full Handle, Bevel Up
EG-2832	2.85mm - 3.2mm, Angled, Full Handle, Bevel Up

Stab Knife **6 per box**



Code	Description
EG-1500	15°, Full Handle
EG-2250	22.5°, Full Handle
EG-3000	30°, Full Handle
EG-4500	45°, Full Handle

Phaco Slit Knife **6 per box**



Code	Description
EG-1022	2.2mm, Angled, Full Handle, Bevel Up
EG-1024	2.4mm, Angled, Full Handle, Bevel Up
EG-5510	2.5mm, Angled, Full Handle, Bevel Up
EG-5568	2.65mm, Angled, Full Handle, Bevel Up
EG-4738	2.75mm, Angled, Full Handle, Bevel Up
EG-4740	2.85mm, Angled, Full Handle, Bevel Up
EG-4741	2.85mm, Straight, Full Handle
EG-5570	3.0mm, Angled, Full Handle, Bevel Up
EG-5571	3.0mm, Straight, Full Handle
EG-5572	3.2mm, Angled, Full Handle, Bevel Up
EG-5573	3.2mm, Straight, Full Handle

IOL Inserter Knife **6 per box**



Code	Description
EG-9350	5.2mm, Angled, Full Handle, Bevel Up
EG-9355	5.5mm, Angled, Full Handle, Bevel Up

Crescent Knife **6 per box**



Code	Description
EG-9402	2.0mm, Angled, Full Handle, Bevel Up

Surgical Blades/Helpers

Needle Blade **6 per box**



Code	Description
EL-5910	3.0mm

Lamellar Blade **6 per box**



EL-6600

Code	Description
EL-6600	60°

Scleral Blade **6 per box**



EL-5700

Code
EL-5700

Miniature Edged Blade **20 per box**



EL-6400



EL-6500



EL-6700

Code	Description
EL-6400	Round Tip
EL-6500	Lance Tip
EL-6700	Curved Tip

Miniature Edged Blade **6 per box**



EL-6900

Code	Description
EL-6900	Sharp All Around

Microsurgical Blade

6 per box



EL-7513



EL-7514

Code	Description
EL-7513	3.0mm, 15°, Blue
EL-7514	3.0mm, 15°, Blue, with Handle
EL-7515	5.0mm, 15°, Green
EL-7516	5.0mm, 15°, Green, with Handle

Precision LASIK Blades

LASIK Blade

10 per pack



ES-3030

Code	Compatible with	Description
ES-3030 +10	Hansatome™ microkeratome	Plus 10
ES-3030 -10	Hansatome™ microkeratome	Minus 10
ES-3030 -20	Hansatome™ microkeratome	Minus 20
ES-3030 -30	Hansatome™ microkeratome	Minus 30
ES-3030 -PLANO	Hansatome™ microkeratome	Plano


LASIK Blade


10 per pack



ES-5050

Code	Compatible with	Description
ES-5050 +10	Zyoptix™ XP microkeratome	Plus 10
ES-5050 +20	Zyoptix™ XP microkeratome	Plus 20
ES-5050 -10	Zyoptix™ XP microkeratome	Minus 10
ES-5050 -20	Zyoptix™ XP microkeratome	Minus 20
ES-5050 -30	Zyoptix™ XP microkeratome	Minus 30
ES-5050 -PLANO	Zyoptix™ XP microkeratome	Plano

LASIK Blade		10 per pack
 <p>ES-7070</p>		
Code	Compatible with	Description
ES-7070 +10	Moria™ M2 microkeratome	Plus 10
ES-7070 -10	Moria™ M2 microkeratome	Minus 10
ES-7070 -20	Moria™ M2 microkeratome	Minus 20
ES-7070 -30	Moria™ M2 microkeratome	Minus 30
ES-7070 -PLANO	Moria™ M2 microkeratome	Plano

LASIK Blade		10 per pack
 <p>ES-9090</p>		
Code	Compatible with	Description
ES-9090 +10	Nidek™ microkeratome	Plus 10
ES-9090 +20	Nidek™ microkeratome	Plus 20
ES-9090 -10	Nidek™ microkeratome	Minus 10
ES-9090 -20	Nidek™ microkeratome	Minus 20
ES-9090 -30	Nidek™ microkeratome	Minus 30
ES-9090 -PLANO	Nidek™ microkeratome	Plano

DCR & LACRIMAL EYETEC[®]

The Eyetec[®] range of ophthalmic DCR & Lacrimal offers a portfolio of products for a variety of procedures including:



DCR & LACRIMAL EYETEC®

Crawford-type Intubation Set:




- Probe length 150mm
- Silicone tube length 320mm
- Retrieval device 105mm length
- Olive tip

Canaliculus Intubation Set:

- Probe length 185mm
- Silicone tube length 300mm

DCR Bodkins Straight and Angled:

- Bodkin length 45mm
- Silicone rod length 480mm



Lacrimal Intubation

Crawford-type Intubation Set



51-918

Code	Description	Units of Sale
51-917	Retrieval device for Crawford-Type Intubation	Box of 5
51-918	Crawford-Type Intubation Set	Box of 1

Canaliculus Intubation Set



51-920

Code	Description	Units of Sale
51-920	Canaliculus Intubation Set	Box of 5

DCR Bodkins

DCR Bodkins Straight and Angled



51-922



51-923


Code	Description	Units of Sale
51-922	DCR Bodkins Straight	Box of 3
51-923	DCR Bodkins Angled	Box of 3

FLUIDS MANAGEMENT

EYETEC®

The Eyetec® brand is suitable for a variety of ophthalmic procedures including cataract, cornea graft and LASIK:

- + PVA and cellulose eye spreads
- + Light shields
- + Eye drains and wicks
- + Instrument wipes and cleaning blocks



FLUIDS MANAGEMENT

EYETEC[®]

PVA and Cellulose ophthalmic sponge products are designed to staunch blood loss and absorb fluids from and around the eye during ophthalmic procedures. The PVA and Cellulose sponge reacts quickly to absorb blood and body fluids.

- Constructed from ultra-smooth micropore PVA sponge
- Ultra-fast wicking action
- Supplied sterile, single use only, declared 5-year shelf life
- Manufactured in the USA




EYETEC[®]

PVA & Cellulose Fluids Management.


KEY PRODUCT FEATURES

- ✓ Ultra-fast wicking
- ✓ Smooth and soft surface
- ✓ Small uniformed pores
- ✓ Less particulate/debris than standard PVA spears
- ✓ Decreased risk of tissue adhesion and post operative bleeding



OTHER PVA

- Rough and uneven surface
- Large irregular pores
- Increased particulate number
- Increased risk of tissue adhesion



EYETEC[®] PVA

- Smooth and soft surface
- Small uniformed pores
- Lower particulate number
- Decreased risk of tissue adhesion and post operative bleeding



PVA Spear vs Standard Spear


Eye Spears and Points

EYETEC® PVA Spears



Code	Description	Units of Sale
40-400	PVA Eye Spears	Pack of 5 box of 90
40-405	PVA Eye Spears	Pack of 5 box of 180
40-406	PVA Points	Pack of 10 box of 180

EYETEC® Cellulose Spears



Code	Description	Units of Sale
40-410	Cellulose Spears	Pack of 5 box of 100
40-415	Cellulose Spears	Pack of 10 box of 200

Corneal Shields

EYETEC® Light Shields



Code	Description	Units of Sale
40-420	PVA Corneal Light Shields 7mm Diameter	Pack of 1 box of 20
40-421	PVA Corneal Light Shields 8mm Diameter	Pack of 1 box of 20
40-423	PVA Corneal Light Shields 10mm Diameter	Pack of 1 box of 20

Eye Shields

EYETEC® LASIK Eye Shields



Code	Description	Units of Sale
40-820	PVA LASIK Eye Shields 9mm Diameter/7mm Edge	Pack of 1 box of 20
40-821	PVA LASIK Eye Shields 8mm Diameter/6mm Edge	Pack of 1 box of 20
40-822	PVA LASIK Eye Shields 4mm Edge/8mm Edge	Pack of 1 box of 20

Eye Drains

EYETEC® Eye Drains




Code	Description	Units of Sale
40-430	PVA Eye Drains 80cc	
40-435	PVA Eye Drains 400cc	

EYETEC® LASIK Eye Drains



Code	Description	Units of Sale
40-830	PVA LASIK Eye Drains 42 x 17 x 11.5mm Opening	Pack of 1 box of 20
40-831	PVA LASIK Eye Drains 40 x 14 x 10.5mm Opening	Pack of 1 box of 20


EYETEC® Eye Wicks



Code	Description	Units of Sale
40-431	Eye Wick	Pack of 2 box of 40
40-470	Wick and Wipe	Pack of 1 box of 20

Instrument Wipes

EYETEC® Instrument Wipe



Code	Description	Units of Sale
40-900	Instrument Wipe	Pack of 1 box of 20

Diamond Cleaning Block

EYETEC® Diamond Cleaning Block



Code	Description	Units of Sale
40-462	Diamond Cleaning Block	Pack of 1 box of 20

PUNCTAL PLUGS

LACRIMEDICS®

Lacrimedics' products for Punctal Occlusion help evaluate, diagnose and treat DED and the Dry Eye component of varying Ocular Surface Diseases.

Absorbable intracanalicular plugs are used for short- and medium-term Occlusion Therapy, and non-absorbable Plugs are used for long-term Occlusion Therapy.:






PUNCTAL PLUGS

LACRIMEDICS®

Punctal Occlusion improves the ocular environment and can help keep patients wearing their Contact Lenses comfortably.

Punctal Occlusion addresses one of the main causes of Dry Eye Disease and may prevent the development of related health conditions.

- Special packaging maintains plug integrity
- Lacrimedics' uses Polydioxanone specifically developed for Punctal Occlusion.
- Unique colour makes them easier to locate
- Dissolvable or Non-dissolvable types
- Specially designed to regulate tear drainage




LACRIMEDICS®

Medium & long term occlusion therapy.

KEY FEATURES

- ✓ Medium-term Occlusion Therapy, lasting for approximately 180 days
- ✓ Long-term Occlusion Therapy for upper or lower Puncta
- ✓ Specially designed to regulate Tear Drainage - Tears to flow through the Punctal opening
- ✓ Equipped with a collapsible tip to ensure ease of insertion, improved patient comfort, and superior retention



DISSOLVABLE VISIPLUG®

Medium-term Occlusion Therapy



LONG TERM COMFORTTIP™

Long-term Occlusion Therapy



Dissolvable Lacrimal Plugs


VisiPlug® Dissolvable Lacrimal Plugs



Code	Description	Pack Configuration	Box Size
1638	Dissolvable VisiPlug® Lacrimal Plugs, 0.4mm	2 Plugs per Pouch with 1 Pouch per Box	2 Units
1639	Dissolvable VisiPlug® Lacrimal Plugs, 0.5mm	2 Plugs per Pouch with 1 Pouch per Box	2 Units


Punctal Occluder

ComfortTip™ Punctal Occluder



Code	Description	Pack Configuration	Box Size
2003	ComfortTip™ Punctal Occluder, X-Small	1 Plug per Tray with 2 Trays per Box	2 Units
2005	ComfortTip™ Punctal Occluder, Small	1 Plug per Tray with 2 Trays per Box	2 Units
2008	ComfortTip™ Punctal Occluder, Medium	1 Plug per Tray with 2 Trays per Box	2 Units
2010	ComfortTip™ Punctal Occluder, Large	1 Plug per Tray with 2 Trays per Box	2 Units

AccuFlo™ Punctal Occluder




Code	Description	Pack Configuration	Box Size
3000	AccuFlo™ Punctal Occluder, Small	1 Plug per Tray with 2 Trays per Box	2 Units
3001	AccuFlo™ Punctal Occluder, Medium	1 Plug per Tray with 2 Trays per Box	2 Units
3002	AccuFlo™ Punctal Occluder, Large	1 Plug per Tray with 2 Trays per Box	2 Units
3003	AccuFlo™ Punctal Occluder, X-Large	1 Plug per Tray with 2 Trays per Box	2 Units

ORBITAL IMPLANTS

CORONET[®]

The Eagle Lab[®] range of ophthalmic knives & blades offers a portfolio of products for a variety of procedures including:




ORBITAL IMPLANTS

CORONET[®]

All CORONET[™] Hydroxyapatite and Alumina Orbital Implants have a 100% interconnected micropore structure, which provides an ideal surface area for absorption of biomolecules, which are essential for cellular attachment and infiltration.

- Ultra-high porosity (80%>)
- Constructed with 100% medical grade synthetic Hydroxyapatite (HA).
- Recommended for PEG procedures for oculoplastic prosthesis.
- Supplied sterile, single use only, one per pack, 5 year declared shelf life.




CORONET®

Hydroxyapatite & Alumina Orbital Implants


KEY FEATURES

- ✓ Fully interconnected, highly porous hydroxyapatite provides a strong yet lightweight implant
- ✓ High macroporosity provides a very lightweight design - over 30% lighter than coral based implants.
- ✓ A non-resorbable, fully synthetic and biocompatible composition provides an ideal scaffold for soft tissue ingrowth. High integral strength - easily drilled for a motility peg.



ALUMINA ORBITAL IMPLANTS

Medium-term Occlusion Therapy




HYDROXYAPATITE ORBITAL IMPLANTS

Medium-term Occlusion Therapy

Orbital Implant

CORONET® Hydroxyapatite Orbital Implant



Code	Implant Size
54-300-14	14 mm
54-300-16	16 mm
54-300-17	17 mm
54-300-18	18 mm
54-300-19	19 mm
54-300-20	20 mm
54-300-21	21 mm
54-300-22	22mm

CORONET® Alumina Orbital Implants



Code	Description
54-305-14	14 mm
54-305-16	16 mm
54-305-17	17 mm
54-305-18	18 mm
54-305-19	19 mm
54-305-20	20 mm
54-305-21	21 mm
54-305-22	22mm

CORNEA GRAFT CORONET®

Coronet® Trephines, Donor Punches, Recipient Trephines and the EndoGlide Cornea Transplant Cartridge system (to name a few), comprise the most comprehensive range of products found anywhere for cornea graft surgery including:



CORNEA GRAFT CORONET®

The most important part of any trephine is the blade. Exclusive CORONET® blade technology is designed to create a sharp trephine.

- Fully patented Cathedral Blade Technology
- Consistent and reliable dual bevel sharpness
- Achieves a vertical cut (donor button & recipient) with minimal tissue compression
- Vertically matched corneas (PK) for wound apposition and tectonic strength (donor button to recipient)
- Made in the UK




CORONET[®]

Cornea Graft Products & Technology.


KEY PRODUCT FEATURES

- ✓ Fully patented Blade Technology
- ✓ Consistent and reliable dual bevel sharpness
- ✓ Achieves a vertical cut (donor button & recipient) with minimal tissue compression
- ✓ Vertically matched corneas for wound apposition and tectonic strength (donor button to recipient)



COMPETITOR TREPHINE X 800 MAGNIFICATION

- Blade edge on the outside only
- Distorts the tissue and creates undercut
- Retains donor graft
- Compromises on quality



CORONET[®] TREPHINE X 800 MAGNIFICATION

- Double bevelled
- Improves wound architecture
- Provides reliable and consistent incisions
- Reduces endothelial cell loss

Vacuum Trephine

CORONET® Recipient Corneal Vacuum Trephine



Code	Sizes (mm)												
51-840*	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.50	

*Blade progression 0.063mm downward per quarter turn and 0.252mm per full rotation

CORONET® PK Clear-Vue™ Recipient Corneal Vacuum Trephine



Code	Sizes (mm)												
51-842*	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.50	

*Blade progression 0.063mm downward per quarter turn and 0.252mm per full rotation

CORONET® DALK Clear-Vue™ Recipient Vacuum Trephine



Code	Sizes (mm)			
51-844	7.50	8.00	8.50	

Trephine Punch

CORONET® Recipient Corneal Vacuum Trephine



Code	Sizes (mm)															
51-835	6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.25	9.50	9.75

CORONET® Donor Corneal Trephine Punch



Code	Sizes (mm)															
51-850	6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.25	9.50	9.75

CORONET® Donor Corneal Trephine Punch



Code	Sizes (mm)						
51-855	7.50	8.00	8.50	9.00	9.50	10.0	

Guarded Trephine

Guarded Corneal Trephines | 300µm



Code	Sizes (mm)																
51-860					7.00		7.50		8.00		8.50						

Guarded Corneal Trephines | 400µm



Code	Sizes (mm)																
51-865					7.00		7.50		8.00		8.50						

Standard Trephine

CORONET® Standard Corneal Trephine



Code	Sizes (mm)																	
51-900	6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.25	9.50	9.75	10.0	

CORONET® Standard Corneal Trephine with Cross Hair



Code	Sizes (mm)																	
51-901	6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.25	9.50	9.75	10.0	

Long Handed Trephine

CORONET® Long Handled Corneal Trephine



Code	Sizes (mm)																	
51-903*	6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75	8.00	8.25	8.50	8.75	9.00	9.20	9.25	9.50	9.75	10.0

*Please replace the asterix with your desired size e.g. "51-903-6.00" for 6mm, or "51-903-7.00" for 7mm

Exactor Trephine

CORONET® Exactor Trephine



Code	Sizes (mm)			
51-905*	15.0	17.0	18.0	20.0

*Please replace the asterix with your desired size e.g. "51-903-6.00" for 6mm, or "51-903-7.00" for 7mm

EndoGlide™

CORONET® EndoGlide™ Ultrathin



Code	Description
51-823	4.5mm Sceral Tunnel Incision, Double Coil Technique

CORONET® DMEK EndoGlide™



Code	Description
51-826	2.65mm Clear Corneal Incision, Tri-Fold Technique

EndoGlide™ Forceps

CORONET® Tan EndoGlide™ Loading Forceps



Code	Description
53-951	Tissue Placement, 23G
53-955	Tissue Placement, 27G

CORONET® Tan EndoGlide™ Loading Forceps



Code	Description
53-952	Tissue Loading, 23G
53-956	Tissue Loading, 27G

Artificial Anterior Chamber

CORONET® Artificial Anterior Chamber

1 per box



Code	Description
51-935	5 Year Shelf Life, Single Use, Sterile

*51-935 is sold sterile, single use with a 5 year shelf life

CORONET® Surgical Trainer

1 per box



Code	Description
51-3000	Surgical Trainer, non-sterile
51-3010	Surgical Trainer Refill Kit (10 per box)

SINGLE-USE INSTRUMENTS

DTR MEDICAL®

The DTR Medical range of ophthalmic Single-Use instruments offers a portfolio of products for a variety of procedures including:

- + Forceps
- + Bipolar Forceps
- + Scissors
- + Castroviejo Needle Holder
- + Barraquer Needle Holder
- + Negus Aspirating Dissector
- + Procedure Packs
- + Double Ended Punctum Dilator
- + Seekers and Probes
- + Blade Handles
- + Kratz-Barraquer Eye Speculum
- + Retractors
- + Kerrison Rongeur




SINGLE-USE INSTRUMENTS

DTR MEDICAL®

Innovia Medical has a range of single-use ophthalmic instruments on offer. The DTR Medical brand is synonymous with reliability, manufacturing excellence, and dependability.

- Disposable surgical instruments for Theatres, Clinics & Primary Care
- Pre-cut silicone tubing and other consumables
- Specialist procedure trays



DTR MEDICAL[®]

Single-use Ophthalmic Instruments.

THE BENEFITS OF SINGLE-USE

- ✓ Matching reusable standards
- ✓ High quality instruments driven by sustainable material, process and design innovation
- ✓ “Supporting visibility, quality and safety in your supply chain”
- ✓ Packaging features the latest machine readable GS1-standard barcoded labels for product traceability
- ✓ Cost-effective alternative to repairs, reprocessing or replacements
- ✓ Suitable for high risk procedures



TIME

MATCHING REUSABLE STANDARDS

Save Time - Reprocessing

First time sharpness

Consistently clean



LIFE

SUITABLE FOR HIGH RISK PROCEDURES

Save Lives - Contamination

Cross-infection is eliminated

Risk is avoided



COST

COST EFFECTIVE

Save Money - Increase Productivity

Theatre and clinic time is reduced

Diminished re-processing, repairs and replacement

Forceps

Typing Forceps

20 Units



Birks Tying Forceps 110mm

Code	Description
BTF2115	Birks Tying Forceps 110mm, 5mm Platform



Straight Tying Forceps 90mm

Code	Description
STF2905	Straight Tying Forceps 90mm, 5mm Platform



Straight Tying Forceps 72mm

Code	Description
STF2725	Straight Tying Forceps 72mm, 5mm Tying Platform

Fine Suture Forceps

20 Units



Fine Suture Forceps Plain

Code	Description
SPF2001-S	Fine Suture Forceps Plain 105mm



Fine Suture Forceps Serrated

Code	Description
SSF2001-S	Fine Suture Forceps Serrated 105mm

Cyst Forceps **20 Units**



Code	Description
TCB2013	Birmingham Tarsal Cyst Forceps 84mm, 13mm Clamp





Code	Description
TCB2016	Birmingham Tarsal Cyst Forceps 84mm, 16mm Clamp



Code	Description
TCL2928	Lambert Tarsal Cyst Forceps 92mm, 8mm Clamp

Pierse Hoskins Forceps **10 Units**




Code	Description
PHF1085	Pierse Hoskins Forceps 85mm, 0.2mm Notched Tips

Crile Artery Forceps

20 Units



Code	Description
CFS2001-D	Crile Artery Forceps 145mm Straight




Code	Description
CFC2001-D	Crile Artery Forceps 145mm Curved on Flat


Adson Forceps




Adson Toothed Forceps



Adson Non-Toothed Forceps



Adson Forceps with Plateau



Fine Tip Adson Toothed Forceps

Adson Fine Tip Toothed Forceps 15mm




Code	Description	Units
ADS5001-S	Adson Non-Toothed Forceps 120mm	50
ATS5001-S	Adson Toothed Forceps 120mm	50
ATP2012-D	Adson Forceps 1x2 with plateau 120mm	20
ATP2015-D	Adson Forceps 1x2 with plateau 150mm	20
ATS2012F	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm	20

Epilation Forceps **20 Units**



Code	Description
EPF2085	Epilation Forceps 85mm, 2 x 3mm Plain Tips

St Martins Forceps **20 Units**



Code	Description
SMF2905	St. Martins Forceps 1:2 Toothed 90mm, 5mm Platform, 0.4mm Tips



Moorfield Forceps **20 Units**



Code	Description
MFF2113	Moorfield Forceps 105mm, 13mm Serration

Mosquito Forceps Curved

20 Units



Code	Description
MQC2012	Micro Mosquito Forceps Curved 120mm
MQC2025-D	Mosquito Forceps Curved 125mm
MQC2125-D	Mosquito Forceps Curved Toothed 1x2 125mm

Mosquito Forceps Straight



50 Units



Code	Description
MQT5001-S	Mosquito Forceps Straight 125mm

Capsulorhexis Forceps

20 Units




Code	Description
CXF2001	Capsulorhexis Forceps, Utrata Tips, 12mm Curved Shafts, Round Handles

Bipolar Forceps

Cawthorne Bipolar Forceps

20 Units



Code	Description
SU760508044	Cawthorne Bipolar Forceps 80 x 0.5mm, 3m, 4mm Plugs
SU761008044	Cawthorne Bipolar Forceps 80 x 1mm, 3m, 4mm Plugs
SU760511044	Cawthorne Bipolar Forceps 110 x 0.5mm, 3m, 4mm Plugs
SU761011044	Cawthorne Bipolar Forceps 110 x 1mm, 3m, 4mm Plugs

Bayonet Bipolar Forceps

20 Units




Code	Description
SU781020044	Bayonet Bipolar Forceps 200 x 1mm, 3m, 4mm Plugs
SU782020044	Bayonet Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs

Scissors

Westcott Scissors

20 Units



Code	Description
WSC2115	Westcott Scissors Curved 115mm, 15mm Blade

Stevens Tenotomy Scissors Straight


20 Units



Code	Description
TSS2001-S	Stevens Tenotomy Scissors 115mm Straight

Stevens Tenotomy Scissors Curved

20 Units



Code	Description
TSC2001-S	Stevens Tenotomy Scissors 115mm Curved

Supercut Stevens Tenotomy Scissors



20 Units



Code	Description
STSS2001	SuperCut Stevens Tenotomy Scissors 115mm Straight
STSC2001	SuperCut Stevens Tenotomy Scissors 115mm Curved

Mayo Scissors

20 Units



Code	Description
MSS2141	Mayo Scissors Straight Bevelled Blade 140mm
MSC2141	Mayo Scissors Curved Bevelled Blade 140mm
MSS2171	Mayo Scissors Straight Bevelled Blade 170mm
MSC2171	Mayo Scissors Curved Bevelled Blade 170mm

Iris Scissors

20 Units



Code	Description
ISS2111-S	Iris Scissors Straight 90mm
ISC2111-S	Iris Scissors Curved 90mm
ISS2110-D	Iris Scissors Straight 110mm
ISC2110-D	Iris Scissors Curved 110mm

Supercut Iris Scissors

20 Units



IMAGE COMING SOON

Code	Description
SISS2111	SuperCut Iris Scissors Straight 90mm
SISC2111	SuperCut Iris Scissors Curved 90mm

Vannas Scissors

20 Units



Code	Description
VSS2081	Vannas Scissors Straight 85mm
VSS2121	Vannas Scissors Straight 125mm
VSC2081	Vannas Scissors Curved 85mm

Metzenbaum Scissors

20 Units



Code	Description
ZSC2121-D	Baby Metzenbaum Scissors Curved 120mm



Code	Description
ZSC2141-D	Metzenbaum Scissors Curved 145mm
ZSC2161-D	Metzenbaum Scissors Curved 160mm
ZSC2181-D	Metzenbaum Scissors Curved 180mm
ZSC2201-D	Metzenbaum Scissors Curved 200mm



Code	Description
ZSS2181-D	Metzenbaum Scissors Straight 180mm
ZSS2201-D	Metzenbaum Scissors Straight 200mm

Supercut Metzenbaum Scissors

20 Units

IMAGESCOMING SOON

Code	Description
ZSC2121	SuperCut Baby Metzenbaum Scissors Curved 120mm

Needle Holders


Castroviejo Needle Holder

20 Units



Code	Description
CVNH2014	Castroviejo Needle Holder with Lock 145mm

Barraquer Needle Holders




Code	Description	Units
BNH1110C	Barraquer Needle Holder Curved 110mm	10
BNH2110C	Barraquer Needle Holder Curved 110mm	20

Aspirating Dissectors

Negus Aspirating Dissectors

5 Units



Code	Description
NDS0511	Negus Aspirating Dissector with One Hole, 2.25mm hole diameter
NDS0521	Negus Aspirating Dissector with Two Holes, 2.25mm hole diameter

Seekers & Probes

Punctum Seeker


20 Units



Code	Description
PSK2010	Punctum Seeker 110mm

Lacrimal Probes

20 Units




Code	Description
LCP2078	Lacrimal Probe, 0.7/0.8mm Diameter Wires
LCP2089	Lacrimal Probe, 0.8/0.9mm Diameter Wires

Blade Handles

Scalpel Blade Handles

20 Units





Code	Description
SBH2031	Scalpel Blade Handle No 3
SBH2041	Scalpel Blade Handle No 4



Eye Speculum

Kratz-Barraquer Eye Speculum

20 Units



Code	Description
KBS2035	Kratz-Barraquer Eye Speculum 35mm (Small)




Code	Description
KBA2040	Kratz-Barraquer Eye Speculum 40mm (Adult)

Retractors

Kilner “Catspaw” Retractor

20 Units



Code	Description
KCP2001	Kilner “Catspaw” Retractor 155mm

Desmarres Retractor

20 Units



Code	Description
DMR2012	Desmarres Retractor 160mm

Kerrison Rongeurs

Kerrison Rongeurs


10 Units



KRUF11029



KRUF11039




Code	Description
KRU1101	Kerrison Rongeur Punch 1mm wide, length 200mm with 40 degree up bite
KRU1201	Kerrison Rongeur Punch 2mm wide, length 200mm with 40 degree up bite
KRU1301	Kerrison Rongeur Punch 3mm wide, length 200mm with 40 degree up bite
KRU1401	Kerrison Rongeur Punch 4mm wide, length 200mm with 40 degree up bite
KRU1501	Kerrison Rongeur Punch 5mm wide, length 200mm with 40 degree up bite
KRUF11029	Kerrison Rongeur Punch 2mm wide, length 100mm with 90 degree up bite
KRUF11039	Kerrison Rongeur Punch 3mm wide, length 100mm with 90 degree up bite

Choppers & Hooks

Phaco Chopper

20 Units



Code	Description
PHC2001	Phaco Chopper, 0.4 x 1.5mm Blade, Rounded Bottom, 11mm, 50° Angled Shaft

Sinskey Hook

20 Units




Code	Description
SIH2041	Sinskey Hook, 11mm 45° Angled Shaft, 0.2mm Diameter x 0.5mm Hook

Nasal Specula

Hartmann Nasal Speculum

20 Units



Code	Description
HNS2011-D	Hartmann Nasal Specula, Small
HNS2021-D	Hartmann Nasal Specula, Medium
HNS2031-D	Hartmann Nasal Specula, Large

Chalazion Currettes

Chalazion Currettes

20 Units



Chalazion Curette 2.5mm

Code	Description
CHC2025	Chalazion Curette 2.5mm



Chalazion Curette 3mm

Code	Description
CHC2030	Chalazion Curette 3mm




Double Ended Chalazion Curette

Code	Description
DHC2012	Double Ended Chalazion Curette 1mm & 2mm

Tubing Products

Sterile Silicone Tube Sets

10 Units



Code	Description
DST2051004	Sterile Silicone DCR Tube 1.02 ID x 2.16 OD x 40mm long
DSTS40506	Sterile Silicone Tube 4.8 ID x 7.9 OD x 600mm long with Spigot
DSTS44506	Sterile Silicone Tube 6.4 ID x 9.5 OD x 600mm long with Spigot

EAGLE LABS® Cannulae

Anaesthesia Administration

Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
100-19	Tri-Port Sub Tenon Cannulae	2x .4mm Side Ports, Open Ended, Angled	25	1	19G	1.0
100-19C	Tri-Port Sub Tenon Cannulae	2x .4mm Side Ports, Open Ended, Curved	25	1	19G	1.0
100-21	Tri-Port Sub Tenon Cannulae	2x .4mm Side Ports, Open Ended, Angled	25	1	21G	0.8
100-21C	Tri-Port Sub Tenon Cannulae	2x .4mm Side Ports, Open Ended, Curved	25	1	21G	0.8
101-23	Atkinson Tip Retrobulbar		38	1 1/2	23G	0.6
101-25	Atkinson Tip Retrobulbar		38	1 1/2	25G	0.5
101-27	Atkinson Tip Retrobulbar		38	1 1/2	27G	0.4
102-23-32	Atkinson Tip Peribulbar		32	1 1/4	23G	0.6
102-25-32	Atkinson Tip Peribulbar		32	1 1/4	25G	0.5
102-27-32	Atkinson Tip Peribulbar		32	1 1/4	27G	0.4
102-25-22	Atkinson Tip Peribulbar		22	7/8	25G	0.5
102-27-22	Atkinson Tip Peribulbar		22	7/8	27G	0.4
103-25	Sharp Point Retrobulbar	Sharp Tip Retrobulbar	38	1 1/2	25G	0.5

Anterior Chamber Maintenance

Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
170-20	Anterior Chamber Maintenaner	Infusion Cannulae with Self Retaining Serrated Tip, 8" Flexible Tubing	4.76	3/16	20G	0.9

Irrigation, Air Injection & Viscoelastic

Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
111-16	Anterior Chamber Washout	25mm Long, Flattened Tip, Angled	25	1	16G	1.4
111-19	Anterior Chamber Washout	25mm Long, Flattened Tip, Angled	25	1	19G	1.1
111-16-BH	Bishop Harmon Cannulae	6mm Bend to Tip	25	1	16G	1.4
111-19-BH	Bishop Harmon Cannulae	8mm Bend to Tip	25	1	19G	1.1
112-25	Anterior Chamber Irrigator	25mm Straight, Round Tip	25	1	25G	0.5
112-27	Anterior Chamber Irrigator	25mm Straight, Round Tip	25	1	27G	0.4
112-30	Anterior Chamber Irrigator	25mm Straight, Round Tip	25	1	30G	0.3
110-25-4	Anterior Chamber Irrigator	4mm Bend to Tip, Angled 30°	25	1	25G	0.5
110-27-4	Anterior Chamber Irrigator	4mm Bend to Tip, Angled 30°	25	1	27G	0.4
110-30-4	Anterior Chamber Irrigator	4mm Bend to Tip, Angled 30°	25	1	30G	0.3
110-20	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	20G	0.9
110-21	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	21G	0.8
110-22	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	22G	0.7
110-23	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	23G	0.6
110-25	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	25G	0.5
110-26	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	26G	0.45
110-27	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	27G	0.4
110-30	Anterior Chamber Irrigator	7mm Bend to Tip, Angled 30°	25	1	30G	0.3
110-27-9	Anterior Chamber Irrigator	9mm Bend to Tip, Angled 30°	25	1	27G	0.4
110-30-9	Anterior Chamber Irrigator	9mm Bend to Tip, Angled 30°		1	30G	0.3
110-25-11	Anterior Chamber Irrigator	11mm Bend to Tip, Angled 30°	25	1	25G	0.5
110-27-11	Anterior Chamber Irrigator	11mm Bend to Tip, Angled 30°	25	1	27G	0.4
110-27-TM	Anterior Chamber Irrigator (Mendez)	3mm Bend to Tip, Angled 45°	10	3/8	27G	0.4
110-26 1/2	Anterior Chamber Irrigator	8mm Bend to Tip, Angled 25°	12	1/2	26G	0.45
110-30SP	Anterior Chamber Irrigator	12mm Bend to Tip, Angled 45°, Beveled	14	9/16	30G	0.3

109-27-25	FLOvisc Viscoelastic Cannulae	20G Stepped Down to 25G, 10mm Angled Tip		1	20G/25G	0.9/0.5
109-27	FLOvisc Viscoelastic Cannulae	20G Stepped Down to 27G, 10mm Angled Tip		1	20G/27G	0.9/0.4
109-22	FLOvisc Viscoelastic Cannulae	22G Stepped Down to 27G, 10mm Angled Tip		1	22G/27G	0.7/0.4
170-20W	Infusion (Wright)	7mm Bend to Beveled Tip, Angled 45°	25	1	20G	0.9
150-23	I and A Cannulae (Simcoe)	Irrigation through Cannulae, 0.4mm Side Port, 30 Aspiration through 13" (33cm) Flexible Tubing	30	13	23G	0.6

Capsulotomy

Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
120-23S	Irrigating Cystotome	Straight Tip	25	1	23G	0.6
120-25S	Irrigating Cystotome	Straight Tip	16	5/8	25G	0.5
122-23S	Irrigating Cystotome	Straight, Reverse Tip	19	3/4	23G	0.6
122-25S	Irrigating Cystotome	Straight, Reverse Tip	16	5/8	25G	0.5
120-25F-8	Irrigating Cystotome	8mm Formed Tip	16	5/8	25G	0.5
120-27F-8	Irrigating Cystotome	8mm Formed Tip	13	1/2	27G	0.4
121-25F-8	Irrigating Cystotome	8mm Formed Tip, Push/Pull	16	5/8	25G	0.5
121-27F-8	Irrigating Cystotome	8mm Formed Tip, Push/Pull	13	1/2	27G	0.4
123-25	Irrigating Cystotome	10mm Angled Tip	16	5/8	25G	0.5
123-27	Irrigating Cystotome	10mm Angled Tip	13	1/2	27G	0.4
120-25F-12	Irrigating Cystotome	12mm Formed Tip	16	5/8	25G	0.5
120-27F-12	Irrigating Cystotome	12mm Formed Tip	13	1/2	27G	0.4
121-25F-12	Irrigating Cystotome	12mm Formed Tip, Push/Pull	16	5/8	25G	0.5
122-25RF	Irrigating Cystotome	12mm Formed Tip, Reverse	16	5/8	25G	0.5
122-27RF	Irrigating Cystotome	12mm Formed Tip, Reverse	13	1/2	27G	0.4
120-01	Irrigating Cystotome	13mm Formed Tip	19	3/4	25G	0.5

Hydrodissection

Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
113-01	Nucleus Hydrodissection (Chang)	1.5mm Flat Beveled Tip, 90° Angled	25	1	27G	0.4
113-23	Nucleus Hydrodissection	7mm Flattened Tip	25	1	23G	0.6
113-25	Nucleus Hydrodissection	7mm Flattened Tip	25	1	25G	0.5
113-27	Nucleus Hydrodissection	7mm Flattened Tip	25	1	27G	0.4
113-30	Nucleus Hydrodissection	7mm Flattened Tip	25	1	30G	0.3
213-25	Nucleus Hydrodissection	7mm Flattened Tip, Vertical Flow	25	1	25G	0.5
213-27	Nucleus Hydrodissection	7mm Flattened Tip, Vertical Flow	25	1	27G	0.4
113-25-8	Nucleus Hydrodissection	8mm Flattened Tip	25	1	25G	0.5
113-25-10	Nucleus Hydrodissection	10mm Flattened Tip	25	1	25G	0.5
113-27-10	Nucleus Hydrodissection	10mm Flattened Tip	25	1	27G	0.4
113-25F-12	Nucleus Hydrodissection	12mm Formed, Flattened Tip	25	1	25G	0.5
113-27F-12	Nucleus Hydrodissection	12mm Formed, Flattened Tip	25	1	27G	0.4
118-25	Nucleus Hydrodissection	Curved, Horizontal Flattened Tip	25	1	25G	0.5
118-25S	Nucleus Hydrodissection	Curved, Vertical Flattened Tip	25	1	25G	0.5
118-27	Nucleus Hydrodissection	Curved, Horizontal Flattened Tip	25	1	27G	0.4
118-27S	Nucleus Hydrodissection	Curved, Vertical Flattened Tip	25	1	27G	0.4
217-27	Cortical Cleaving Hydrodissection (Grabow)	11mm Bend to Round Tip, Angled 45°	25	1	27G	0.4
218-27	Cortical Cleaving Hydrodissection	11mm Bend to Flattened Tip, Angled 45°	25	1	27G	0.4
313-30	Hydrodissection (Tuebingen)	75° Curved, Round Tip	25	1	30G	0.3

Hydrodelineation						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
114-30	Nucleus Hydrolineation	Curved, Beveled Tip	25	1	30G	0.3
Nucleus Removal						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
140-25-1	Irrigating Vectis	5mm Wide (1 ea.), 0.3mm Open Port			25G	0.5
140-25-2	Irrigating Vectis	5mm Wide (2 ea.), 0.3mm Open Port			25G	0.5
140-25-3	Irrigating Vectis	5mm Wide (3 ea.), 0.3mm Open Port			25G	0.5
140-27-1	Irrigating Vectis	5mm Wide (1 ea.), 0.3mm Open Port			27G	0.4
Capsule Polishing						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
124-23	Capsule Polishers (Lu)	8mm Bend, 30mm Top Port, Rounded Tip, 25 Curved, Roughened Olive Tip, Open Port	25	1	23G	0.6
125-23	Capsule Polishers (Terry)	7mm Angled, Silicone Tip	25	1	23G	0.6
125-27	Capsule Polishers (Terry)	7mm Angled, Silicone Tip	25	1	27G	0.4
125-24	Capsule Polishers (Kratz)	8mm Angled, Silicone Tip	25	1	24G	0.55
126-23	Capsule Polishers	10mm Angled, Olive Tip	25	1	23G	0.6
126-25	Capsule Polishers	10mm Angled, Olive Tip	25	1	25G	0.5
Cortex Removal						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
108-23	Simcoe Cortex Extractor	Curved, 0.3mm Inside Port	25	1	23G	0.6
105-275	Aspirating J-Hook	2.0mm x 2.0mm Straight Hook	25	1	27G	0.4
106-275	Aspirating Cortex Extractor	1.5mm x 1.5mm Straight Hook	25	1	27G	0.4
191-23	Cortex Aspiration	0.4mm Side Port, Male Adapter, Curved	25	1	23G	0.6
192-23	Bimanual I/A Cortex Removal System	0.4mm Side Port, Male/Female Adapters, Curved			21G/23G	0.8/0.6
Corneal Abrasion						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
165-25	Stromal Puncture Needle	19mm Bend to Tip, 30° Angled	25	1	25G	0.5
Fluid Administration						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
169-30L-6	Ocular Irrigation Tube	Extension Tube with Male/Female Luer Connections				
169-30L-10	Ocular Irrigation Tube	Extension Tube with Male/Female Luer Connections				
Glaucoma						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
177-33	Visco Canalostomy	5mm Beveled Tip, Angled			32G	
178-27	Glaucoma Revision Pick (Lederer)	10mm Bend to Tip, Angled 30°	13	1/2	27G	0.4
Oculoplastic						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
160-21	Lacrimal Cannulae	Curved, Open End	38	1 1/2	21G	0.8
160-25	Lacrimal Cannulae	Curved, Open End	38	1 1/2	25G	0.5
160-26	Lacrimal Cannulae	Curved, Open End	38	1 1/2	26G	0.45
161-23	Lacrimal Cannulae	Curved, 0.3mm Side Port	38	1 1/2	23G	0.6

161-25	Lacrimal Cannulae	Curved, 0.3mm Side Port	38	11 1/2	25G	0.5
163-21	Lacrimal Cannulae	Straight, 0.3mm Side Port	38	11 1/2	21G	0.8
163-25	Lacrimal Cannulae	Straight, 0.3mm Side Port	38	11 1/2	25G	0.5
162-23	Lacrimal Cannulae	25mm Straight with 13mm 23G Tip	25	11 1/2	19G/23G	1.0
Refractive (LASIK)						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
605-25	LASIK Irrigating Cannulae (Guell)	9mm Formed Bend to Tip, 6 Ports Open Ended	25	1	25G	0.5
604-25	LASIK Irrigating Cannulae (Banaji)	12mm Formed, 6x ports, Closed End	25	1	25G	0.5
606-25	LASIK Irrigating Cannulae (Burrato)	9mm Formed Bend, 2 Ports Open Ended	25	1	25G	0.5
600-16	LASIK Flap Irrigator (Vidaurri)	9mm Angled Bend to Tip, Dual Multi-Port Irrigation Tips, 2 Ports (ea.)	25	1	25G	0.5
600-06	LASIK Flap Irrigator (Vidaurri)	13mm Angled Bend to Tip, Dual Multi-Port Irrigation Tips, 4 Ports (ea.)	25	1	25G	0.5
603-25	LASIK Irrigating Cannulae (Lindstrom)	7mm Angled, 3x 0.3mm Ports	25	1	25G	0.5
600-25	LASIK Irrigating Cannulae (Nichamin)	12mm Formed, Flattened Tip, Open Irrigation	25	1	25G	0.5
600-27	LASIK Irrigating Cannulae (Nichamin)	12mm Formed, Flattened Tip, Open Irrigation	25	1	27G	0.4
610-23	LASIK Flap Spatula/Irrigator (Hersh)	9mm Angled Bend to Tip	25	1	23G	0.6
610-25	LASIK Flap Spatula/Irrigator (Hersh)	9mm Angled Bend to Tip	25	1	25G	0.5
613-26	LASIK Irrigating Spatula (Slade)	9mm Formed Bend to Tip	25	1	26G	0.45
614-25	LASIK Irrigating Cannulae (Soloman)	Formed Spatula Tip, 3 Ports	25	1	25G	0.5
617-27	LASIK Irrigating Cannulae	10mm Formed Bend to Tip, Flattened, 4 Ports	25	1	27G	0.4
Vitreoretinal						
Code	Name	Description	Length (mm)	Length (inch)	Gauge	Dia (mm)
166-25	Subretinal Fluid Cannulae (Gizzard)	Silicone Sleeve Extended 6mm Beyond Tip, Straight	37	1 7/16	20G	0.9

Eagle Labs® Knives & Blades

Surgical Knives

Code	Name	Description
EG-5010	Double Bevel Sideport Knife	1.0mm, Angled, Full Handle, Double Bevel
EG-5012	Double Bevel Sideport Knife	1.2mm, Angled, Full Handle, Double Bevel
EG-6821	Spoon Knife	3.0mm Angled, Full Handle, Double Bevel
EG-5555	MVR Knife	19GA, Full Handle
EG-5560	MVR Knife	20GA, Full Handle
EG-5563	MVR Knife	23GA, Full Handle
EG-5555	Clear Corneal Knife	2.2mm, Angled, Full Handle, Double Bevel
EG-5663	Clear Corneal Knife	2.4mm, Angled, Full Handle, Double Bevel
EG-5663	Clear Corneal Knife	2.75mm, Angled, Full Handle, Double Bevel
EG-1130	Clear Corneal Knife	3.0mm, Angled, Full Handle, Double Bevel
EG-1132	Clear Corneal Knife	3.2mm, Angled, Full Handle, Double Bevel
EG-1327	Clear Corneal Knife, Parallel Sides	2.75mm, Angled, Full Handle, Double Bevel
EG-2022	Trapezoid Clear Corneal/ Implant Knife	2.0mm - 2.2mm, Angled, Full Handle, Bevel Up
EG-2832	Trapezoid Clear Corneal/ Implant Knife	2.85mm - 3.2mm, Angled, Full Handle, Bevel Up
EG-1500	Stab Knife	15°, Full Handle
EG-2250	Stab Knife	22.5°, Full Handle
EG-3000	Stab Knife	30°, Full Handle
EG-4500	Stab Knife	45°, Full Handle
EG-1022	Phaco Slit Knife	2.2mm, Angled, Full Handle, Bevel Up
EG-1024	Phaco Slit Knife	2.4mm, Angled, Full Handle, Bevel Up
EG-5510	Phaco Slit Knife	2.5mm, Angled, Full Handle, Bevel Up
EG-5568	Phaco Slit Knife	2.65mm, Angled, Full Handle, Bevel Up
EG-4738	Phaco Slit Knife	2.75mm, Angled, Full Handle, Bevel Up
EG-4740	Phaco Slit Knife	2.85mm, Angled, Full Handle, Bevel Up
EG-4741	Phaco Slit Knife	2.85mm, Straight, Full Handle
EG-5570	Phaco Slit Knife	3.0mm, Angled, Full Handle, Bevel Up
EG-5571	Phaco Slit Knife	3.0mm, Straight, Full Handle
EG-5572	Phaco Slit Knife	3.2mm, Angled, Full Handle, Bevel Up
EG-5573	Phaco Slit Knife	3.2mm, Straight, Full Handle
EG-9350	IOL Inserter Knife	5.2mm, Angled, Full Handle, Bevel Up
EG-9355	IOL Inserter Knife	5.2mm, Angled, Full Handle, Bevel Up
EG-9402	Crescent Knife	2.0mm, Angled, Full Handle, Bevel Up

Surgical Blades/holders		
Code	Name	Description
EL-5910	Needle Blade	3.0mm
EL-6600	Lamellar Blade	60°
EL-5700	Scleral Blade	
EL-6400	Miniature Edged Blade	Round Tip
EL-6500	Miniature Edged Blade	Lance Tip
EL-6700	Miniature Edged Blade	Curved Tip
EL-6900	Miniature Edged Blade	Sharp All Around
EL-7513	Microsurgical Blade	3.0mm, 15°, Blue
EL-7514	Microsurgical Blade	3.0mm, 15°, Blue, with Handle
EL-7515	Microsurgical Blade	5.0mm, 15°, Green
EL-7516	Microsurgical Blade	5.0mm, 15°, Green, with Handle
Precision LASIK Blades		
Code	Name	Description
ES-3030 +10	LASIK Blade	Plus 10, Compatible with Hansatome™ microkeratome
ES-3030 -10	LASIK Blade	Minus 10, Compatible with Hansatome™ microkeratome
ES-3030 -20	LASIK Blade	Minus 20, Compatible with Hansatome™ microkeratome
ES-3030 -30	LASIK Blade	Minus 30, Compatible with Hansatome™ microkeratome
ES-3030 -PLANO	LASIK Blade	Plano, Compatible with Hansatome™ microkeratome
ES-5050 +10	LASIK Blade	Plus 10, Compatible with Zyoptix™ XP microkeratome
ES-5050 +20	LASIK Blade	Plus 20, Compatible with Zyoptix™ XP microkeratome
ES-5050 -10	LASIK Blade	Minus 10, Compatible with Zyoptix™ XP microkeratome
ES-5050 -20	LASIK Blade	Minus 20, Compatible with Zyoptix™ XP microkeratome
ES-5050 -30	LASIK Blade	Minus 30, Compatible with Zyoptix™ XP microkeratome
ES-5050 -PLANO	LASIK Blade	Plano, Compatible with Zyoptix™ XP microkeratome

Precision LASIK Blades		
Code	Name	Description
ES-7070 +10	LASIK Blade	Plus 10, Moria™ M2 microkeratome
ES-7070 -10	LASIK Blade	Minus 10, Moria™ M2 microkeratome
ES-7070 -20	LASIK Blade	Minus 20, Moria™ M2 microkeratome
ES-7070 -30	LASIK Blade	Minus 30, Moria™ M2 microkeratome
ES-7070 -PLANO	LASIK Blade	Plano, Moria™ M2 microkeratome
ES-9090 +10	LASIK Blade	Plus 10, Nidek™ microkeratome
ES-9090 +20	LASIK Blade	Plus 20, Nidek™ microkeratome
ES-9090 -10	LASIK Blade	Minus 10, Nidek™ microkeratome
ES-9090 -20	LASIK Blade	Minus 20, Nidek™ microkeratome
ES-9090 -30	LASIK Blade	Minus 30, Nidek™ microkeratome
ES-9090 -PLANO	LASIK Blade	Plano, Nidek™ microkeratome

EYETEC® DCR & Lacrimal

Lacrimal Intubation

Code	Name	Description
51-917	Crawford-type Intubation Set	Retrieval device for Crawford-Type Intubation, Box of 5
51-918	Crawford-type Intubation Set	Crawford-Type Intubation Set, Box of 1

Canaliculus Intubation Set

Code	Name	Description
51-920	Canaliculus Intubation Set	Canaliculus Intubation Set, Box of 1

Code	Name	Description
51-922	DCR Bodkins Straight	DCR Bodkins Straight, Box of 3
51-923	DCR Bodkins Angled	DCR Bodkins Angled, Box of 3

EYETEC® PVA & Cellulose

Eye Spears and Points

Code	Name	Description
40-400	EYETEC® PVA Spears	PVA Eye Spears, Pack of 5 box of 90
40-405	EYETEC® PVA Spears	PVA Eye Spears, Pack of 5 box of 180
40-406	EYETEC® PVA Spears	PVA Points, Pack of 10 box of 180
40-410	EYETEC® Cellulose Spears	Cellulose Spears, Pack of 5 box of 100
40-411	EYETEC® Cellulose Spears	Cellulose Points, Pack of 5 box of 100
40-415	EYETEC® Cellulose Spears	Cellulose Spears, Pack of 10 box of 200
40-416	EYETEC® Cellulose Spears	Cellulose Points, Pack of 10 box of 250

Corneal Shields

Code	Name	Description
40-420	EYETEC® Light Shields	PVA Corneal Light Shields 7mm Diameter, Pack of 1 box of 20
40-421	EYETEC® Light Shields	PVA Corneal Light Shields 8mm Diameter, Pack of 1 box of 20
40-423	EYETEC® Light Shields	PVA Corneal Light Shields 10mm Diameter, Pack of 1 box of 20

Eye Shields

Code	Name	Description
40-820	EYETEC® LASIK Eye Shields	PVA LASIK Eye Shields 9mm Diameter/7mm Edge, Pack of 1 box of 20
40-821	EYETEC® LASIK Eye Shields	PVA LASIK Eye Shields 8mm Diameter/6mm Edge, Pack of 1 box of 20
40-822	EYETEC® LASIK Eye Shields	PVA LASIK Eye Shields 4mm Edge/8mm Edge, Pack of 1 box of 20

Eye Drains

Code	Name	Description
40-430	EYETEC® Eye Drains	PVA Eye Drains 80cc
40-435	EYETEC® Eye Drains	PVA Eye Drains 400cc
40-830	EYETEC® LASIK Eye Drains	PVA LASIK Eye Drains 42 x 17 x 11.5mm Opening, Pack of 1 box of 20
40-831	EYETEC® LASIK Eye Drains	PVA LASIK Eye Drains 40 x 14 x 10.5mm Opening, Pack of 1 box of 20
40-431	EYETEC® Eye Wicks	Eye Wick, Pack of 2 box of 40
40-470	EYETEC® Eye Wicks	Wick and Wipe, Pack of 1 box of 20

Instrument Wipes

Code	Name	Description
40-900	EYETEC® Instrument Wipe	Instrument Wipe, Pack of 1 box of 20

Diamond Cleaning Blocks

Code	Name	Description
40-462	EYETEC® Diamond Cleaning Block	Diamond Cleaning Block, Pack of 1 box of 20

LACRIMEDICS® Punctal Plugs

VisiPlug®

Code	Name	Description
1638	VisiPlug® Dissolvable Lacrimal Plugs	Dissolvable VisiPlug® Lacrimal Plugs, 0.4mm
1639	VisiPlug® Dissolvable Lacrimal Plugs	Dissolvable VisiPlug® Lacrimal Plugs, 0.5mm

ComfortTip™

Code	Name	Description
2003	ComfortTip™ Punctal Occluder	ComfortTip™ Punctal Occluder, X-Small
2005	ComfortTip™ Punctal Occluder	ComfortTip™ Punctal Occluder, Small
2008	ComfortTip™ Punctal Occluder	ComfortTip™ Punctal Occluder, Medium
2010	ComfortTip™ Punctal Occluder	ComfortTip™ Punctal Occluder, Large

AccuFlo™

Code	Name	Description
3000	AccuFlo™ Punctal Occluder	AccuFlo™ Punctal Occluder, Small
3001	AccuFlo™ Punctal Occluder	AccuFlo™ Punctal Occluder, Medium
3002	AccuFlo™ Punctal Occluder	AccuFlo™ Punctal Occluder, Large
3003	AccuFlo™ Punctal Occluder	AccuFlo™ Punctal Occluder, X-Large

CORONET® Orbital Implants

Dissolvable Lacrimal Plugs

Code	Name	Implant Size
54-300-14	CORONET® Hydroxyapatite Orbital Implant	14 mm
54-300-16	CORONET® Hydroxyapatite Orbital Implant	16 mm
54-300-17	CORONET® Hydroxyapatite Orbital Implant	17 mm
54-300-18	CORONET® Hydroxyapatite Orbital Implant	18 mm
54-300-19	CORONET® Hydroxyapatite Orbital Implant	19 mm
54-300-20	CORONET® Hydroxyapatite Orbital Implant	20 mm
54-300-21	CORONET® Hydroxyapatite Orbital Implant	21 mm
54-300-22	CORONET® Hydroxyapatite Orbital Implant	22 mm
54-305-14	CORONET® Alumina Orbital Implants	14 mm
54-305-16	CORONET® Alumina Orbital Implants	16 mm
54-305-17	CORONET® Alumina Orbital Implants	17 mm
54-305-18	CORONET® Alumina Orbital Implants	18 mm
54-305-19	CORONET® Alumina Orbital Implants	19 mm
54-305-20	CORONET® Alumina Orbital Implants	20 mm
54-305-21	CORONET® Alumina Orbital Implants	21 mm
54-305-22	CORONET® Alumina Orbital Implants	22 mm

CORONET® Corneal Trephine

Vacuum Trephine

Code	Name	Sizes (mm)
51-840	CORONET® Recipient Corneal Vacuum Trephine	6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.50
51-842	CORONET® PK Clear-Vue™ Recipient Corneal Vacuum Trephine	6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.50
51-844	CORONET® DALK Clear-Vue™ Recipient Vacuum Trephine	7.50 8.00 8.50

Trephine Punch

Code	Name	Sizes (mm)
51-835	CORONET® Recipient Corneal Vacuum Trephine	6.00 6.25 6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.25 9.50 9.75
51-850	CORONET® Donor Corneal Trephine Punch	6.00 6.25 6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.25 9.50 9.75
51-855	CORONET® Donor Corneal Trephine Punch	7.50 8.00 8.50 9.00 9.50 10.0

Guarded Trephine

Code	Name	Sizes (mm)
51-860	Guarded Corneal Trephines 300µm	7.00 7.50 8.00 8.50
51-865	Guarded Corneal Trephines 400µm	7.00 7.50 8.00 8.50

Standard Trephine

Code	Name	Sizes (mm)
51-900	CORONET® Standard Corneal Trephine	6.00 6.25 6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.25 9.50 9.75 10.0
51-901	CORONET® Standard Corneal Trephine with Cross Hair	6.00 6.25 6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.25 9.50 9.75 10.0

Long Handed Trephine

Code	Name	Sizes (mm)
51-903	CORONET® Long Handed Corneal Trephine	6.00 6.25 6.50 6.75 7.00 7.25 7.50 7.75 8.00 8.25 8.50 8.75 9.00 9.25 9.50 9.75 10.0

Extractor Trephine

Code	Name	Sizes (mm)
51-905	CORONET® Extractor Trephine	15.0 17.0 18.0 20.0

EndoGlide™

Code	Name	Description
51-823	CORONET® EndoGlide™ Ultrathin	4.5mm Sceral Tunnel Incision, Double Coil Technique
51-826	CORONET® DMEK EndoGlide™	2.65mm Clear Corneal Incision, Tri-Fold Technique

EndoGlide™ Forceps

Code	Name	Description
53-951	CORONET® Tan EndoGlide™ Loading Forceps	Tissue Placement, 23G
53-955	CORONET® Tan EndoGlide™ Loading Forceps	Tissue Placement, 27G
53-952	CORONET® Tan EndoGlide™ Loading Forceps	Tissue Placement, 23G
53-956	CORONET® Tan EndoGlide™ Loading Forceps	Tissue Placement, 27G

Artificial Anterior Chamber

Code	Name	Description
51-935	CORONET® Artificial Anterior Chamber	5 Year Shelf Life, Single Use, Sterile, 1 per box
51-300	CORONET® Surgical Trainer	Surgical Trainer, non-sterile, 1 per box
51-3010	CORONET® Surgical Trainer	Surgical Trainer Refill Kit, 10 per box

Tubing Products		
Code	Name	Description
DST2051004	Sterile Silicone Tube Sets	Sterile Silicone DCR Tube 1.02 ID x 2.16 OD x 40mm long
DSTS40506	Sterile Silicone Tube Sets	Sterile Silicone Tube 4.8 ID x 7.9 OD x 600mm long with Spigot
DSTS44506	Sterile Silicone Tube Sets	Sterile Silicone Tube 6.4 ID x 9.5 OD x 600mm long with Spigot

DTR MEDICAL® Single-Use

Forceps

Code	Name	Description
BTF2115	Typing Forceps	Birks Typing Forceps 110mm, 5mm Platform
STF2905	Typing Forceps	Straight Typing Forceps 90mm, 5mm Platform
STF2725	Typing Forceps	Straight Typing Forceps 72mm, 5mm Tying Platform
SPF2001-S	Fine Suture Forceps	Fine Suture Forceps Plain 105mm
STF2725	Fine Suture Forceps	Straight Tying Forceps 72mm, 5mm Tying Platform
TCB2013	Cyst Forceps	Birmingham Tarsal Cyst Forceps 84mm, 13mm Clamp
TCB2016	Cyst Forceps	Birmingham Tarsal Cyst Forceps 84mm, 16mm Clamp
TCL2928	Cyst Forceps	Lambert Tarsal Cyst Forceps 92mm, 8mm Clamp
PHF1085	Pierse Hoskins Forceps	Pierse Hoskins Forceps 85mm, 0.2mm Notched Tips
CFS2001-D	Crile Artery Forceps	Crile Artery Forceps 145mm Straight
CFC2001-D	Crile Artery Forceps	Crile Artery Forceps 145mm Curved on Flat
ADS5001-S	Adson Forceps	Adson Non-Toothed Forceps 120mm
ATS5001-S	Adson Forceps	Adson Toothed Forceps 120mm
ATP2012-D	Adson Forceps	Adson Forceps 1x2 with plateau 120mm
ATP2015-D	Adson Forceps	Adson Forceps 1x2 with plateau 150mm
ATS2012F	Adson Forceps	Fine Tip Adson Toothed Forceps 1x2 with plateau 120mm
EPF2085	Epilation Forceps	Epilation Forceps 85mm, 2 x 3mm Plain Tips
SMF2905	St Martins Forceps	St. Martins Forceps 1:2 Toothed 90mm, 5mm Platform, 0.4mm Tips
MFF2113	Moorfield Forceps	Moorfield Forceps 105mm, 13mm Serration
MQC2012	Mosquito Forceps Curved	Mosquito Forceps Curved 120mm
MQC2025-D	Mosquito Forceps Curved	Mosquito Forceps Curved 125mm
MQC2125-D	Mosquito Forceps Curved	Mosquito Forceps Curved Toothed 1x2 125mm
MQT5001-S	Mosquito Forceps Straight	Mosquito Forceps Straight 125mm
CXF2001	Capsulorhexis Forceps	Capsulorhexis Forceps, Utrata Tips, 12mm Curved Shafts, Round Handles

Bipolar Forceps

Code	Name	Description
SU760508044	Cawthorne Bipolar Forceps	Cawthorne Bipolar Forceps 80 x 0.5mm, 3m, 4mm Plugs
SU761008044	Cawthorne Bipolar Forceps	Cawthorne Bipolar Forceps 80 x 1mm, 3m, 4mm Plugs
SU760511044	Cawthorne Bipolar Forceps	Cawthorne Bipolar Forceps 110 x 0.5mm, 3m, 4mm Plugs
SU761011044	Cawthorne Bipolar Forceps	Cawthorne Bipolar Forceps 110 x 1mm, 3m, 4mm Plugs
SU781020044	Bayonet Bipolar Forceps	Bayonet Bipolar Forceps 200 x 1mm, 3m, 4mm Plugs
SU782020044	Bayonet Bipolar Forceps	Bayonet Bipolar Forceps 200 x 2mm, 3m, 4mm Plugs

Scissors

Code	Name	Description
WSC2115	Westcott Scissors	Westcott Scissors Curved 115mm, 15mm Blade
TSS2001-S	Stevens Tenotomy Scissors Straight	Stevens Tenotomy Scissors 115mm Straight
TSC2001-S	Stevens Tenotomy Scissors Curved	Stevens Tenotomy Scissors 115mm Curved
STSS2001	Supercut Stevens Tenotomy Scissors 20	SuperCut Stevens Tenotomy Scissors 115mm Straight
STSC2001	Supercut Stevens Tenotomy Scissors 20	SuperCut Stevens Tenotomy Scissors 115mm Curved

Scissors		
Code	Name	Description
MSS2141	Mayo Scissors	Mayo Scissors Straight Bevelled Blade 140mm
MSC2141	Mayo Scissors	Mayo Scissors Curved Bevelled Blade 140mm
MSS2171	Mayo Scissors	Mayo Scissors Straight Bevelled Blade 170mm
MSC2171	Mayo Scissors	Mayo Scissors Curved Bevelled Blade 170mm
ISS2111-S	Iris Scissors	Iris Scissors Straight 90mm
ISC2111-S	Iris Scissors	Iris Scissors Curved 90mm
ISS2110-D	Iris Scissors	Iris Scissors Straight 110mm
ISC2110-D	Iris Scissors	Iris Scissors Curved 110mm
SISS2111	Supercut Iris Scissors	SuperCut Iris Scissors Straight 90mm
SISC2111	Supercut Iris Scissors	SuperCut Iris Scissors Curved 90mm
VSS2081	Vannas Scissors	Vannas Scissors Straight 85mm
VSS2121	Vannas Scissors	Vannas Scissors Straight 125mm
VSC2081	Vannas Scissors	Vannas Scissors Curved 85mm
ZSC2121-D	Metzenbaum Scissors	Baby Metzenbaum Scissors Curved 120mm
ZSC2141-D	Metzenbaum Scissors	Metzenbaum Scissors Curved 145mm
ZSC2161-D	Metzenbaum Scissors	Metzenbaum Scissors Curved 160mm
ZSC2181-D	Metzenbaum Scissors	Metzenbaum Scissors Curved 180mm
ZSC2201-D	Metzenbaum Scissors	Metzenbaum Scissors Curved 200mm
ZSS2181-D	Metzenbaum Scissors	Metzenbaum Scissors Straight 180mm
ZSC2121	Supercut Metzenbaum Scissors	SuperCut Baby Metzenbaum Scissors Curved 120mm
Needle Holders		
Code	Name	Description
CVNH2014	Castroviejo Needle Holder	Castroviejo Needle Holder with Lock 145mm
BNH1110C	Barraquer Needle Holders	Barraquer Needle Holder Curved 110mm, 10 Units
BNH2110C	Barraquer Needle Holders	Barraquer Needle Holder Curved 110mm, 20 Units
Aspirating Dissectors		
Code	Name	Description
NDS0511	Negus Aspirating Dissectors	Negus Aspirating Dissector with One Hole, 2.25mm hole diameter
NDS0521	Negus Aspirating Dissectors	Negus Aspirating Dissector with Two Holes, 2.25mm hole diameter
Seekers & Probes		
Code	Name	Description
PSK2010	Punctum Seeker	Punctum Seeker 110mm
LCP2078	Lacrimal Probes	Lacrimal Probe, 0.7/0.8mm Diameter Wires
LCP2089	Lacrimal Probes	Lacrimal Probe, 0.8/0.9mm Diameter Wires
Blade Handles		
Code	Name	Description
SBH2031	Scalpel Blade Handles	Scalpel Blade Handle No 3
SBH2041	Scalpel Blade Handles	Scalpel Blade Handle No 4
Eye Speculum		
Code	Name	Description
KBS2035	Kratz-Barraquer Eye Speculum	Kratz-Barraquer Eye Speculum 35mm (Small)
KBA2040	Kratz-Barraquer Eye Speculum	Kratz-Barraquer Eye Speculum 40mm (Adult)

Retractors		
Code	Name	Description
KCP2001	Kilner "Catspaw" Retractor	Kilner "Catspaw" Retractor 155mm
DMR2012	Desmarres Retractor	Desmarres Retractor 160mm
Punch		
Code	Name	Description
KRU1101	Kerrison Rongeurs	Kerrison Rongeur Punch 1mm wide, length 200mm with 40 degree up bite
KRU1201	Kerrison Rongeurs	Kerrison Rongeur Punch 2mm wide, length 200mm with 40 degree up bite
KRU1301	Kerrison Rongeurs	Kerrison Rongeur Punch 3mm wide, length 200mm with 40 degree up bite
KRU1401	Kerrison Rongeurs	Kerrison Rongeur Punch 4mm wide, length 200mm with 40 degree up bite
KRU1501	Kerrison Rongeurs	Kerrison Rongeur Punch 5mm wide, length 200mm with 40 degree up bite
KRUF11029	Kerrison Rongeurs	Kerrison Rongeur Punch 2mm wide, length 100mm with 90 degree up bite
KRUF11039	Kerrison Rongeurs	Kerrison Rongeur Punch 3mm wide, length 100mm with 90 degree up bite
Choppers & Hooks		
Code	Name	Description
PHC2001	Phaco Chopper	Phaco Chopper, 0.4 x 1.5mm Blade, Rounded Bottom, 11mm, 50° Angled Shaft
SIH2041	Sinsky Hook	Sinsky Hook, 11mm 45° Angled Shaft, 0.2mm Diameter x 0.5mm Hook
Nasal Specula		
Code	Name	Description
HNS2011-D	Hartmann Nasal Speculum	Hartmann Nasal Specula, Small
HNS2021-D	Hartmann Nasal Speculum	Hartmann Nasal Specula, Medium
HNS2031-D	Hartmann Nasal Speculum	Hartmann Nasal Specula, Large
Chalazion Curettes		
Code	Name	Description
CHC2025	Chalazion Curettes	Chalazion Curette 2.5mm
CHC2030	Chalazion Curettes	Chalazion Curette 3mm
DHC2012	Chalazion Curettes	Double Ended Chalazion Curette 1mm & 2mm
Tubing Products		
Code	Name	Description
DST2051004	Sterile Silicone Tube Sets	Sterile Silicone DCR Tube 1.02 ID x 2.16 OD x 40mm long
DSTS40506	Sterile Silicone Tube Sets	Sterile Silicone Tube 4.8 ID x 7.9 OD x 600mm long with Spigot
DSTS44506	Sterile Silicone Tube Sets	Sterile Silicone Tube 6.4 ID x 9.5 OD x 600mm long with Spigot

USA Sales:

Summit Medical

An Innovia Medical Company.

815 Vikings Parkway, Suite 100. St. Paul, MN 55121, USA

☎ 1.888.229.2875

✉ customerservice@innoviamedical.com

Rest of World Sales:

Network Medical Products

An Innovia Medical Company.

Coronet House, Kearsley Road, Ripon, North Yorkshire, HG4 2SG, UK

☎ +44 (0) 1765 609555

✉ networkmedicalinfo@innoviamedical.com

Innovia Medical Ophthalmic Products

Innovia Medical.

815 Vikings Parkway, Suite 100. St. Paul, MN 55121, USA

International +44 (0) 1765 609555 | US 1.888.229.2875

www.innoviamedical.com