

Shippert

aesthetics



By Summit Medical LLC



About Innovia Medical

Innovia Medical is a family of specialty surgical companies that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com

4 TISSU-TRANS

Filtron
Filtron Stand

6 ASPIRATOR ACCESSORIES

Waste Canisters
Aspirator Filters

7 TUBING

Infiltration
Aspiration

8 FAT GRAFTING GUIDE**9 CANNULAS**

Cannula Cleaning
How to Order
Hole Patterns
Accessories
Infusion Handle & Cannula
Disposable Cannulas

18 WOUND DRAINAGE

Kits

19 COMPRESSION GARMENTS

Size Chart
Breasts
Face

22 RHINOPLASTY

Nasal Packing
Epistaxis Balloons
Internal Nasal Splints
External Nasal Splints
Cross Reference Guide

SUBMIT PURCHASE ORDERS TO:



**AN
INNOVIA
MEDICAL
COMPANY**

Summit Medical LLC
815 Vikings Parkway, Suite 100
St. Paul, MN 55121 | USA

P: 1-888-229-2875 | +1 651-789-3939
F: 1-888-229-1941 | +1 651-789-3979
www.summitmedicalusa.com

Our products are available through distributors and Group Purchasing Organization contracts. Medline, Henry Schein, McKesson, Cardinal Health and Owens & Minor are some of the participating companies. Please call customer service with any questions.

FILTRON®

The FILTRON is a closed, sterile filtering system that filters and collects viable fat, removing waste, non-viable cells and excess fluid during fat harvest. The FILTRON does not change the normal flow of the physician's technique and does not require decant time. The FILTRON unit is placed in-line between the harvest cannula and the waste canister/suction source. Suction pressure is used at 15.0" Hg / 50.8kPa and fat is drawn into the mesh filter within the canister. Waste (blood, epinephrine, xylocaine) and non-viable cells are pulled through the filter and are suctioned out to the waste canister. When the desired volume is collected, the transfer tube is attached and the FILTRON unit is inverted. Virtually Pure Fat* is then drawn directly into a luer lock syringe and is ready to reinject.

Each unit includes six feet of harvest tubing and a transfer tube with luer lock tip.

PART NO.	VOLUME	QTY/BOX
3-TT-FILTRON 100	100cc	1
3-TT-FILTRON 250	250cc	1
3-TT-FILTRON 500	500cc	1
3-TT-FILTRON 1000	1000cc	1
3-TT-FILTRON 2000	2000cc	1

FAT HARVESTED

WASTE FILTERED AND
SUCTIONED OUT



VIRTUALLY PURE FAT*
RETAINED IN FILTER

- Closed, Single-Use System
- Sterile, Disposable
- Separates Tissue and Fluids by Filtration
- Filters and Removes Waste, Non-Viable Cells and Excess Fluid During Harvest
- No Decant Time Needed
- Filter Technology
- Eliminates Centrifuge
- Fat immediately ready for reinjection
- Compatible with Most Aspirators and Instruments

**Demonstrated in a study to produce virtually pure fat containing trace amounts (1%) of residual oil and fluid.*

**2000cc****1000cc****500cc****250cc****100cc**

AVAILABLE IN FIVE VOLUMES


VOLUME IS TOTAL FAT READY FOR RE-INJECTION WITH THE WASTE FRACTION REMOVED

FILTRON Stand

Canister stand has four universal brackets and has five legs for stability. FILTRON devices can be hung upright for harvest or upside down for fat removal. The extended height FILTRON stand raises the sterile field.

Tip: Hang a sterile drape or mayo stand cover over the stand prior to hanging a FILTRON to maintain sterility.


PART NO.	DESCRIPTION	HEIGHT	QTY/BOX
3-T-STAND-DLX	Standard Height	34"	1
3-T-STAND-TALL	Extended Height	52"	1



FILTRON IV Pole Bracket

This reusable, autoclavable universal bracket with clamp fits any IV pole. Devices can be hung upright for harvest or upside down for fat removal.

PART NO.	QTY/BOX
3-TT-BRACK-HK	1



FILTRON Holder for IV Pole

Reusable/Autoclavable holder connects to the FILTRON and designed to fit any IV pole. Holder comes in two different size to accommodate all sizes of the FILTRON. Easy rotation after fat is collected in the canister and ready for reinjection.

PART NO.	DESCRIPTION	QTY/BOX
TT-HOLDER SM	Fits FILTRON 100 and FILTRON 250	1
TT-HOLDER LG	Fits FILTRON 500, FILTRON 1000, & FILTRON 2000	1



WASTE CANISTERS

Implosion resistant waste collection canisters are disposable, liner-free and have sealable ports. Canisters offer easy view of aspirate and have a high pathogen filter.

PART NO.	VOLUME	QTY/BOX
LS7020*	2000cc	10
LS7030	3000cc	10

**most popular*



Waste Canister Holder for IV Pole

Autoclavable holder connects the LS7020 and LS7030 Waste Canisters to any IV Pole.

- Allows for varying placement of the waste containers to maximize surgical space
- Reusable

PART NO.	QTY/BOX
TT-HOLDER NR	1

ASPIRATOR FILTERS



Replacement Absolute Biofilter

Replacement Absolute Biofilter prevents aerosolized pathogens from your exhaust from seeping into the air. Biofilters have a 15 hour replacement interval and each is 99.999% DOP rated for 0.02 micron particles.

PART NO.	QTY/BOX
FL3000-05	5
FL3000M-05*	5

**works with Microaire aspirator*

INFILTRATION

Infiltration Tubing with Pump Segment

Sterile auto-fuse infiltration tubing has a luer connector and a spike on opposite ends. Disposable tubing is compatible with all peristaltic head infiltration pumps.

PART NO.	DESCRIPTION	OUTER DIAMETER	QTY/BOX
AFTD	Single Spike	$\frac{3}{8}$ " (.375")	10
AFTDY	Double Spike	$\frac{3}{8}$ " (.375")	10



SINGLE SPIKE



DOUBLE SPIKE

ASPIRATION

High Pressure Tubing for Liposuction

Sterile PVC tubing is latex-free and has thick walls for high pressure suction. Disposable tubing is compatible with all available aspirators and has excellent drapability.

PART NO.	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY/BOX
LS7010	$\frac{3}{8}$ " (.375")	$\frac{5}{8}$ " (.625")	8 ft	10
LS7012	$\frac{3}{8}$ " (.375")	$\frac{3}{4}$ " (.625")	10 ft	10



PAL Microaire® Tapered Tubing

Sterile tapered tubing is individually packaged, 12 feet in length, and disposable. Pair tubing with the luer adapter to use with LL cannula and PAL hand piece.

PART NO.	DESCRIPTION	LENGTH	QTY/BOX
PAL TUBING	Tubing	12 ft	10
3-PAL-CPLR	Luer Adapter	N/A	1

FAT GRAFTING GUIDE

				<100cc	Eyelids	Lips	Face	Body (Small)	Body (Large)	Breast	Buttocks
Infusion Cannulas	3-C-BAN-16x15		Luer Lock 6 Hole Spray 16gaugex15cm	•	•	•	•				
	3-C-GAN-16x15										
	3-C-BAN-14x30		Luer Lock 6 Hole Spray 14gaugex30cm					•	•	•	•
	3-C-GAN-14x30										
Harvest Cannulas	3-C-24H-3x15		24 hole 3mmx15cm	•	•	•	•				
	3-H-BA-24H-3x15 (handle)										
	3-H-BA-4H-3x30 (handle)		24 hole 3mmx30cm					•	•	•	•
Collection Units	3-TT-FILTRON 100		100cc		•	•	•				
	3-TT-FILTRON 250		250cc				•	•			
	3-TT-FILTRON 500		500cc				•	•	•	•	
	3-TT-FILTRON 1000		1000cc						•	•	
	3-TT-FILTRON 2000		2000cc						•	•	•
Injection Cannulas	3-C-1H2-18x8		Luer Lock 1 Hole, Type II 18gaugex8cm	•	•	•					
	3-C-1H2-16x8		Luer Lock 1 Hole, Type II 16gaugex8cm	•			•				
	3-C-1H2-14x8		Luer Lock 1 Hole, Type II 14gaugex8cm	•				•		•	
	3-C-1H2-3x15		Luer Lock 1 Hole, Type II 3mmx15cm						•	•	•

This is a guide only. Treatment plan is determined by healthcare professional.

CANNULA CLEANING BRUSHES

Non-sterile stainless steel rod with medical grade nylon bristles.

PART NO.	SIZE	QTY/BOX
3-C-BR2	2mm x 42cm (0.07" x 16.5")	1
3-C-BR3	3mm x 48cm (0.11" x 19")	1
3-C-BR4	4mm x 48cm (0.15" x 19")	1
3-C-BR5	5mm x 48cm (0.2" x 19")	1
3-C-BR6	6mm x 48cm (0.23" x 19")	1
3-C-BR8	8mm x 48cm (0.31" x 19")	1
3-C-BR10	10mm x 48cm (0.39" x 19")	1



HOW TO ORDER

A multitude of hole patterns, diameters and length combinations are available. Use the example below to help build your cannula. Replace **XX** with desired hole pattern, diameter, and length.

LUER LOCK

3 - C - XX - XX - XX
 Part Number Hole Pattern Diameter Length

HANDLE

3 - H - BA - XX - XX - XX
 Part Number Hole Pattern Diameter Length

Diameter

18gauge

16gauge

14gauge

2mm

3mm

4mm

5mm

6mm

Length

8cm

10cm

12cm

15cm

20cm

25cm

30cm

35cm

40cm

45cm

CANNULA AND HANDLES

Blue Aluminum Handle with Cannula

3-H-BA-Hole Pattern-DiameterxLength


1" Autoclavable Blue Aluminum Handle with Cannula of choice.



Luer Lock Handle

3-H-BAL

1" Autoclavable Blue Aluminum Handle with Luer Lock Hub.



Luer Lock Cannulas

3-C-Hole Pattern-DiameterxLength

Autoclavable stainless steel Luer Lock Cannula.



Toomey 60cc Cannulas

3-C-T60-Hole Pattern-DiameterxLength

Autoclavable Stainless Steel 60cc hub fits 60cc Toomey Syringes and Handles.



Toomey Handle

3-H-TH60

Autoclavable Black Delrin 60cc Toomey Handle.

CANNULA HOLE PATTERNS

INFILTRATION

BAN 6 Hole Body Spray

GAN 6 Hole Garden Spray

MH 6 Multi Hole Spray

HARVEST

4H Large Volume 4 Hole Harvester

24H Small Volume 24 Hole Harvester

STANDARD

MC Mercedes

TP Triport

SP Spatula

BB Becker Basket

REINJECTION

1H1 One Hole Type I

1H2 One Hole Type II (Coleman)

1HV One Hole V Type

Contact customerservice@innoviamedical.com for custom cannulas.
Please allow 10-14 business days for cannula orders.
No returns on cannula orders.


ACCESSORIES



Disposable Luer to Luer Transfer

Sterile female Luer Lock to Luer Lock transfer system is made of Polypropylene and has large finger grips to protect against contamination during tissue transfer or emulsification.


PART NO.	QTY/BOX
3-S-13901	50



Reusable Luer to Luer Transfer

Reusable stainless steel female Luer Lock to Luer Lock transfer system. Used to transfer tissue or emulsify between syringes.

PART NO.	SIZE	QTY/BOX
3-T-TLL2	2mm	1
3-T-TLL3	3mm	1
3-T-TLL4	4mm	1



Reusable Toomey to Luer Lock Syringe Adapter

Reusable stainless steel 60cc Toomey Monoject syringe hub to Luer Lock Syringe transfers tissue between syringes and is autoclavable.


PART NO.	QTY/BOX
3-S-T60-LL-SS	1

ACCESSORIES

Monoject™ Syringes


Monoject™ Syringes are sterile, individually packaged, single-use syringes without a needle. They have a polypropylene barrel and plunger rod with graduated barrel markings.

PART NO.	DESCRIPTION	VOLUME	QTY/BOX
3-S-T60-20	Toomey Tip	60cc	20



Syringe Racks


Syringe Racks are made of a unique polypropylene that is autoclavable up to 121°C (250°F) and maintain stability while submerged in water. Racks are easy to assemble and stackable.



Dimension of assembled:
 24.6 L x 10.4 W x 6.4 H (mm)
 9.7 L x 4.09 W x 2.52 H (in)

PART NO.	SIZE	HOLE	HOLE DIAMETER	FOR SYRINGE	QTY/BOX
3-S-R5	5cc	90	13mm	5cc, 3cc & 1cc BD	1
3-S-R10	10cc	60	16mm	10cc BD	1
3-S-R30	30cc	24	25mm	30cc BD, 20cc Monoject	1
3-S-R60	60cc	21	30mm	60cc BD, 35cc & 60cc Monoject Syringe	1


ACCESSORIES



Horseshoe Locks

Locks are used to maintain negative pressure in syringes during manual tissue aspiration. Horseshoe locks are aluminum and autoclavable.

PART NO.	DESCRIPTION	FOR SYRINGE	QTY/BOX
3-S-HL10	10cc	10cc BD Syringe	1
3-S-HL10M	10cc	10cc Monoject Syringe	1
3-S-HL20	20cc	20cc BD and Terumo	1
3-S-HL30	30cc	30cc BD Syringe	1
3-S-HL60	60cc	60cc BD or Monoject Syringe	1



Clicker Locks

Locks are used to maintain negative pressure in syringes during manual tissue aspiration. Clicker locks are aluminum with stainless steel springs and are autoclavable.


PART NO.	DESCRIPTION	FOR SYRINGE	QTY/BOX
3-S-CL10	10cc	10cc BD Syringe	1
3-S-CL10-M	10cc	10cc Monoject Syringe	1
3-S-CL20	20cc	20cc BD	1
3-S-CL30	30cc	300cc BD	1
3-S-CL35	35cc	35cc BD and Monoject Syringe	1
3-S-CL60	60cc	60cc BD and Monoject Syringe	1

ACCESSORIES

Syringe Power Adaptor 60cc

Syringe Power adaptor connects a 60cc syringe to tubing.


PART NO.	QTY/BOX
SPA60CC	1



Luer Lock Syringe Tip Caps

Sterile red polypropylene luer lock caps for standard luer lock syringes are individually packaged and disposable

PART NO.	QTY/BOX
3-S-65001	100




INFUSION HANDLES


Infiltration Handles

Autoclavable blue aluminum stainless steel handles with Luer Lock fittings on both ends.

PART NO.	SIZE
3-H-IHSL	Large 1" diameter
3-H-IHSS	Small ¾" diameter



3-H-IHSL




3-H-IHSS

INFUSION CANNULAS



GAN 6 Hole Garden Spray



BAN 6 Hole Body Spray




MH 6 Multi Hole Spray

Tumescent Infiltration Cannulas

Autoclavable stainless steel Luer Lock cannulas for infusion and infiltration.

PART NO.	DESCRIPTION
3-C-GAN	6 Hole Garden Spray
3-C-BAN	6 Hole Body Spiral
3-C-MH	6 Multi Hole Spray



Infiltration IV Spike

Connects an IV bag to infiltration tubing through spike connection.

PART NO.	QTY/BOX
IVSPIKE	1

DISPOSABLE CANNULAS

Luer Lock Cannulas

Disposable sterile hub with standard Luer Lock connection and Lipocyte friendly hole patterns used to remove fat from collection devices. 10/box.

One Hole Type II


CF301590	1 Hole II Injector, 3mmx15cm
CF201590	1 Hole II Injector, 2mmx15mm
CF10790	1 Hole II Injector, 1mmx7mm
CF16890	1 Hole II Injector, 16gaugex8cm
CF12790	1 Hole II Injector, 18gaugex7cm

Open-Ended

CPOE3300	OP, 3mmx30cm
-----------------	--------------

Multi-Hole

CP301524H	24 Hole, 3mmx15cm
CP141508H	8 Hole, 2.1mmx15cm
CP141506HS	6 Hole, Spiral, 2.1mmx15cm
CP30154H	4 Hole, 3mmx15cm



Handle and Cannula

A white disposable, lightweight sterile white handle with thumb groove. 6/box.

3 Hole

LC32303	3 Hole, 3mmx23cm
LC42803	3 Hole, 4mmx28cm
LC52803	3 Hole, 5mmx28cm

24 Hole

LC328024H	24 Hole, 3mmx28cm
LC428024H	24 Hole, 4mmx28cm

Description

3 Hole, 3mmx23cm
3 Hole, 4mmx28cm
3 Hole, 5mmx28cm

Description

24 Hole, 3mmx28cm
24 Hole, 4mmx28cm



3 Hole




24 Hole

WOUND DRAINAGE

SeromaCath® Wound Drainage Kits

Sterile wound drainage systems for post-operative seromas following mastectomy, axillary dissection, groin dissection, latissimus dorsi flap, TRAM flap, liposuction, abdominoplasty, breast implants and other procedures. Kit comes with connecting drainage tube for a bulb suction reservoir; bulb sold separately.

PART NO.	KIT CONTENTS	QTY/BOX
GR-1275	2.75" (7cm) Blunt for Seromas surrounding Breast implants 2" (5cm) 14gauge Teflon catheter/stylet needle/plug with two side holes 2.75" (7cm) 17gauge blunt metal catheter with two side holes Connecting drainage tube for a bulb suction reservoir	5 kits
GR-1035	3.5" (8.9cm) for Deep Seromas common with abdominoplasty 3.5" (8.9cm) 14gauge Teflon catheter/stylet needle/plug with six side holes Connecting drainage tube for a bulb suction reservoir	5 kits
GR-1002	2" (5cm) for Common Seromas 2" (5cm) 14gauge Teflon catheter/stylet needle/plug with six side holes Connecting drainage tube for a bulb suction reservoir	5 kits



WOUND DRAINAGE BULBS ARE NOT INCLUDED IN THE DRAINAGE KITS

VERONIQUE SIZE CHART

Female Size Chart

Our compression garments are latex-free and come in a variety of sizes and colors for both males and females. Garment stock is limited to availability. We recommend placing your order at least one week prior to the date the garment is needed. We do not offer returns on garment orders.

SIZE	WAIST	HIP	BUST	BREAST	UNDER BICEP CIRCUMFERENCE
XS	21-31"	31-33"	27-29"	25-26"	9.50-10.25"
Small	24-26"	34-36"	30-32"	27-30"	10.50-11.25"
Medium	27-29"	37-39"	33-36"	31-34"	11.50-12.25"
Large	30-32"	40-42"	37-40"	35-38"	12.50-13.25"
XL	33-35"	43-45"	41-44"	39-42"	13.50-14.25"
2XL	36-39"	45-48"	45-48"	43-46"	14.50-15.25"
3XL	40-43"	49-52"	49-52"	47-49"	15.75-16.75"
4XL	44-47"	53-56"	53-56"	50-52"	17.00-18.00"

NO RETURNS ON GARMENT ORDERS


*Garments are not guaranteed in stock. Please order one week prior to procedure.
All Veronique Garments are LATEX-FREE.*




BREASTS

Breast/Body Wraps

All of our breast wraps are made of nylon and spandex, offer adjustable downward compression, and have adjustable front hooks and loop clasps. They are one size fits all.



V-510-3
3" Breast Wrap



V-510-3U
3" Breast Wrap with
Cotton Lining




V-510-6
6" Breast Wrap



V-510-9
9" Breast Wrap


**NO RETURNS ON
GARMENT ORDERS**

Specialty Breast Wraps



V-510-J
3" Breast Wrap with Lace

Breast wrap with lace is made of nylon and spandex, offers adjustable downward compression, and has adjustable front hooks and loop clasps. They are one size fits all.



V-510-X
Versatile Stabilizing Breast Wrap

Versatile stabilizing breast wrap is made of nylon and spandex, offers customized band positioning with, and has adjustable front hooks and loop clasps. They are one size fits all.

FACE

Universal Facial Bands

Our universal facial bands have adjustable contouring face and neck bands and hook and loop clasps for a soft and durable stretch. They come in white and champagne and are available in S-XL.



V-210-2
2" Universal Facial Band



V-2100B
Occipital Universal Facial Band



V-210
4" Universal Facial Band



V-210S
Universal Split Ear Facial Band

Universal split ear facial band has adjustable and detachable anchor straps, split ear design, and hook and loop clasps for a soft and durable stretch. It is one size fits all.



V-210U
Universal Facial Band with Cotton Lining

Universal facial band has adjustable and detachable anchor straps and hook and loop clasps for a soft and durable stretch. It has a soft cotton lining and is one size fits all.



V-210K
Universal Facial Band with Cold/Hot Compress

Universal facial band with adjustable and detachable anchor straps as well as detachable cotton knit pouches for cold or hot compresses. Hook and loop clasps create a soft and durable stretch. It is one size fits all.



V-210N
Universal Facial Band

Universal facial band has adjustable contouring face and neck bands and hook and loop clasps for a soft and durable stretch. It is available in S-XL.

NO RETURNS ON GARMENT ORDERS


NASAL AND SINUS PACKING

Our nasal and sinus packing products are the perfect product for facial aesthetic procedures like rhinoplasties and septoplasties. They are sterile, expand to various sizes, and are specifically designed for comfortable insertion and fast removal.

AbsorbENT™ Nasal and Sinus Packs

Single-use nasal and sinus packs are made of sterile Polyvinyl Alcohol (PVA) foam that provides gentle but firm compression to nasal mucosa. They are lint and fiber free and have interconnected open cells that are soft and pliable for epistaxis management.

Packs expand quickly and evenly and have a string attached for retention and easy removal. Dimensions listed are for the hydrated state while inserted.




PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
RH-7530-10	Ultra Thin	3.5 x 1.2 x .6	10
RH-7532-10	Ultra Thin	3.5 x 1.2 x .9	10
RH-2315-20	Small	5.5 x 1.5 x 2	20
RH-7401-10	Small	4.5 x 1.5 x 2	10
RH-7408-10	Small with Airway	4.5 x 1.5 x 2	10
RH-7404-10	Large	8 x 1.5 x 2	10
RH-7405-10	Large with Airway	8 x 1.5 x 2	10
RH-7409-10	Large Anatomical	8 x 1.5 x 3	10
RH-7410-10	Large Anatomical with Airway	8 x 1.5 x 3	10
RH-7411-10	Large Anatomical - Thin	8 x 1 x 3	10
RH-7406-10	Large - Thin	8 x 1 x 2	10



Compressed



Hydrated




**Anatomical with
Airway**

NASAL AND SINUS PACKING

Rhino Rocket® Nasal Packs with Applicator

Rhino Rocket Nasal Packs are made of sterile PVA Expandacell® foam that provides epistaxis management after trauma, nasal or sinus surgery. Rhino Rockets provide gentle but firm mucosal compression with a tamponade of anterior or middle segment epistaxis. They take approximately 30 seconds to position and insert and swell to up to six times their compressed diameter when wet. Attached string for retention and easy removal. Dimensions listed are for hydrated state while inserted.



Aids in Insertion

Patented Applicator

Protects Foam from Moisture

Minimizes Blood Containment Concerns

Slimline

PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
11S-S0300-08AS	Small	3 x 1 x 2	8
11S-S0500-08AS	Medium	5 x 1 x 2	8
11S-S0800-08AS	Large	8 x 1 x 2	8
11S-S1000-08AS	Extra Large	10 x 1 x 2	8

Wideline

PART NO.	SIZE	DIMENSIONS LxWxH (CM)	QTY/BOX
11W-W0300-08AS	Small	3 x 1.5 x 2	8
11W-W0500-08AS	Medium	5 x 1.5 x 2	8
11W-W0800-08AS	Large	8 x 1.5 x 2	8
11W-W1000-08AS	Extra Large	10 x 1.5 x 2	8

Starter Pack


PART NO.	SIZE	DIMENSIONS	QTY/BOX
11S-STARTERRHINO	Variety pack for new users. One of each size.		8

EPISTAXIS BALLOONS

Epistaxis® Catheters

Double balloon catheter made of medical-grade silicone that controls intranasal bleeding. Single-use catheter is completely sterile and is available with or without a syringe.

PART NO.	SIZE	DESCRIPTION	QTY/BOX
RH-4200-01	97mm	No Syringe	1
RH-4200-10	97mm	No Syringe	10
RH-4220-01	97mm	With 20cc Syringe	1
RH-4400-01	102mm	No Syringe	1
RH-4400-10	102mm	No Syringe	10
RH-4420-01	102mm	With 20cc Syringe	1



INTERNAL NASAL SPLINTS

Breeze™ Nasal Airway Splint

The Breeze Nasal Airway Splint is a single-use product made of medical grade silicone that is commonly referred to as a “Doyle Splint.” Its unique airway geometry creates a pressure gradient force that increases airflow. It has suture indicators located at the anterior tip to aid in suture placement. The splint has a lubricious coating straight out of the package that aids insertion and removal.


PART NO.	DESCRIPTION	QTY/BOX
SP-79300	With Airway, Sterile	1 Pair
SP-79300-05	With Airway, Sterile	5 Pairs
SP-79300-50	With Airway, Sterile	50 Pairs
SP-79200-B	With Airway, Non-Sterile	50 Pairs
SP-79400	Without Airway, Sterile	1 Pair



Nasal Septal Button

Sterile, single-use Nasal Septal Buttons are made of medical grade silicone and can be easily trimmed to ensure proper fit using the 1 cm measurement rings. The Nasal Septal Button's one-piece design with 2 flange provides maximum security and cleanliness for nasal septal perforations. 6mm posts on all sizes.

PART NO.	DESCRIPTION	DIAMETER AND THICKNESS	QTY/BOX
SP-78100	Small	3cm x 4mm	1
SP-78105	Medium	5cm x 4mm	1
SP-78107	Large	7cm x 4mm	1



Bi-Valve Nasal Splint

Sterile, single-use Bi-Valve Nasal Splints are made of clear medical grade fluoroplastic. Their non-stick surface reduces adhesion and/or hematoma, has a medial slit that facilitates insertion and removal and holes for easy suturing.

Standard Size: 26.1 mm x 42.6 mm

Oversize: 39.11 mm x 52 mm

PART NO.	DESCRIPTION	QTY/BOX
SP-27000	Standard Thin (.25mm)	6 Pairs
SP-27005	Standard Thick (.5mm)	6 Pairs
SP-27010	Oversized Thin (.25mm)	6 Pairs
SP-27015	Oversized Thick (.5mm)	6 Pairs




DENVER SPLINTS®


The Denver Nasal Splint stabilizes the nose after surgery reducing edema after nasal surgery or trauma. The moldable aluminum retains nasal bones in position and multiple layers provide cushioned pressure on subcutaneous tissue. The dorsal pad reduces dead space and the possibility of polybeak. The Denver Splints come in four sizes: mini, petite, small/medium, large and are sold standalone or in kits.

Denver Splint® Series 1500

The Denver Splint Series 1500 is made of soft moldable aluminum and Velcro™ and can be worn for 7-10 days.



PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-1500-05KM	Kit, Mini	50mm x 28mm	5
10-1500-05KP	Kit, Petite	50mm x 30mm	5
10-1500-05KS	Kit, Small/Medium	65mm x 35mm	5
10-1500-05KL	Kit, Large	69mm x 43mm	5
10-1500-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-1500-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-1500-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-1500-10NKL	Non-Kit, Large	69mm x 43mm	10




- Adhesive
- Moldable Aluminum
- No Water Mess
- Inconspicuous

DENVER SPLINTS®

Denver Splint® Series 2000 - One Piece

The Denver Splint Series 2000 is made of soft moldable aluminum and 1/8" of ventilating foam and can be worn for 5-7 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-2000-05KM	Kit, Mini	50mm x 28mm	5
10-2000-05KP	Kit, Petite	50mm x 30mm	5
10-2000-05KS	Kit, Small/Medium	65mm x 35mm	5
10-2000-05KL	Kit, Large	69mm x 43mm	5
10-2000-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-2000-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-2000-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-2000-10NKL	Non-Kit, Large	69mm x 43mm	10



Denver Splint® Series 5000

The Denver Splint Series 5000 is made of soft moldable aluminum and Velcro™ and has an attached dorsal pad. It can be worn for 7-10 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-5000-05KM	Kit, Mini	50mm x 28mm	5
10-5000-05KP	Kit, Petite	50mm x 30mm	5
10-5000-05KS	Kit, Small/Medium	65mm x 35mm	5
10-5000-05KL	Kit, Large	69mm x 43mm	5
10-5000-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-5000-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-5000-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-5000-10NKL	Non-Kit, Large	69mm x 43mm	10


Denver Splint® Series 1800

The Denver Splint Series 1800 is made of soft moldable aluminum and Velfoam. It can be worn for 7-10 days.

PART NO.	DESCRIPTION	DIMENSIONS	QTY/BOX
10-1800-05KM	Kit, Mini	50mm x 28mm	5
10-1800-05KP	Kit, Petite	50mm x 30mm	5
10-1800-05KS	Kit, Small/Medium	65mm x 35mm	5
10-1800-05KL	Kit, Large	69mm x 43mm	5
10-1800-10NKM	Non-Kit, Mini	50mm x 28mm	10
10-1800-10NKP	Non-Kit, Petite	50mm x 30mm	10
10-1800-10NKS	Non-Kit, Small/Medium	65mm x 35mm	10
10-1800-10NKL	Non-Kit, Large	69mm x 43mm	10

KIT INCLUDES:

- 1 Denver Splint®
- 2 Pre-cut Protecto Strips, sterile
- 1 Alcohol wipe, sterile
- 1 Skin prep, sterile
- 1 Dorsal pad



NON-KIT INCLUDES:

- 1 Denver Splint®
- 1 Dorsal pad

THERMOPLASTIC NASAL SPLINTS




Pre-Cut Denver Splint® Series 4000

- Rigid, opaque, lightweight material
- Soften with hot water and mold over taped nose

PART NO.	SIZE	DESCRIPTION	QTY/BOX
10-4021-05KL	Large	Kit, Ivory 1/16"	5
10-4021-05KM	Mini	Kit, Ivory 1/16"	5
10-4021-05KP	Petite	Kit, Ivory 1/16"	5
10-4021-05KS	Small	Kit, Ivory 1/16"	5
10-4021-10NKL	Large	Non-Kit, Ivory 1/16"	10
10-4021-10NKM	Mini	Non-Kit, Ivory 1/16"	10
10-4021-10NKP	Petite	Non-Kit, Ivory 1/16"	10
10-4021-10NKS	Small	Non-Kit, Ivory 1/16"	10
10-4022-05KL	Large	Kit, Blush 1/16"	5
10-4022-05KM	Mini	Kit, Blush 1/16"	5
10-4022-05KP	Petite	Kit, Blush 1/16"	5
10-4022-05KS	Small	Kit, Blush 1/16"	5
10-4022-10NKL	Large	Non-Kit, Blush 1/16"	10
10-4022-10NKM	Mini	Non-Kit, Blush 1/16"	10
10-4022-10NKP	Petite	Non-Kit, Blush 1/16"	10
10-4022-10NKS	Small	Non-Kit, Blush 1/16"	10
10-4031-05KL	Large	Kit, Ivory 3/32"	5
10-4031-05KM	Mini	Kit, Ivory 3/32"	5
10-4031-05KP	Petite	Kit, Ivory 3/32"	5
10-4031-05KS	Small	Kit, Ivory 3/32"	5
10-4031-10NKL	Large	Non-Kit, Ivory 3/32"	10
10-4031-10NKM	Mini	Non-Kit, Ivory 3/32"	10
10-4031-10NKP	Petite	Non-Kit, Ivory 3/32"	10
10-4031-10NKS	Small	Non-Kit, Ivory 3/32"	10

KIT INCLUDES:

- 1 Denver Splint®
- 2 Pre-cut Protecto Strips, sterile
- 1 Alcohol wipe, sterile
- 1 Skin prep, sterile
- 1 Dorsal pad



NON-KIT INCLUDES:

- 1 Denver Splint®
- 1 Dorsal pad

INDICATE SIZE



THERMOPLASTIC NASAL SPLINTS

Thermoplast Splints


Thermoplast splints are made of rigid but lightweight material that softens with hot water allowing them to be affixed to the patient. Splints are available in pre-cut sizes (3"x3") or in sheets to cut your own shape (9"x12").

3x3" Squares

PART NO.	DESCRIPTION	QTY/BOX
AQ-F-05-21	Ivory 1/16"	5
AQ-F-10-21	Ivory 1/16"	10
AQ-F-05-31	Ivory 3/32"	5
AQ-F-10-31	Ivory 3/32"	10
AQ-F-05-22	Blush 1/16"	5
AQ-F-10-22	Blush 1/16"	10

9x12" Sheets

PART NO.	DESCRIPTION	QTY/BOX
AQ-9x12-21	Ivory 1/16"	1
AQ-9x12-31	Ivory 3/32"	1
AQ-9x12-22	Blush 1/16"	1



SINUS PACK AND SPLINTING CROSS REFERENCE

All nasal and sinus packs are made from sterile, polyvinyl alcohol (PVA) foam. This chart offers easy cross referencing and ordering.

OUR PART NO.	COMPETITORS PART NO.	DESCRIPTION
NASAL M-TYPE		
RH-7401-10	440400/400400	M-Type Small
RH-7408-10	440409/400409	M-Type Small with Airway
RH-7404-10	440402/400402	M-Type Large
RH-7405-10	440410/400410	M-Type Large with Airway
NASAL P-TYPE		
RH-2310-10	440406 (Pope)/400406	P-Type XL
RH-2315-20	400421 (Weimert)	P-Type Small
NASAL D-TYPE		
RH-7409-10	440424 (Doyle)/400424	D-Type Wideline Anatomical
RH-7411-10	440413 (Doyle)/400413	D-Type Slimline Anatomical
SINUS NASAL K & T TYPES		
RH-7530-10	400530 (Kennedy)	K-Type, Ultra-Slim
RH-7532-10	400428 (Toffel)	T-Type, Slim
EAR PACKS		
OT-5201-50	400141	Ear Speedi-Wick 15mm
OT-5202-50	400146	Ear Speedi-Wick 25mm
30-305	425008	Ear Fenestrated 15mm
OT-5204-10	400355	Ear Post-Op Small
OT-5203-10	400162	Ear Post-Op Large
DENVER SPLINTS		
10-2000-05KL	1528136	Denver Splint, 2000 Kit Large
10-2000-05KS	1528126	Denver Splint, 2000 Kit Sm/Med
10-2000-05KP	1528116	Denver Splint, 2000 Kit Petite
10-2000-05KM	1528116	Denver Splint, 2000 Kit Mini
10-2000-10NKL	1528137/1528138	Denver Splint, 2000 Non-Kit Large
10-2000-10NKS	1528120/1528121	Denver Splint, 2000 Non-Kit Sm/Med
10-2000-05NKP	1528110/1528111	Denver Splint, 2000 Non-Kit Petite
10-2000-05NKM	1528110/1528111	Denver Splint, 2000 Non-Kit Mini
THERMOPLAST SPLINTS		
10-4021-05KL	1529120	Denver Splint, 4000 Kit Large
10-4021-05KS	1529110	Denver Splint, 4000 Kit Sm/Med
10-4021-05KP	1529100	Denver Splint, 4000 Kit Petite
10-4021-10NKS	1529010	Denver Splint, 4000 Non-Kit Sm/Med
10-4021-10NKP	1529000	Denver Splint, 4000 Non-Kit Petite
10-4022-10NKL	1529020	Denver Splint, 4000 Non-Kit Large

**All of our packs have an attached string*

Shippert Aesthetics
Aesthetic & Reconstructive

Summit Medical LLC
An Innovia Medical Company

815 Vikings Parkway
Suite 100
St. Paul, MN 55121 | USA

1.888.229.2875
www.innoviamedical.com/aesthetic-reconstructive

Summit Medical LLC is an FDA registered facility, ISO 13485:2016 certified and MDSAP certified.

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.

Specifications and designs are subject to change without any notice or obligation on the part of the manufacturer.

Denver Splint, Tissu-Trans and Rhino Rocket are all registered trademarks of Shippert Medical. AbsorbENT and Breeze are all trademarks of Summit Medical.

All brands are trademarks or registered trademarks of their respective owners.

©2022 Innovia Medical. All rights reserved. M-81 | RevC | OCT2022

